

eXTReme™ Seal

Extra strong adhesive minimizes evaporation without residue

Application:

PCR and robotics


Temperature Range:

-20°C to +120°C

CATALOG #	DESCRIPTION	QUANTITY
XTR-100	eXTReme Seal	100

eXTReme Seal features an extra-strong thick adhesive to minimize evaporation, yet removes easily without leaving residue on the plate. Features chamfered corners to eliminate plate overhang for robotics applications and tight storage conditions.

- Certified DNase- & RNase- free
- Ideal for use with raised-rim and non-skirted PCR plates
- Available in a larger cut for flat-top semi-skirted and full-skirted PCR plates: eXTReme Seal, Large



"For a cut compatible with flat-top PCR plates, check out eXTReme Seal, Large"