

eXTreme™ Seal, Large

Extra strong adhesive minimizes evaporation without residue

Application:

PCR

Temperature Range:



-20°C to +120°C

CATALOG #	DESCRIPTION	QUANTITY
XTR-LG-100	eXTreme Seal, Large	100

eXTreme Seal, Large features an extra strong, thick adhesive to minimize evaporation, yet removes easily without leaving residue on the plate; ideal when plates need to be sealed multiple times. The larger cut is designed for greater plate coverage for a more complete seal around the edges to reduce evaporation rates.

- Certified DNase- & RNase- free
- Ideal for use with flat-top, semi-skirted and full-skirted PCR plates as well as deep-well plates
- Large cut allows for easier positioning of the film to the plate via edge of plate alignment vs. centering

Available in a smaller cut for raised-rim and non-skirted PCR plates: eXTreme Seal



"For a smaller cut compatible with raised-rim & non-skirted plates, check out eXTreme Seal"