


PR1MA™ Pipette Tips

Technical Information



Pipette Tips shown at actual size.

Tip Dimensions

Volume	Length (mm)	Top - Inner Dia (mm)
10 µL	31.1	3.8
10 µL XL	45.8	3.8
20 µL Barrier	50.8	5.3
100 µL Barrier	50.8	5.3
200 µL	50.8	5.3
200 µL XL	91.1	6.1
300 µL	55.9	5.8
1000 µL	76.1	7.5
1250 µL	90.2	7.4

PR1MA™ Pipette Tips

PR1MA Pipette Tip – Description & Packaging

PR1MA Part No	Volume	Description	Sterile?	Contents per CS	Purity Statement
PR-10BK	10 µL	Bulk Tips		1,000 tips	Certified RNase, DNase, DNA and Endotoxin Free
PR-10RFL	10 µL	Reloads		10 reloads of 96 tips	Certified RNase, DNase, DNA and Endotoxin Free
PR-10RK-NS	10 µL	Racked Tips		10 racks of 96 tips	Certified RNase, DNase, DNA and Endotoxin Free
PR-10RK-S	10 µL	Racked Sterile Tips	Y	10 racks of 96 tips	Certified RNase, DNase, DNA, Endotoxin and Pyrogen-Free
PR-10RK-FL	10 µL Barrier	Racked Sterile Barrier Tips	Y	10 racks of 96 tips	Certified RNase, DNase, DNA, Endotoxin and Pyrogen-Free
PR-10XLBK	10 µL XL	Bulk Tips		1,000 tips	Certified RNase, DNase, DNA and Endotoxin Free
PR-10XRFL	10 µL XL	Reloads		10 reloads of 96 tips	Certified RNase, DNase, DNA and Endotoxin Free
PR-10XRK-NS	10 µL XL	Racked Tips		10 racks of 96 tips	Certified RNase, DNase, DNA and Endotoxin Free
PR-10XRK-S	10 µL XL	Racked Sterile Tips	Y	10 racks of 96 tips	Certified RNase, DNase, DNA, Endotoxin and Pyrogen-Free
PR-10XRK-FL	10 µL XL Barrier	Racked Sterile Barrier Tips	Y	10 racks of 96 tips	Certified RNase, DNase, DNA, Endotoxin and Pyrogen-Free
PR-200BK	200 µL	Bulk Tips		1,000 tips	Certified RNase, DNase, DNA and Endotoxin Free
PR-200RFL	200 µL	Reloads		10 reloads of 96 tips	Certified RNase, DNase, DNA and Endotoxin Free
PR-200RK-NS	200 µL	Racked Tips		10 racks of 96 tips	Certified RNase, DNase, DNA and Endotoxin Free
PR-20RK-FL	20 µL Barrier	Racked Sterile Barrier Tips	Y	10 racks of 96 tips	Certified RNase, DNase, DNA, Endotoxin and Pyrogen-Free
PR-100RK-FL	100 µL Barrier	Racked Sterile Barrier Tips	Y	10 racks of 96 tips	Certified RNase, DNase, DNA, Endotoxin and Pyrogen-Free
PR-200RK-S	200 µL	Racked Sterile Tips	Y	10 racks of 96 tips	Certified RNase, DNase, DNA, Endotoxin and Pyrogen-Free
PR-200RK-FL	200 µL Barrier	Racked Sterile Barrier Tips	Y	10 racks of 96 tips	Certified RNase, DNase, DNA, Endotoxin and Pyrogen-Free
PR-200XRK-NS	200 µL XL	Racked Tips		10 racks of 96 tips	Certified RNase, DNase, DNA and Endotoxin Free
PR-200XRK-S	200 µL XL	Racked Sterile Tips	Y	10 racks of 96 tips	Certified RNase, DNase, DNA, Endotoxin and Pyrogen-Free
PR-200XRK-FL	200 µL XL Barrier	Racked Sterile Barrier Tips	Y	10 racks of 96 tips	Certified RNase, DNase, DNA, Endotoxin and Pyrogen-Free
PR-300BK	300 µL	Bulk Tips		1,000 tips	Certified RNase, DNase, DNA and Endotoxin Free
PR-300RFL	300 µL	Reloads		10 reloads of 96 tips	Certified RNase, DNase, DNA and Endotoxin Free
PR-300RK-NS	300 µL	Racked Tips		10 racks of 96 tips	Certified RNase, DNase, DNA and Endotoxin Free
PR-300RK-S	300 µL	Racked Sterile Tips	Y	10 racks of 96 tips	Certified RNase, DNase, DNA, Endotoxin and Pyrogen-Free
PR-300RK-FL	300 µL Barrier	Racked Sterile Barrier Tips	Y	10 racks of 96 tips	Certified RNase, DNase, DNA, Endotoxin and Pyrogen-Free
PR-1000BK	1000 µL	Bulk Tips		1,000 tips	Certified RNase, DNase, DNA and Endotoxin Free
PR-1000RK-NS	1000 µL	Racked Tips		8 racks of 96 tips	Certified RNase, DNase, DNA and Endotoxin Free
PR-1000RK-S	1000 µL	Racked Sterile Tips	Y	8 racks of 96 tips	Certified RNase, DNase, DNA, Endotoxin and Pyrogen-Free
PR-1000RK-FL	1000 µL Barrier	Racked Sterile Barrier Tips	Y	8 racks of 96 tips	Certified RNase, DNase, DNA, Endotoxin and Pyrogen-Free
PR-1250BK	1250 µL	Bulk Tips		1,000 tips	Certified RNase, DNase, DNA and Endotoxin Free
PR-1250RFL	1250 µL	Reloads		8 reloads of 96 tips	Certified RNase, DNase, DNA and Endotoxin Free
PR-1250RK-NS	1250 µL	Racked Tips		8 racks of 96 tips	Certified RNase, DNase, DNA and Endotoxin Free
PR-1250RK-S	1250 µL	Racked Sterile Tips	Y	8 racks of 96 tips	Certified RNase, DNase, DNA, Endotoxin and Pyrogen-Free
PR-1250RK-FL	1250 µL Barrier	Racked Sterile Barrier Tips	Y	8 racks of 96 tips	Certified RNase, DNase, DNA, Endotoxin and Pyrogen-Free