

MidSCI™

2025

QUICK SHOP

guide



TABLE OF CONTENTS

3 •• Liquid Handling	50 •• Hoods
11 •• Centrifuge Tubes	51 •• Spectrophotometers
14 •• Labels	52 •• Incubators
18 •• Storage Racks	54 •• PCR Equipment
19 •• PCR/qPCR Supplies	55 •• Swabs & Saliva Kits
24 •• PCR/qPCR Reagents	56 •• Electrophoresis
27 •• Chemicals	58 •• Transilluminators
28 •• Cold Storage Supplies	59 •• Imagers
31 •• Reservoirs	60 •• Homogenizers
32 •• Petri Dishes & Spreaders	61 •• Baths
33 •• Erlenmeyer Flasks	62 •• Cold Storage
34 •• Filtration	67 •• Water System
35 •• Tissue Cell Culture	68 •• Balances
40 •• Mixing	69 •• Microscopes
46 •• Centrifuges	70 •• General Lab Supplies

Stay connected with MIDSCI for the latest updates, industry insights, and exciting content!

Follow us on LinkedIn, Facebook, and YouTube to join our growing community.



[linkedin.com/company/midsci](https://www.linkedin.com/company/midsci)



[FB.me/midsci](https://www.facebook.com/midsci)



[youtube.com/@midsci5374](https://www.youtube.com/@midsci5374)

Don't miss out on valuable information, tips, and news? Sign up for our emails.



[midsci.com/contact-us](https://www.midsci.com/contact-us)



If you can't find it here, give us a call at 800-227-9997 or email us at tech@midsci.com

MIDSCI assumes no liability for printing errors, mistakes, or product changes. All prices subject to change without notice.


All terms and conditions apply. Visit our website for more information: <https://www.midsci.com/terms>

V.01.2025

LIQUID HANDLING

PR1MA™ Low-Retention Pipette Tips

- All PR1MA™ tips are low-retention!
- Don't lose your sample in the tip
- Easier identification with color coded refills & clear racks
- Fits all standard pipettors
- Certified free of detectable RNase, DNase, PCR inhibitors, endotoxins and tested pyrogen free.



- Barrier tips have extra working volume, so liquid doesn't touch the filter
- Great stackable racks
- 10 µL XL tips have 1 µL graduations
- Clear & readable graduations

PR1MA™ Low-Retention Barrier Pipette Tips



ITEM #	DESCRIPTION	QTY	PRICE
PR-10RK-FL	10 µL, Sterile	96/rk, 10 rks/cs	\$68.17
PR-10XLRK-FL	10 µL XL, Sterile	96/rk, 10 rks/cs	\$67.97
PR-20RK-FL	20 µL, Sterile	96/rk, 10 rks/cs	\$67.97
PR-100RK-FL	100 µL, Sterile	96/rk, 10 rks/cs	\$67.97
PR-200RK-FL	200 µL, Sterile	96/rk, 10 rks/cs	\$67.97
PR-200XLRK-FL	200 µL XL 5 in. Long, Sterile	96/rk, 8 rks/cs	\$85.78
PR-300RK-FL	300 µL, Sterile	96/rk, 10 rks/cs	\$67.97
PR-1000RK-FL	1,000 µL, Sterile	96/rk, 8 rks/cs	\$59.40
PR-1250RK-FL	1,250 µL, Sterile	96/rk, 8 rks/cs	\$56.25



LIQUID HANDLING

PR1MA™ Sterile Racked Pipette Tips



ITEM #	DESCRIPTION	QTY	PRICE
PR-10RK-S	10 µL, Sterile	96/rk, 10 rks/cs	\$42.33
PR-10XLRK-S	10 µL XL, Sterile	96/rk, 10 rks/cs	\$42.33
PR-200RK-S	200 µL, Sterile	96/rk, 10 rks/cs	\$44.40
PR-200XLRK-S	200 µL XL 5 in. Long, Sterile	96/rk, 8 rks/cs	\$62.21
PR-300RK-S	300 µL, Sterile	96/rk, 10 rks/cs	\$54.63
PR-1000RK-S	1,000 µL, Sterile	96/rk, 8 rks/cs	\$54.94
PR-1250RK-S	1,250 µL, Sterile	96/rk, 8 rks/cs	\$46.96



PR1MA™ Non-Sterile Racked Pipette Tips

ITEM #	DESCRIPTION	QTY	PRICE
PR-10RK-NS	10 µL, Non-Sterile	96/rk, 10 rks/cs	\$37.05
PR-10XLRK-NS	10 µL XL, Non-Sterile	96/rk, 10 rks/cs	\$37.05
PR-200RK-NS	200 µL, Non-Sterile	96/rk, 10 rks/cs	\$39.18
PR-200XLRK-NS	200 µL XL 5 in. Long, Non-Sterile	96/rk, 8 rks/cs	\$56.10
PR-300RK-NS	300 µL, Non-Sterile	96/rk, 10 rks/cs	\$49.01
PR-1000RK-NS	1,000 µL, Non-Sterile	96/rk, 8 rks/cs	\$49.27
PR-1250RK-NS	1,250 µL, Non-Sterile	96/rk, 8 rks/cs	\$43.69

PR1MA™ Pipette Tip Reloads

ITEM #	DESCRIPTION	QTY	PRICE
PR-10RFL	10 µL, Non-Sterile	96/rk, 10 rks/cs	\$28.27
PR-10XLRFL	10 µL XL, Non-Sterile	96/rk, 10 rks/cs	\$28.27
PR-200RFL	200 µL, Non-Sterile	96/rk, 10 rks/cs	\$26.99
PR-300RFL	300 µL, Non-Sterile	96/rk, 10 rks/cs	\$30.03
PR-1250RFL	1,250 µL, Non-Sterile	96/rk, 8 rks/cs	\$37.16

Racks only, no boxes included



LIQUID HANDLING



PR1MA™ Bulk Pipette Tips



ITEM #	DESCRIPTION	QTY	PRICE
PR-10BK	10 µL, Non-Sterile	1,000/pk	\$17.45
PR-10XLBK	10 µL XL, Non-Sterile	1,000/pk	\$17.45
PR-200BK	200 µL, Non-Sterile	1,000/pk	\$13.57
PR-200XLBK	200 µL XL 5 in. Long, Non-Sterile	1,000/pk	\$44.54
PR-300BK	300 µL, Non-Sterile	1,000/pk	\$21.88
PR-1000BK	1,000 µL, Non-Sterile	1,000/pk	\$22.50
PR-1250BK	1,250 µL, Non-Sterile	1,000/pk	\$21.43

- Round or flat tip options
- Short body 200 µL available in bulk packaging
- Patented press-to-open racks provide easy, one-handed opening
- Universal fit to accommodate most pipettor brands
- Multichannel-friendly with a standard 96-tip arrangement
- Certified free of detectable RNase, DNase, DNA, PCR inhibitors, and tested pyrogen-free. Sterile product is certified endotoxin-free





Gel-Loading Tips



ITEM #	DESCRIPTION	QTY	PRICE
4817-S0	10 µL, 0.57 mm, Round, Racked, Sterile	96/rk, 10 rks/cs	\$187.43
4817-SF	10 µL, 0.57 mm, Round, Filtered, Racked, Sterile	96/rk, 10 rks/cs	\$224.99
4837-00	200 µL, 0.57 mm, Round, Racked, Non-Sterile	96/rk, 10 rks/cs	\$171.30
4837-S0	200 µL, 0.57 mm, Round, Sterile, Racked	96/rk, 10 rks/cs	\$187.43
4847-AF	20 µL, 0.57 mm, Short-Body, Rnd (MBP/QSP Style), Racked, Sterile	96/rk, 10 rks/cs	\$224.99
4847-BF	100 µL, 0.57 mm, Short-Body, Rnd (MBP/QSP Style), Racked, Sterile	96/rk, 10 rks/cs	\$224.99
4847-00	200 µL, 0.57 mm, Short-Body, Rnd (MBP/QSP Style), Racked, Non-Sterile	96/rk, 10 rks/cs	\$127.93
4847-S0	200 µL, 0.57 mm, Short-Body, Rnd (MBP/QSP Style), Racked, Sterile	2 rks/pk, 5 pks/cs	\$152.20
4840-00	200 µL, 0.57 mm, Short-Body, Bulk (MBP/QSP Style), Bulk, Non-Sterile	1,000/bg, 10 bgs/cs	\$708.40

LIQUID HANDLING

- Compatible with Rainin® single- and multi-channel LiteTouch® pipettors
- A light force is required for adding or ejecting tips, making them easy to use
- LTS tips are made from medical-grade virgin polypropylene
- Racks and grids are made from our exclusive post-consumer recycled polypropylene, which reduces the use of virgin polypropylene by up to 61%



PR1MA™ LTS® Compatible Pipette Tips

ITEM #	DESCRIPTION	QTY	PRICE
PR-LTS20RFL	20 µL LTS-Style Tips, Reloads	96/pk, 10pks/cs	\$32.03
PR-LTS20RFL-S	20 µL LTS-Style Tips, Reloads, Pre-Sterile	96/pk, 10pks/cs	\$32.03
PR-LTS20RK-NS	20 µL LTS-Style Tips, Racked	96/pk, 10pks/cs	\$41.98
PR-LTS20RK-S	20 µL LTS-Style Tips, Racked, Sterile	96/pk, 10pks/cs	\$47.96
PR-LTS20RK-FL	20 µL LTS-Style Filter Tips, Racked, Sterile	96/pk, 10pks/cs	\$77.01
PR-LTS200RFL	200 µL LTS-Style Tips, Reloads	96/pk, 10pks/cs	\$30.58
PR-LTS200RFL-S	200 µL LTS-Style Tips, Reloads, Pre-Sterile	96/pk, 10pks/cs	\$33.92
PR-LTS200RK-NS	200 µL LTS-Style Tips, Racked	96/pk, 10pks/cs	\$44.39
PR-LTS200RK-S	200 µL LTS-Style Tips, Racked, Sterile	96/pk, 10pks/cs	\$50.31
PR-LTS200RK-FL	200 µL LTS-Style Filter Tips, Racked, Sterile	96/pk, 10pks/cs	\$77.01
PR-LTS300RFL	300 µL LTS-Style Tips, Reloads	96/pk, 10pks/cs	\$32.03
PR-LTS300RFL-S	300 µL LTS-Style Tips, Reloads, Pre-Sterile	96/pk, 10pks/cs	\$35.38
PR-LTS300RK-NS	300 µL LTS-Style Tips, Racked	96/pk, 10pks/cs	\$44.39
PR-LTS300RK-S	300 µL LTS-Style Tips, Racked, Sterile	96/pk, 10pks/cs	\$50.31
PR-LTS300RK-FL	300 µL LTS-Style Filter Tips, Racked, Sterile	96/pk, 10pks/cs	\$77.01
PR-LTS1000RFL	1,000 µL XL LTS-Style Tips, Reloads	96/pk, 8pks/cs	\$42.10
PR-LTS1000RFL-S	1,000 µL XL LTS-Style Tips, Reloads, Pre-Sterile	96/pk, 8pks/cs	\$48.79
PR-LTS1000RK-NS	1,000 µL XL LTS-Style Tips, Racked	96/pk, 8pks/cs	\$49.49
PR-LTS1000RK-S	1,000 µL XL LTS-Style Tips, Racked, Sterile	96/pk, 8pks/cs	\$53.21
PR-LTS1000RK-FL	1,000 µL XL LTS-Style Filter Tips, Racked, Sterile	96/pk, 8pks/cs	\$63.73

Certified free of detectable RNase, DNase, DNA, PCR inhibitors, and tested pyrogen-free. Sterile product is certified endotoxin-free.

Rainin, LTS and LiteTouch are registered trademarks of Mettler-Toledo Rainin, LLC.

LIQUID HANDLING

AlphaPette™ by PR1MA™ Single-Channel Pipettors

- Universal shaft will accept most pipette tips
- Lower ejection force
- Autoclavable
- Distinguishing color sets it apart for dedicated use
- Lifetime warranty (*except calibration and abuse*)

ITEM #	DESCRIPTION	QTY	PRICE
A-2	0.1 to 2 µL	1/ea	\$280.83
A-10	0.5 to 10 µL	1/ea	\$280.83
A-20	2 to 20 µL	1/ea	\$280.83
A-100	10 to 100 µL	1/ea	\$280.83
A-200	20 to 200 µL	1/ea	\$280.83
A-1000	100 to 1,000 µL	1/ea	\$280.83
A-5000	1,000 to 5,000 µL	1/ea	\$297.07
P4405	Multi-Channel Pipette Racks, 5-Place	1/ea	\$93.26
MID4280-CR	Carousel Pipettor Stand, 6-Place	1/ea	\$101.21



First calibration is free
with included certificate!

- Universal shaft will accept most pipette tips
- Lower ejection force
- 3-year warranty



Alphapette™ by PR1MA™ Multi-Channel Pipettors

ITEM #	DESCRIPTION	QTY	PRICE
AP4508-10	8-Channel, 1 to 10 µL	1/ea	\$710.01
AP4508-50	8-Channel, 5 to 50 µL	1/ea	\$710.01
AP4508-200	8-Channel, 20 to 200 µL	1/ea	\$710.01
AP4508-300	8-Channel, 50 to 300 µL	1/ea	\$710.01
AP4512-10	12-Channel, 1 to 10 µL	1/ea	\$864.95
AP4512-50	12-Channel, 5 to 50 µL	1/ea	\$864.95
AP4512-200	12-Channel, 20 to 200 µL	1/ea	\$864.95
AP4512-300	12-Channel, 50 to 300 µL	1/ea	\$864.95

Halo™ Single-Channel Pipettors

ITEM #	DESCRIPTION	QTY	PRICE
H6800-2	0.1 to 2.0µL	1/ea	\$255.97
H6800-10	1 to 10 µL	1/ea	\$255.97
H6800-20	2 to 20 µL	1/ea	\$255.97
H6800-50	5 to 50 µl	1/ea	\$255.97
H6800-100	10 to 100 µL	1/ea	\$255.97
H6800-200	20 to 200 µL	1/ea	\$255.97
H6800-1M	100 to 1,000 µL	1/ea	\$255.97
H6800-1MXL	200 to 1,250µL	1/ea	\$255.97
H6800-5M	0.5 to 5 mL	1/ea	\$255.97
H6800-10M	1 to 10 mL	1/ea	\$255.97
H6800-CR	Carousel for Halo Pipettors	1/ea	\$92.01
H6800-SK	Halo Pipettor Starter Kit	1/ea	\$718.57



Universal
shaft will
accept most
pipette tips!

LIQUID HANDLING

- Universal shaft will accept most pipette tips
- Multifunctional with six operational modes
- 1-year warranty



Excel™ Electronic Pipettors

ITEM #	DESCRIPTION	QTY	PRICE
P3600L-10	Single-Channel, 0.5 to 10 µL	1/ea	\$406.20
P3600L-20	Single-Channel, 2 to 20 µL	1/ea	\$406.20
P3600L-200	Single-Channel, 10 to 200 µL	1/ea	\$406.20
P3600L-1200	Single-Channel, 100 to 1,200 µL	1/ea	\$406.20
P3608L-10	8-Channel, 0.5 to 10 µL	1/ea	\$829.44
P3608L-20	8-Channel, 2 to 20 µL	1/ea	\$829.44
P3608L-200	8-Channel, 10 to 200 µL	1/ea	\$829.44
P3608L-1200	8-Channel, 100 to 1200 µL	1/ea	\$940.27
P3612L-10	12-Channel, 0.5 to 10 µL	1/ea	\$968.18
P3612L-20	12-Channel, 2 to 20 µL	1/ea	\$968.18
P3612L-200	12-Channel, 10 to 200 µL	1/ea	\$968.18
P3612L-1200	12-Channel, 100 to 1,200 µL	1/ea	\$1,106.31
P3600L-BAT-1	Lithium Ion Battery	1/ea	\$48.11

The PR1MA™ Starter Kit offers a complete liquid handling solution, including four precision pipettes that cover a wide volume range from 0.5 µl to 1,000 µl, making it ideal for a variety of applications.

It also features three racks of PR1MA Low-Retention Pipette Tips, designed to minimize sample loss and ensure accurate, reproducible results.

For convenience, the kit includes a six-position carousel stand for organized storage and a PR1MA Backpack for easy transport.

Additionally, the PR1MA Starter Kit provides over \$250 in savings compared to purchasing the individual components separately.



PR1MA™ Starter Kit

ITEM #	DESCRIPTION	QTY	PRICE
PR-LQH01-SK	PR1MA™ Starter Kit, Includes: Backpack, A-10, A-20, A-200, A-1000, 3rks PR1MA Tips, 6-pl Carousel*	1/ea	\$997.00

*No mix and match or switching pipettor or tip sizes. PR1MA™ Low Retention Graduated Pipette Tips are certified free of detectable RNase, DNase, DNA, PCR inhibitors, and tested pyrogen-free.

LIQUID HANDLING

- Made of high clarity, 100% virgin polystyrene
- Sterilized by E-beam, SAL=10⁻⁶ and individually wrapped in paper/plastic packaging
- Non-pyrogenic, DNase/Rnase free
- Edged color-coded band for identification
- Polyolefin fiber filter to reduce contamination

New & Improved!

PR1MA™ Serological Pipettes



ITEM #	DESCRIPTION	QTY	PRICE
PR-SERO-1	1 mL, Individ. Wrapped, Sterile	800/cs	\$120.44
PR-SERO-2	2 mL, Individ. Wrapped, Sterile	600/cs	\$136.46
PR-SERO-5	5 mL, Individ. Wrapped, Sterile	200/cs	\$74.19
PR-SERO-10	10 mL, Individ. Wrapped, Sterile	200/cs	\$80.93
PR-SERO-25	25 mL, Individ. Wrapped, Sterile	200/cs	\$99.12
PR-SERO-50	50 mL, Individ. Wrapped, Sterile	90/cs	\$99.78

Also available in sterile, bulk packaging!

Great for high volume pipetting!



- Manufactured with the highest grade polystyrene, extremely break resistant and robust
- Reservoir offers large capacity volumes in a compact, ergonomic design

TPP® Serological Pipettes

ITEM #	DESCRIPTION	QTY	PRICE
TP94525	25/40 mL, Individ. Wrapped, Sterile	60/cs	\$146.52
TP94550	50/70 mL, Individ. Wrapped, Sterile	50/cs	\$170.62

- **Turbo-Speed Functionality:** The Turbo mode delivers an impressive 25% increase in pipetting speed, ideal for handling large-volume pipettes with reservoirs
- Fits all commercially available glass and plastic serological pipettes ranging from 1 to 100 mL
- Adjustable thumb wheel enables you to set the desired pipetting speed

TPP® Turbo-Fix Pipet Aid

ITEM #	DESCRIPTION	QTY	PRICE
TP94701	TPP® Turbo-Fix Pipet Aid	1/ea	\$1,027.57



LIQUID HANDLING

- Featuring the QuickStand™ Retractable Kickstand, supports pipette at an angle and folds away when not in use
- Rechargeable battery with 8 hours of continuous use on a single charge
- Ergonomic design, lightweight and comfortable for extended use
- Variable speed control, adjustable speed for precise pipetting
- Pipette compatibility, fits pipettes from 1 mL to 100 mL



PR1MA™ Pipette Controllers

ITEM #	DESCRIPTION	QTY	PRICE
MPIP-1	PR1MA™ Pipette Controller, 115V	1/ea	\$326.80
MID6950	PR1MA™ Reach Long Neck, 115V	1/ea	\$379.85

- Ergonomic Comfort: Designed to reduce hand and arm fatigue during extended use
- Precise Control: Multiple speed settings for accurate liquid transfer
- Rechargeable Power: Can be charged during use for continuous operation
- Safety Features: Double filters prevent contamination; includes extra filters and a stand
- Quiet Performance: Advanced, low-noise pump for efficient operation



Drummond Pipet-Aid Pipette Controllers

ITEM #	DESCRIPTION	QTY	PRICE
4-098-101	XP, White and Yellow, 110V	1/ea	\$385.00
4-098-501	XP2, White and Yellow, 110V	1/ea	\$385.00
4-098-105	XL, White and Yellow, 110V	1/ea	\$385.00


- Ergonomic, lightweight design
- Continuously adjustable pipetting speed via thumb-wheel
- Can be charged while in use
- Overfilling and contamination prevented by hydrophobic 0.45 mm filter
- Fits 1 to 100 mL serological pipettes
- 1-year warranty

PIPETBOY acu 2 Pipet Aid

ITEM #	DESCRIPTION	QTY	PRICE
155-000	1 to 100 mL, 0.45 µm Filter, 100 - 240 VAC, Grey	1/ea	\$512.50
155-015	1 to 100 mL, 0.45 µm Filter, 100 - 240 VAC, Green	1/ea	\$512.50
155-016	1 to 100 mL, 0.45 µm Filter, 100 - 240 VAC, Red	1/ea	\$512.50
155-017	1 to 100 mL, 0.45 µm Filter, 100 - 240 VAC, Blue	1/ea	\$512.50
155-019	1 to 100 mL, 0.45 µm Filter, 100 - 240 VAC, Grey	1/ea	\$512.50



CENTRIFUGE TUBES



PR1MA™ Premium Microcentrifuge Tubes

ITEM #	DESCRIPTION	QTY	PRICE
PR-MCT6	0.5 mL, Snap Cap, Boil Proof, Clear, Non-Sterile	1,000/pk	\$33.07
PR-MCT6-LR	0.5 mL, Low-Retention, Snap Cap, Boil Proof, Clear, Non-Sterile	500/pk	\$24.49
PR-MCT6-S	0.5 mL, Snap Cap, Boil Proof, Clear, Sterile	500/pk	\$25.49
PR-MCT15	1.5 mL, Snap Cap, Boil Proof, Clear, Non-Sterile	500/pk	\$19.61
PR-MCT15-LR	1.5 mL, Low-Retention, Snap Cap, Boil Proof, Clear, Non-Sterile	250/pk	\$13.99
PR-MCT15-S	1.5 mL, Snap Cap, Boil Proof, Clear, Sterile	250/pk	\$14.72
PR-MCT17	1.5 mL, Large Snap Cap, Boil Resist, Clear, Non-Sterile	500/pk	\$19.61
PR-MCT2	2 mL, Snap Cap, Boil Proof, Clear, Non-Sterile	500/pk	\$17.75
PR-MCT2-LR	2 mL, Low-Retention, Snap Cap, Boil Proof, Clear, Non-Sterile	250/pk	\$17.75
PR-MCT2-S	2 mL, Snap Cap, Boil Proof, Clear, Sterile	250/pk	\$14.93

For colors add:

-B (blue), -G (green), -O (orange), -R (red), -P (purple), -Y (yellow), -RB (assorted), or -RBF (assorted neon) to the Item #
Sterile available for clear tubes only. Low-retention available for clear, non-sterile tubes only.

LARGE CAP: BOIL-RESISTANT, EASY TO OPEN AND CLOSE, AND A LARGE LABELING SURFACE

SMALL CAP: BOIL-PROOF, IDEAL FOR RIGOROUS HIGH-TEMP APPLICATIONS

OPTICAL CLARITY ENSURES MAXIMUM SAMPLE VISIBILITY

RATED TO 20,000 XG

GRADUATIONS EVERY 100 µL

AUTOCCLAVABLE

LARGE WRITING SURFACE ON CAP TOP



Lots more to choose from on our website!

*Certified free of detectable RNase, DNase, DNA, PCR inhibitors and tested pyrogen-free.
 Sterile product is also certified endotoxin-free.*

CENTRIFUGE TUBES

- Ideal for aliquots, histology samples, & so much more
- Racks, inserts for your centrifuge, and all the accessories you need
- Graduations at every 0.2 mL and 0.5 mL

*Great alternative
when 15 mL is too much
& 2 mL is not enough*



5 mL Microcentrifuge Tubes

ITEM #	DESCRIPTION	QTY	PRICE
C1005-SC5-B	5 mL [PP], Screw Cap, Non-Sterile	100/bg, 5 bgs/cs	\$101.77
C1005-SC5-R	5 mL [PP], Screw Cap, Racked, Sterile	50/rk, 10 rks/cs	\$115.50

See page 46 for the perfect centrifuge!



- Fits most 50 mL rotors - Rim eliminates need for adapter
- Printed graduations - From 5 mL to 2 mL
- Compact low-profile design - Minimizes shipping, storage and waste
- Gamma sterilized - RNase, DNase, & Pyrogen free

25 mL Centrifuge Tubes

ITEM #	DESCRIPTION	QTY	PRICE
C2625-B	25 mL [PP], Screw Cap, Sterile, Bagged	25/bg, 8 bgs/cs	\$104.03
C2625-R	25 mL [PP], Screw Cap, Sterile, Racked	25/pk, 8 pks/cs	\$107.58

- Features a 'skirted' bottom that allows the tubes to stand upright without a separate test tube rack
- Sterile, RNase, DNase, pyrogen, and cytotoxin free
- Compatible with standard rotors
- Frosted marking surface

30 & 50 mL Self-Standing Centrifuge Tubes

ITEM #	DESCRIPTION	QTY	PRICE
MID-2630	30 mL [PP], Flat Screw Cap, Sterile, Bagged	25/bg, 20 bgs/cs	\$173.78
C2605	50 mL [PP], Flat Screw Cap, Sterile, Bagged	25/bg, 20 bgs/cs	\$211.61

**NEW! 15 and 50 mL
Eco-Racked, now available!**



CENTRIFUGE TUBES

New & Improved!



- Printed graduations with white marking area
- Polypropylene tube and high density polyethylene cap
- High centrifugal strength up to 12,000 xg
- Flat markable cap with leakproof seal
- Sterilized by E-beam, SAL=10⁻⁶
- Non-pyrogenic, DNase/Rnase free

PR1MA™ Conical Centrifuge Tubes

ITEM #	DESCRIPTION	QTY	PRICE
C15B	15 mL [PP], Screw Cap, Bagged, Sterile	500/cs	\$155.82
C15R	15 mL [PP], Screw Cap, Racked, Sterile	500/cs	\$142.62
C50B	50 mL [PP], Screw Cap, Bagged, Sterile	500/cs	\$187.67
C50R	50 mL [PP], Screw Cap, Racked, Sterile	500/cs	\$187.03

- Outstanding for high centrifugation or ultra low freezing
- Liquid & temp resistant label on leak-proof screw cap
- Specs:
[PP] -190°C to 100°C & 15,500 xg
[PS] -30°C to +70°C & 1,700 xg
- White graduations down to 0.1 mL & 0.5 mL respectively



TPP® Centrifuge Tubes

ITEM #	DESCRIPTION	QTY	PRICE
TP91016	13 mL, Round [PP], Sterile	40/pk, 20 pks/cs	\$779.23
TP91019	13 mL, Skirted Conical [PP], Sterile	40/pk, 20 pks/cs	\$779.91
TP91017	15 mL, Round [PP], Sterile	40/pk, 20 pks/cs	\$808.16
TP91015	15 mL, Conical, [PP], Sterile	40/pk, 20 pks/cs	\$678.13
TP91115	15 mL, Conical, [PS], Sterile	40/pk, 20 pks/cs	\$827.31
TP91050	50 mL, Conical, [PP], Sterile	20/pk, 18 pks/cs	\$365.50
TP91051	50 mL, Skirted Conical, [PP], Sterile	20/pk, 16 pks/cs	\$352.25
TP91056	50 mL, Round, [PP], Sterile	20/pk, 18 pks/cs	\$238.93



PR1MA™ Screw Cap Tubes & Cap Assemblies



2 mL conical

2 mL skirted

ITEM #	DESCRIPTION	QUANTITY	PRICE
PR-SC5AC1	0.5 mL, Conical, Clear, Sterile*	50/bg, 10 bgs/cs	\$109.68
PR-SC5AS1	0.5 mL, Skirted, Clear, Sterile*	50/bg, 10 bgs/cs	\$109.68
PR-SC15AC1G	1.5 mL, Conical, Graduated, Clear, Sterile*	50/bg, 10 bgs/cs	\$102.05
PR-SC15AS1	1.5 mL, Skirted, Clear, Sterile*	50/bg, 10 bgs/cs	\$102.05
PR-SC2AC1G	2 mL, Conical, Graduated, Clear, Sterile*	50/bg, 10 bgs/cs	\$102.05
PR-SC2AS1G	2 mL, Skirted, Graduated, Clear, Sterile*	50/bg, 10 bgs/cs	\$102.05

*Screw cap included (PR-SCC)

LABELS

PRIMA-DOTS Rolls & Sheets

NEW!



ROLLS

ITEM #	DESCRIPTION	QTY	PRICE
PR-DR95-W	3/8 in. Dots for 0.5-1.5 mL Tubes, White	1,000/pk	\$23.00
PR-DR95-RB	3/8 in. Dots for 0.5-1.5 mL Tubes, Asst	5,000/pk	\$60.00
PR-DR13-W	1/2 in. Dots for 1.5-2 mL Tubes, White	1,000/pk	\$25.00
PR-DR13-RB	1/2 in. Dots for 1.5-2 mL Tubes, Asst	5,000/pk	\$80.00
PR-DR19-W	3/4 in. Dots for 15-50 mL Tubes, White	1,000/pk	\$30.00
PR-DR25-W	1 in. Dots for 50 mL Tubes, White	500/pk	\$30.00

SHEETS

PR-DS95-W	3/8 in. Dots for 0.5-1.5 mL Tubes, White	192 labels/sheet, 20 sheets/pk	\$65.00
PR-DS95-RB	3/8 in. Dots for 0.5-1.5 mL Tubes, Asst	960 labels/color, 20 sheets/pk	\$99.00
PR-DS13-W	1/2 in. Dots for 1.5-2 mL Tubes, White	192 labels/sheet, 20 sheets/pk	\$65.00
PR-DS13-RB	1/2 in. Dots for 1.5-2 mL Tubes, Asst	960 labels/color, 20 sheets/pk	\$80.00



PRIMA-TAGS Rolls & Sheets

NEW!

ITEM #	DESCRIPTION	QTY	PRICE
PR-TR217-W	0.83 x 0.28 in. for 0.5 mL Tubes, White	1,500/pk	\$25.00
PR-TR2413-W	0.95 x 0.5 in. for 0.5 mL Tubes, White	1,000/pk	\$25.00
PR-TR3313-W	1.3 x 0.5 in. for 1.5/2 mL Tubes, White	1,000/pk	\$25.00
PR-TR3819-W	1.5 x 0.75 in. for General Use, White	1,000/pk	\$25.00

ROLLS

PR-TS238-W	0.92 x 0.32 in. for 0.2 mL Tubes, White	119 labels/sheet, 20 sheets/pk	\$66.00
PR-TS2413-W	0.95 x 0.5 in. for 0.5 mL Tubes, White	119 labels/sheet, 20 sheets/pk	\$65.00
PR-TS3313-W	1.3 x 0.5 in. for 1.5-2 mL Tubes, White	85 labels/sheet, 20 sheets/pk	\$65.00
PR-TS386-W	1.5 x 0.24 in. for Microplates, White	156 labels/sheet, 20 sheets/pk	\$65.00
PR-TS3819-W	1.5 x 0.75 in. for General Use, White	60 labels/sheet, 20 sheets/pk	\$65.00
PR-TS4319-W	1.69 x 0.75 in. for Vials, White	52 labels/sheet, 20 sheets/pk	\$65.00
PR-TS6725-W	2.64 x 0.98 in. for Racks/Boxes, White	30 labels/sheet, 20 sheets/pk	\$65.00
PR-TS7651-W	2.99 x 2.01 in. for Racks/Boxes, White	12 Labels/sheet, 20 sheets/pk	\$65.00
PR-TS2525-W	0.98 x 0.98 in. for Microscope Slides, White	63 labels/sheet, 20 sheets/pk	\$65.00

SHEETS




Looking for other colors or varieties?

Scan the code for more options!

LABELS

NEW!

PR1MA™-CRYODOTS Rolls & Sheets



ROLLS

ITEM #	DESCRIPTION	QTY	PRICE
PR-CDR95-W	3/8 in. Dots for 0.5-1.5 mL Tubes, White	1,000/pk	\$25.00
PR-CDR95-RB	3/8 in. Dots for 0.5-1.5 mL Tubes, Asst	5,000/pk	\$120.00
PR-CDR13-W	1/2 in. Dots for 1.5-2 mL Tubes, White	1,000/pk	\$29.00
PR-CDR13-RB	1/2 in. Dots for 1.5-2 mL Tubes, Asst	5,000/pk	\$80.00
PR-CDR19-W	3/4 in. Dots for 15-50 mL Tubes, White	1,000/pk	\$35.00
PR-CDR25-W	1 in. Dots for 50 mL Tubes, White	500/pk	\$35.00

SHEETS

PR-CDS95-W	3/8 in. Dots for 0.5-1.5 mL Tubes, White	192 labels/sheet, 20 sheets/pk	\$65.00
PR-CDS13-W	1/2 in. Dots for 1.5-2 mL Tubes, White	192 labels/sheet, 20 sheets/pk	\$65.00

Labels Temp Range -196°C to 90°C



Looking for other colors or varieties?

Scan the code for more options!

PR1MA™-CRYOTAGS Rolls & Sheets

NEW!

ITEM #	DESCRIPTION	QTY	PRICE
PR-CTR217-W	0.83 x 0.28 in. for 0.5 mL Tubes, White	1,500/pk	\$29.00
PR-CTR2413-W	0.95 x 0.5 in. for 0.5 mL Tubes, White	1,000/pk	\$29.00
PR-CTR3313-W	1.3 x 0.5 in. for 1.5/2 mL Tubes, White	1,000/pk	\$29.00
PR-CTR3319-W	1.5 x 0.75 in. for General Use, White	1,000/pk	\$31.00
PR-CTSPCR-W	0.92 x 0.32 in. for 0.2 mL Tubes, White	119 labels/sheet, 20 sheets/pk	\$66.00
PR-CTS2413-W	0.95 x 0.5 in. for 0.5 mL Tubes, White	119 labels/sheet, 20 sheets/pk	\$65.00
PR-CTS3313-W	1.3 x 0.5 in. for 1.5-2 mL Tubes, White	85 labels/sheet, 20 sheets/pk	\$65.00
PR-CTS386-W	1.5 x 0.24 in. for Microplates, White	156 labels/sheet, 20 sheets/pk	\$65.00
PR-CTS3819-W	1.5 x 0.75 in. for General Use, White	60 labels/sheet, 20 sheets/pk	\$65.00
PR-CTS4319-W	1.69 x 0.75 in. for Vials, White	52 labels/sheet, 20 sheets/pk	\$65.00
PR-CTS6725-W	2.64 x 0.98 in. for Racks/Boxes, White	30 labels/sheet, 20 sheets/pk	\$65.00
PR-CTS7651-W	2.99 x 2.01 in. for Racks/Boxes, White	12 Labels/sheet, 20 sheets/pk	\$65.00
PR-CTS2525-W	0.98 x 0.98 in. for Microscope Slides, White	63 labels/sheet, 20 sheets/pk	\$65.00

ROLLS

SHEETS

Labels Temp Range -196°C to 90°C

ORDER TODAY!

midsci.com • custserv@midsci.com • 800.227.9997 • f: 636.225.9998

LABELS

- Tough Spots® are circular, PVC labels designed for handwriting with most lab markers

- Available in a variety of colors and sizes to suit different labeling needs

- Ideal for use in laboratory environments requiring durable, easy-to-write labels

Tough-Spots® Rolls & Sheets



-40°C to 121°C

-196°C to 90°C

ITEM #	DESCRIPTION	QTY	PRICE
SPRL-1000	3/8 in. for 1.5-2.0 mL Tubes/Vials, White	5,000/pk	\$111.25
SPRL-1100	3/8 in. for 1.5-2.0 mL Tubes/Vials, Asst	5,000/pk	\$111.25
SPRL-2000	1/2 in. for 2.0 mL Tubes/Vials, White	5,000/pk	\$111.25
SPRL-2100	1/2 in. for 2.0 mL Tubes/Vials, Assorted	5,000/pk	\$111.25
SPOT-1000	3/8 in. for 0.5 mL Tubes, White	192 labels/sheet, 20 sheets/pk	\$67.00
SPOT-2000	1/2 in. for 1.5 mL and 2.0 mL Tubes, White	192 labels/sheet, 20 sheets/pk	\$67.00
SPOT-3000	3/4 in. for 15 mL Tubes, White	108 labels/sheet, 20 sheets/pk	\$67.00
SPOT-4000	1 in. for 50 mL Tubes, White	63 labels/sheet, 20 sheets/pk	\$67.00
SPOT-5000	7/16 in. for 1.5 mL and 2.0 mL Tubes, White	234 labels/sheet, 20 sheets/pk	\$67.00

ROLLS

SHEETS

- Tough-Tags® are durable, chemically inert laboratory labels that adhere well to microtubes, test tubes, storage boxes, and microscope slides

- They are resistant to organic solvents, caustic agents, and are ideal for use with samples stored in moderate freezer temperatures

- Available in two formats: one for handwriting and one for laser printing

Tough-Tags® Rolls & Sheets

ITEM #	DESCRIPTION	QTY	PRICE
TTSC-1000	0.94 x 0.5 in. for 0.5-1.5 mL Tubes, Clear	1,000/pk	\$33.50
TTSW-1000	0.94 x 0.5 in. for 0.5-1.5 mL Tubes, White	1,000/pk	\$28.00
TTSH-1000	0.94 x 0.5 in. for "Cut to Size", White	25 sheets/pk	\$68.00
TTSW-2240	0.94 x 0.5 in. for 0.5-1.5 mL Tubes, White	119 labels/sheet, 25 sheets/pk	\$72.75

SHEETS — ROLLS

-40°C to 121°C



Looking for other colors or varieties?
Scan the code for more options!

LABELS

- Cryo-Babies® and Cryo-Tags® are designed for long-term cryogenic storage and conventional freezer conditions
- Labels feature solvent-resistant adhesive, ensuring they adhere to a wide range of surfaces, including cryogenic vials, microtubes, test tubes, glass vials, beakers, and metal racks
- Heat and cold resistant, they withstand flash freezing in liquid nitrogen and dry heat up to 80°C



Cryo-Babies On A Roll

ITEM #	DESCRIPTION	QTY	PRICE
CRY-0500	0.94 x 0.5 in. Babies for Cryo Vials, White	1,000/pk	\$34.50
CRY-1000	1.28 x 0.5 in. Babies for Cryo Vials, White	1,000/pk	\$34.50
CRY-2000	1.5 x 0.75 in. Babies for Cryo Vials, White	1,000/pk	\$34.50

-196°C to 150°C

ROLLS



Laser Cryo Tags & Babies

ITEM #	DESCRIPTION	QTY	PRICE
LCRY-1100	1.69 x 0.75 in. Tags for Large Tubes/Boxes/Racks, White	52 labels/sheet, 20 sheets/pk	\$67.00
LCRY-1200	1.50 x 0.75 in. Tags for Large Tubes/Vials, White	60 labels/sheet, 20 sheets/pk	\$67.00
LCRY-1258	2.625 x 1.0 in. Tags for Large Tubes/Boxes/Racks, White	30 labels/sheet, 20 sheets/pk	\$67.00
LCRY-1300	1.125 x 1.75 in. Tags for Large Tubes/Boxes/Racks, White	42 labels/sheet, 20 sheets/pk	\$67.00
LCRY-1400	1.0 x 1.0 in. Babies for Microscope Slides, White	80 labels/sheet, 20 sheets/pk	\$67.00
LCRY-1700	1.28 x 0.50 in. Babies for Cryogenic Vials, White	85 labels/sheet, 20 sheets/pk	\$67.00
LCRY-2380	0.94 x 0.50 in. Babies for Cryogenic Vials, White	119 labels/sheet, 20 sheets/pk	\$67.00
LCRY-2500	2.50 x 0.50 in. Tags for Large Tubes/Boxes/Racks, White	60 labels/sheet, 20 sheets/pk	\$67.00
LCRY-3020	3.0 x 2.0 in. Tags for Racks/Bottles/Boxes, White	12 labels/sheet, 20 sheets/pk	\$67.00


SHEETS

-196°C to 150°C




STORAGE RACKS

IsoFreeze® PCR Color-Change Racks



ITEM #	DESCRIPTION	QTY	PRICE
5640-T4	96-Well, Purple to Pink, Cold, With lid	2/pk	\$94.83
5640-T6	96-Well, Green to Yellow, Cold, With lid	2/pk	\$94.83
5650-T4	96-Well, Purple to Pink, Cold, No lid	2/pk	\$99.94
5650-T6	96-Well, Green to Yellow, Cold, No lid	2/pk	\$99.94



4-Way Flipper® Rack

ITEM #	DESCRIPTION	QTY	PRICE
FL-4N	0.5 to 50 mL 4-way flipper rack, Natural	5/pk	\$38.17
FL-4A	0.5 to 50 mL 4-way flipper rack, Assorted	5/pk	\$41.57



See more storage rack options online!

Utility Racks



ITEM #	DESCRIPTION	QTY	PRICE
WRI-4A	8-Place PCR WorkUp Rack & Workstation w/ Lid, Assorted Fluorescent	5/pk	\$128.71
WR-96A	96-Place PCR WorkUp Rack, w/ Lid, Assorted Fluorescent	5/pk	\$39.86
WUR-96A	96-Place Stacking PCR WorkUp Rack, Assorted Standard	5/pk	\$42.99
HS29025A	80-Place Storage Rack, Natural	5/pk	\$17.94
HS29025G	80-Place Storage Rack, Assorted Fluorescent	5/pk	\$20.10



HS29025G



WR-96A

- Multi-tube compatibility: 12 x 0.5 mL; 9 x 1.5-2.0 mL; 6 x 15 mL; 2 x 50 mL
- Each Flipper can be independently oriented
- 4-way racks can be used in the tray or independently



4-Way Flipper® Racks with Tray

ITEM #	DESCRIPTION	QTY	PRICE
5400-29	4-way Flipper Rack w/ Tray, Assorted Fluorescent	4 rks/tray	\$60.80
5400-39	4-way Flipper Rack w/ Tray, Assorted Standard	4 rks/tray	\$57.69



PCR/qPCR SUPPLIES



PR1MA™ PCR Strips with Attached Caps

ITEM #	DESCRIPTION	QTY	PRICE
PR-PCR28ACD	0.2 mL x 8 Strip PCR Tubes, Clear w/ Individ. Attached Dome Cap	120/cs	\$69.90
PR-PCR28ACF	0.2 mL x 8 Strip PCR Tubes, Clear w/ Individ. Attached Flat Cap	120/cs	\$69.90
PR-PCR28ACFF	0.2 mL x 8 Strip "Easy Split" PCR Tubes, Clear w/ Individ. Attached Flat Cap	120/cs	\$69.90
PR-PCR28ACZD	0.2 mL x 8 Strip PCR Tubes, Clear w/ Hinge Attached Dome Cap	125/cs	\$82.96
PR-PCR28ACZF	0.2 mL x 8 Strip PCR Tubes, Clear w/ Hinge Attached Flat Cap	125/cs	\$72.68

PR1MA™ PCR Strips with Caps



ITEM #	DESCRIPTION	QTY	PRICE
PR-PCR18CF	0.1 mL x 8 PCR Strip Tubes & Flat Caps	125/cs	\$77.97
PR-PCR28CF	0.2 mL x 8 PCR Strip Tubes & Flat Caps	125/cs	\$77.97
PR-PCR28CD	0.2 mL x 8 PCR Strip Tubes & Domed Caps	125/cs	\$77.97
PR-PCR212CF	0.2 mL x 12 PCR Strip Tubes & Flat Caps	80/cs	\$92.47
PR-PCR212CD	0.2 mL x 12 PCR Strip Tubes & Domed Caps	80/cs	\$92.47

Need white or colors?
Give us a call!



Strip Caps work with our
PR1MA™ PCR Plates!



PR1MA™ PCR Strips & Caps

ITEM #	DESCRIPTION	QTY	PRICE
PR-PCR18	0.1 mL x 8 PCR Strip Tubes	125/cs	\$71.47
PR-PCR28	0.2 mL x 8 PCR Strip Tubes	125/cs	\$72.40
PR-PCRC8F	PCR Strip of 8 Flat Caps	125/cs	\$18.88
PR-PCRC8D	PCR Strip of 8 Domed Caps	125/cs	\$19.93
PR-PCR212	0.2 mL x 12 PCR Strip Tubes	80/cs	\$82.61
PR-PCRC12F	PCR Strip of 12 Flat Caps	80/cs	\$19.12
PR-PCRC12D	PCR Strip of 12 Domed Caps	80/cs	\$20.16

Strip and Caps Sold Separately

All tubes and caps on this page are certified free of detectable RNase, DNase, PCR inhibitors, & pyrogen free.

PCR/qPCR SUPPLIES

- Individual PCR tubes feature evaporation-resistant, easy-to-open-and-close lids that prevent sample loss
- Ultra-thin, uniform wells ensure optimal heat transfer and reaction efficiency
- Dome, frosted flat, or optically-clear caps for qPCR (real time)
- Conveniently packaged in reclosable bags to minimize risk of contamination
- Available in ultra-clear, assorted, or specific colors

PR1MA™ Individual PCR Tubes

ITEM #	DESCRIPTION	QTY	PRICE
PR-PCR2F	0.2 mL Individual PCR Tube w/Clear Flat Cap	1,000/cs	\$39.38
PR-PCR2FF	0.2 mL Individual PCR Tube w/Frosted Flat Cap	1,000/cs	\$39.38
PR-PCR2D	0.2 mL Individual PCR Tube w/Domed Cap	1,000/cs	\$46.00

Tubes and caps are certified free of detectable RNase, DNase, PCR inhibitors, & pyrogen free.



Not sure what will fit your equipment? **We have a solution for that!**
Check out our **PR1MA™** Plates and Tubes Selection Tool



- Uniformed wells and ultra-thin walls provide optimal heat transfer and reaction efficiency
- Fits most thermal cyclers
- High contrast black alphanumeric for easy well identification

PR1MA™ Full-Skirted Plates



ITEM #	DESCRIPTION	QTY	PRICE
PR-PCR3196	0.1 mL x 96 Full-Skirted Low Profile PCR Plate A12 Cut	10/cs	\$25.18
PR-PCR33841	384 Full-Skirted ABI Type 1-Notch PCR Plate A24 Cut	10/cs	\$36.75
PR-PCR33842	384 Full-Skirted Roche Type 2-Notch PCR Plate A24/P24 Cut	10/cs	\$30.89



Need white or colors?
Give us a call!

PCR/qPCR SUPPLIES

PR1MA™ Semi-Skirted Plates



ITEM #	DESCRIPTION	QTY	PRICE
PR-PCR2196F	0.1 mL x 96 Semi-Skirted Fast Type PCR Plate A1 Cut	10/cs	\$25.18
PR-PCR2196LC	0.1 mL x 96 Semi-Skirted LightCycler Type PCR Plate H12 Cut	10/cs	\$25.18
PR-PCR2296	0.2 mL x 96 Semi-Skirted ABI Type PCR Plate A12 Cut	10/cs	\$30.30
PR-PCR2296R	0.2 mL x 96 Semi-Skirted ABI Type Raised Rim PCR Plate A12 Cut	10/cs	\$25.18
PR-PCR2296S	0.2 mL x 96 Semi-Skirted Straight Side PCR Plate A12 Cut	10/cs	\$25.18



PR-PCR2196F



PR-PCR2196LC



PR-PCR2296



PR-PCR2296R



PR-PCR2296S

- High-contrast, black alphanumerics for easy well identification


- Ultra thin for uniformity & high reaction efficiency
- Well ridges reduce sample evaporation

PR1MA™ Non-Skirted Plates




**Broad thermal
cycler compatibility**


ITEM #	DESCRIPTION	QTY	PRICE
PR-PCR1196	0.1 mL x 96 Non-Skirted PCR Plate, H12 Cut	20/cs	\$60.59
PR-PCR1296	0.2 mL x 96 Non-Skirted PCR Plate, A12 Cut	10/cs	\$25.18
PR-PCR1296-CUT	0.2 mL x 96 Cuttable Non-Skirted PCR Plate, A12 Cut	10/cs	\$29.70
PR-PCR1C296	0.2 mL x 96 Non-Skirted Chimney Top PCR Plate, H12 Cut	10/cs	\$29.72



PR-PCR1196



PR-PCR1296 &
PR-PCR1296-CUT



PR-PCR1C296

PCR/qPCR SUPPLIES

Discover the next generation of PCR plates!

- Polycarbonate frame is rigid and flat to facilitate a firm grip by robotic arms and ensure reliable placement
- Polypropylene wells provide optimal thermal transfer and higher reaction efficiency
- Reduced outer well evaporation for more consistent data across all wells
- White well options enhance the fluorescence signal during qPCR
- Wells are slightly raised to ensure evaporation-proof sealing with film, foil, or strip caps
- ANSI/SBS footprint and well spacing



PRIMA™-Shell sets a new standard in performance and convenience.

PRIMA™-Shell PCR Plates



ITEM #	DESCRIPTION	QTY	PRICE
PR-PCR3742-C	0.1 mL x 96 Full-Skirt, Clear, Clear Wells	10/pk	\$58.77
PR-PCR3742-RB**	0.1 mL x 96 Full-Skirt, Assorted, Clear Wells	10/pk	\$51.10
PR-PCR3742-CW	0.1 mL x 96 Full-Skirt, Clear, White Wells	10/pk	\$58.77
PR-PCR3742BC-C	0.1 mL x 96 Full-Skirt, Clear, Clear Wells, Barcoded	10/pk	\$58.03
PR-PCR3750-C	0.2 mL x 96 Semi-Skirt, Clear, Clear Wells	10/pk	\$51.10
PR-PCR3750-RB**	0.2 mL x 96 Semi-Skirt, Assorted, Clear Wells	10/pk	\$51.10
PR-PCR3735-C	384 Full-Skirt, Clear, A24 Cut, Clear Wells	10/pk	\$75.39
PR-PCR3735-RB**	384 Full-Skirt, Assorted, A24 Cut, Clear Wells	10/pk	\$71.80
PR-PCR3735BC-C	384 Full-Skirt, Clear, A24 Cut, Clear Wells, Barcoded	10/pk	\$78.73
PR-PCR3730-C	384 Full-Skirt, Clear, A24 & P24 Cut, Clear Wells	10/pk	\$71.80
PR-PCR3730BC-C	384 Full-Skirt, Clear, A24 & P24 Cut, Clear Wells Barcoded	10/pk	\$78.73

All plates are made with rigid polycarbonate frames and polypropylene wells Certified free of detectable RNase, DNase, DNA, PCR inhibitors, and tested pyrogen-free.

For colors add: -B (blue), -G (green), -O (orange), -R (red), -V (violet), -Y (yellow) to the Item #
 **Assorted colors include Blue, Green, Red, Violet, and Yellow. Available with barcode. Custom barcoding is available.

Trust MIDSCI™ for PCR consumables that deliver accurate and reproducible results, empowering you to make groundbreaking discoveries.

PCR/qPCR SUPPLIES


PR1MA™ Sealing Mat

ITEM #	DESCRIPTION	QTY	PRICE
PR-PCRMAT	4.67 x 2.9 in. Silicone Pierceable Sealing Mat for 96-Well Plates	5/cs	\$52.99



- Application: PCR
- Temperature Range: -20°C to +120°C*
40°C to +120°C**
- Certified DNase & RNase free
- Ideal for use with flat-top, semi-skirted and full-skirted PCR plates

PCR Transparent Sealing Film



ITEM #	DESCRIPTION	QTY	PRICE
XTR2-100	eXTReme™ 2 Seal* Extra Strong Sealing Film, 2 mil Clear Polyester, Indiv. Sheets, Non-Sterile	100/bx	\$87.65
XTR2-LG-100	eXTReme™ 2 Seal* Extra Strong Sealing Film, 2 mil Clear Polyester, LG Indiv. Sheets, Non-Sterile	100/bx	\$101.10
100-THER-PLT	ThermalSeal® Sealing Film**, Easy to Remove Film for PCR & General Storage, Indiv. Sheets, Non-Sterile	100/bx	\$112.48
STR-THER-PLT	ThermalSeal® Sealing Film**, Easy to Remove Film for PCR & General Storage, Indiv. Sheets, Sterile	100/bx	\$147.89
TS2-100	ThermalSeal® 2™** Soft, Extra Thick, Conformable Adhesive, Indiv. Sheets, Non-Sterile	100/bx	\$123.76
TSA-100	ThermalSeal® Advanced***, Soft, Extra Thick, Conformable Adhesive, SM Indiv. Sheets, Non-Sterile	100/bx	\$124.09

ThermalSeal® RTS™ & RT2™ qPCR Sealing Film

ITEM #	DESCRIPTION	QTY	PRICE
TS-RT2-100	-40°C to +120°C, Optically Transparent, Highly Conformable Adhesive, Indiv. Sheets, Non-Sterile	100/bx	\$223.80
TS-RT2RR-100	-40°C to +120°C, Optically Transparent, Highly Conformable Adhesive, SM Indiv. Sheets, Non-Sterile, for Raised Rim Plates	100/bx	\$223.80



- Certified DNase & RNase free
- Excellent chemical resistance (e.g. DMSO)
- Applications: PCR, Long-term storage, Light sensitive assays, Robotics



AlumaSeal® II™ PCR Aluminum Film

ITEM #	DESCRIPTION	QTY	PRICE
AF-100	-80°C to +120°C, Easy-to-Pierce Aluminum Foil Seal, Indiv. Sheets, Non-Sterile	100/bx	\$138.59
AFS-25	-80°C to +120°C, Easy-to-Pierce Aluminum Foil Seal, Indiv. Sheets, Sterile	50/bx	\$115.47

PCR/qPCR REAGENTS

Ideal for routine PCR applications as well as genotyping, colony PCR and fast PCR

PR1MA™ Taq



ITEM #	DESCRIPTION	QTY	PRICE
PR1000-500	Taq	500 units	\$108.84
PR1000-1000	Taq	1,000 units	\$175.56
PR1000-5000	Taq	5,000 units	\$902.28
PR1000-10000	Taq	10,000 units	\$1,494.74
PR1000-R-500	Taq Red Dye	500 units	\$144.69
PR1001-500	Taq Master Mix	500 reactions	\$231.68
PR1001-R-200	Taq Master Mix Red Dye	200 reactions	\$92.77
PR1001-R-500	Taq Master Mix Red Dye	500 reactions	\$233.72

PR1MA™ High Fidelity Taq



Ideal for cloning, mutagenesis and microarrays

ITEM #	DESCRIPTION	QTY	PRICE
PR1001-HF-500	HF Taq	500 units	\$409.31
PR1001-HF-200	HF Master Mix	200 reactions	\$184.28
PR1001-HFHS-200	HFHS Master Mix	200 reactions	\$181.75

Ideal for multiplex PCR and amplification of GC-rich DNA

PR1MA™ Hot Start Taq



ITEM #	DESCRIPTION	QTY	PRICE
PR1000-HS-500	HS Taq	500 units	\$226.30
PR1000-HS-2500	HS Taq	2,500 units	\$1,141.44
PR1000-HS-10000	HS Taq	10,000 units	\$3,723.26
PR1001-HS-200	HS Taq Master Mix	200 reactions	\$169.36

PR1MA™ dNTPs



ITEM #	DESCRIPTION	QTY	PRICE
PR3040-M-2	40 mM dNTP Mix	10 mM/5 µmol each, 1 mL	\$108.71
PR3040-M-8	40 mM dNTP Mix	10 mM/20 µmol each, 4 mL	\$380.48
PR3101-S-1	100 mM dNTP set dATP, dCTP, dGTP, dTTP	250 µL/25 µmol each	\$148.34
PR3101-S-4	100 mM dNTP set dATP, dCTP, dGTP, dTTP	1 mL/100 µmol each	\$531.71
PR3101-S-20	100 mM dNTP set dATP, dCTP, dGTP, dTTP	20 x 250 µL/25 µmol each	\$3,587.37



24 ORDER TODAY!

midsci.com • custserv@midsci.com • 800.227.9997 • f: 636.225.9998

PCR/qPCR REAGENTS

- Superior sensitivity and fast cycling with exceptional results
- Ideal for low copy number templates
- Early Ct values and detection across a broad dynamic range
- Ready to use 2x mastermix
- Includes PR1MA™ Hot Start Taq Polymerase for greater specificity and accuracy
- Inert, blue tracking dye assists with accurate pipetting
- Optimized detection system specific, no change in protocol or results but significantly more cost effective

Tell us what qPCR system you have, and we can guide you to which enzyme is made for YOUR machine!



PR1MA™ qMAX Green or Probe Mix

ITEM #	DESCRIPTION	QTY	PRICE
PR2000-H-1000	Green High ROX	1,000 reactions	\$655.36
PR2000-L-1000	Green Low ROX	1,000 reactions	\$655.36
PR2000-N-1000	Green No ROX	1,000 reactions	\$655.36
PR2001-H-1000	Probe High ROX	1,000 reactions	\$655.36
PR2001-L-1000	Probe Low ROX	1,000 reactions	\$655.36
PR2001-N-1000	Probe No ROX	1,000 reactions	\$655.36



Need 500 or 5000 reactions?
Please visit midsci.com for other quantities

cDNA Synthesis Kit is a 2-tube format for easy reaction set up and is ideal for 4 pg to 0.5 µg of input RNA. One tube includes our exceptionally stable Reverse Transcriptase combined with a potent RNase inhibitor, and the other tube contains a 5X reaction buffer with an optimal mixture of anchored oligo (dT) primers and random hexamers to produce a non-biased population of cDNA.

1st Strand cDNA Synthesis Flex Kit includes four components: a high capacity Reverse Transcriptase, an optimized 5X buffer, and separate solutions of oligo (dT) primers, and random hexamer primers. This multi-component format is ideal for 10 pg to 2.0 µg of input RNA and allows for greater flexibility in assay design.



PR1MA™ qMAX cDNA Synthesis Kit



ITEM #	DESCRIPTION	QTY	PRICE
PR2100-C-S	cDNA Synthesis Sample	5 reactions	\$19.75
PR2100-C-25	cDNA Synthesis Kit	25 reactions	\$128.59
PR2100-C-100	cDNA Synthesis Kit	100 reactions	\$452.60
PR2100-C-250	cDNA Synthesis Kit	250 reactions	\$1,073.29
PR2110-S	First Strand cDNA Synthesis Flex Sample Kit	5 reactions	\$12.40
PR2110-50	First Strand cDNA Synthesis Flex Kit	50 reactions	\$110.10


Compare to *IScript™*, *Superscript™ II/III/IV*, *VILO*, *Thermo™ First Strand* - let us know what kit you use and we can give you a money saving swap!

PCR/qPCR REAGENTS

- Ready-to-use formulation includes loading buffer and tracking dye
- 500 µL suitable for 100 lanes (5 µL per lane)

- Higher intensity reference bands
- Ultra pure production allows economical ambient shipping
- 0.1 mg/mL concentration

Additional loading buffer is not required!



PR1MA™ DNA Ladders

ITEM #	DESCRIPTION	QTY	PRICE
PR4010-100	100bp DNA Ladder	500 µL/100 lanes	\$42.34
PR4010-200	100bp DNA Ladder	2 x 500 µL/200 lanes	\$108.71
PR4010-500	100bp DNA Ladder	5 x 500 µL/500 lanes	\$189.53
PR4100-200	1kb DNA Ladder	2 x 500 µL/200 lanes	\$108.71
PR4100-600	1kb DNA Ladder	6 x 500 µL/600 lanes	\$266.34
PR4005-100	50bp DNA Ladder	500 µL/100 lanes	\$42.34
PR4005-500	50bp DNA Ladder	5 x 500 µL/500 lanes	\$189.53

- Replaces hazardous Ethidium Bromide (EtBr)
- Better sensitivity than EtBr
- Detect as little as 0.1ng of DNA
- Excitation by UV or Blue Light

PR1MA™ SmartGlow™ DNA Safe Stain

ITEM #	DESCRIPTION	QTY	PRICE
BE4500-LD	Loading Dye, Blue & UV light	1 mL	\$126.20
BE4500-PS	Pre-Stain, Blue & UV light	1 mL	\$126.20



- Yields high resolution sharp DNA bands with high clarity and low background
- Optimized gel strength for ease of gel preparation
- High purity and multi-purpose

PR1MA™ Agarose

ITEM #	DESCRIPTION	QTY	PRICE
MIDSCI-500	General Purpose	500 grams	\$240.07

CHEMICALS

Chemicals are packed
at ISO/GMP Facility



Find more chemicals and
information online!

PRIMA™ General Chemicals

ITEM #	DESCRIPTION	QTY	PRICE
KCA20080-100	100 g, Agarose, Low EEO	1	\$195.01
KCA20600-1000	1 kg, Ammonium Acetate	1	\$60.46
KCB32050-1000	1 kg, Boric Acid	1	\$51.90
KCA22645-4000	4 L, Ethanol, Reagent Grade	1	\$85.60
KCP50220-0.1	100 mg, Proteinase K	1	\$140.97

PRIMA™ Antibiotics

ITEM #	DESCRIPTION	QTY	PRICE
KCA40040-10	10 g, Ampicillin, Sodium Salt	1	\$73.83
KCB12200-0.025	25 mg, Blasticidin S HCl - Powder	1	\$101.38
KCC46000-5	5 g, Carbenicillin Disodium Salt	1	\$171.47
KCC51000-5	5 g, Cefotaxime	1	\$279.00
KCP40160-5	5 g, Polymixin B Sulfate	1	\$255.73

PRIMA™ Buffers

ITEM #	DESCRIPTION	QTY	PRICE
KCC34500-1000	1 kg, Citric Acid Anhydrous	1	\$60.46
KCG36050-1000	1 kg, Glycine, Free Base	1	\$73.30
KCH75030-1000	1 kg, HEPES, Free Acid	1	\$396.97
KCT60040-1000	1 kg, Tris Base Ultra Pure	1	\$94.70
KCT60050-1000	1 kg, TRIS	1	\$132.41

PRIMA™ Growth Media

ITEM #	DESCRIPTION	QTY	PRICE
KCA20020-500	500 g, Agar, High Gel Strength, Powder	1	\$196.35
KCG22020-1	1 L, Glycerol	1	\$51.90
KCL24066-1000	1 kg, LB Broth, Granulated [Lennox LB Broth]	1	\$148.73
KCL24022-1000	1 kg, Luria Agar Granulated [Miller's LB Agar]	1	\$232.73
KCP20250-1000	1 kg, Peptone, Granulated	1	\$148.73
KCT60060-1000	1 kg, Tryptone, Powder	1	\$140.97

COLD STORAGE SUPPLIES

- Manufactured from high-clarity, ultra-pure polypropylene (meets FDA 21 CFR 177.1520 & USP Class VI)
- Tubes comply with the In Vitro Devices Directive (IVD)

- Contrasting black 100 µL graduations are bold and easy to read
- Large white writing area on all tubes
- Optimized for vapor-phase liquid nitrogen storage
- Available with round base or self-standing star base
- Resealable packaging, tubes are packaged in resealable stand-up pouches that minimize contamination and provide easy access

PR1MA™ Cryogenic Tubes

ITEM #	DESCRIPTION	QTY	PRICE
PR-CVET10-SB	1.0 mL, Cryo Tubes, Ext. Thread, Star Base, Sterile, Clear	100/pk	\$26.97
PR-CVET18-RB	1.8 mL, Cryo Tubes, Ext. Thread, Round Base, Sterile, Clear	100/pk	\$27.56
PR-CVET18-SB	1.8 mL, Cryo Tubes, Ext. Thread, Star Base, Sterile, Clear	100/pk	\$27.56
PR-CVET45-SB	4.5 mL, Cryo Tubes, Ext. Thread, Star Base, Sterile, Clear	100/pk	\$35.66
PR-CVIT10-SB	1.0 mL, Cryo Tubes, Int. Thread, Star Base, Sterile, Clear	100/pk	\$30.66
PR-CVIT18-RB	1.8 mL, Cryo Tubes, Int. Thread, Round Base, Sterile, Clear	100/pk	\$30.66
PR-CVIT18-SB	1.8 mL, Cryo Tubes, Int. Thread, Star Base, Sterile, Clear	100/pk	\$30.66
PR-CVIT36-RB	3.6 mL, Cryo Tubes, Int. Thread, Round Base, Sterile, Clear	100/pk	\$40.22
PR-CVIT36-SB	3.6 mL, Cryo Tubes, Int. Thread, Star Base, Sterile, Clear	100/pk	\$40.22
PR-CVIT45-RB	4.5 mL, Cryo Tubes, Int. Thread, Round Base, Sterile, Clear	100/pk	\$41.13
PR-CVIT45-SB	4.5 mL, Cryo Tubes, Int. Thread, Star Base, Sterile, Clear	100/pk	\$41.13
PR-CVCI	Cap Inserts Cryo Tubes, Assorted	500/pk	\$31.65

Note: not suitable for immersion in liquid nitrogen




- Unique cap & vial design eliminates silicone gasket pops
- Self standing

- Externally threaded for reduced contamination
- Color coded cap inserts available

TPP® CryoTubes

ITEM #	DESCRIPTION	QTY	PRICE
TP89012-100	1.2 mL [PP], External Thread, Self Standing, Natural, Sterile	100/pk	\$85.79
TP89020-100	2 mL [PP], External Thread, Self Standing, Natural, Sterile	100/pk	\$88.16
TP89050-100	4.5 mL [PP], External Thread, Self Standing, Natural, Sterile	100/pk	\$109.35
TP99020	Cap Inserts, Six Assorted Colors	600/cs	\$95.23



COLD STORAGE SUPPLIES

- Used for upright freezer or chest freezer racks
- Temperature range: -196°C to +121°C
- Drain holes on bottom and writing area on lid
- Standard 5.25 x 5.25 x 2 inch
- Available cell dividers: 144, 100, 81, 64, 49, 36, 25, or 16



PR1MA™ Cardboard Freezer Boxes

ITEM #	DESCRIPTION	QTY	PRICE
SB2CC-81	Standard White Cardboard 5.25 x 5.25 x 2 in. Box w/ 81 Dividers	1/pk	\$5.00
SB2CC-100	Standard White Cardboard 5.25 x 5.25 x 2 in. Box w/ 100 Dividers	1/pk	\$5.00

Need other options?
Just let your trusted MIDSCI™ Lab Consultant know what you need.

- Stores tubes up to 2.0 mL
- Suitable for temperatures as low as -90°C
- Rugged polypropylene
- High-contrast alphanumeric indexing on lid

81-Well Freezer Storage Boxes with Lift-Off Lids

ITEM #	DESCRIPTION	QTY	PRICE
FB-81	81-Well 5.25 x 5.25 x 2 in. [PP], w/ Lift-Off Lid, Clear	5/pk	\$37.46
FB-81C	81-Well 5.25 x 5.25 x 2 in. [PP], w/ Lift-Off Lid, Assorted	5/pk	\$40.65



- Fit all the most popular cryo and cold storage vials
- Fit into standard freezer racks
- Temperature range: -196°C to 121°C
- Drain holes for venting and safety



Polycarbonate Freezer Box

ITEM #	DESCRIPTION	QTY	PRICE
MPB2-25	25-Well, 3 x 3 x 2.06 in. [PC], White	1/ea	\$8.35
PB2-81	81-Well, 5.25 x 5.25 x 2 in. [PC], White	1/ea	\$10.99
PB4-81	81-Well, 5.25 x 5.25 x 3.25 in. [PC], White	1/ea	\$13.29
PB2-100	100-Well, 5.25 x 5.25 x 2 in. [PC], White	1/ea	\$10.99

COLD STORAGE SUPPLIES

Upright Freezer Racks

ITEM #	RACK DIMS D x H x W	BOX HEIGHT	BOX CONFIG (DxH)	TYPE	PRICE
UF-332	16.88 x 6.1 x 5.5 in.	2 in.	3 x 3	Universal	\$75.97
UF-342	16.88 x 8.88 x 5.5 in.	2 in.	3 x 4	Universal	\$97.52
UF-432	22.19 x 6.1 x 5.5 in.	2 in.	4 x 3	Universal	\$92.24
UF-442	26.75 x 8.88 x 5.5 in.	2 in.	4 x 4	Universal	\$108.97
UF-452	22.19 x 5.5 x 11 in.	2 in.	4 x 5	Universal	\$141.09
UF-552	27.31 x 11 x 5.5 in.	2 in.	5 x 5	Universal	\$289.00
UFE-332	16.81 x 6.63 x 5.5 in.	2 in.	3 x 3	Eco	\$68.39
UFE-342	16.81 x 8.88 x 5.5 in.	2 in.	3 x 4	Eco	\$78.42
UFE-432	22.13 x 6.63 x 5.5 in.	2 in.	4 x 3	Eco	\$83.02
UFE-442	22.13 x 8.88 x 5.5 in.	2 in.	4 x 4	Eco	\$122.88
UFE-452	22.13 x 11 x 5.5 in.	2 in.	4 x 5	Eco	\$139.60
UFE-552	27.19 x 11 x 5.5 in.	2 in.	5 x 5	Eco	\$175.07



UF-332



UFE-332



UFD-332



UFD-542

Upright Drawer Freezer Racks

ITEM #	D x H x W	BOX HEIGHT	BOX CONFIG (DxH)	PRICE
UFD-332	16.75 x 7.13 x 5.5 in.	2 in.	3 x 3	\$100.04
UFD-342	16.75 x 9.44 x 5.5 in.	2 in.	3 x 4	\$134.38
UFD-432	22.31 x 7.13 x 5.5 in.	2 in.	4 x 3	\$134.38
UFD-442	22.31 x 9.44 x 5.5 in.	2 in.	4 x 4	\$159.70
UFD-452	22.44 x 11.81 x 5.5 in.	2 in.	4 x 5	\$189.23
UFD-552	27 x 11.81 x 5.5 in.	2 in.	5 x 5	\$283.10

Vertical Freezer Racks

ITEM #	D x H x W	BOX HEIGHT	BOX CONFIG	PRICE
CF-9-2	20.06 x 5.63 x 5.5 in.	2 in.	9	\$69.18
CF-10-2	22.31 x 5.63 x 5.5 in.	2 in.	10	\$76.98
CF-11-2	24.5 x 5.63 x 5.5 in.	2 in.	11	\$84.49
CF-12-2	26.63 x 5.63 x 5.5 in.	2 in.	12	\$92.24



CF-9-2



CF-12-2

Not sure what rack will maximize storage in your freezer? - No Problem.

Visit our Free Freezer Rack Configurator for more information



RESERVOIRS


Single Chamber Reagent Reservoirs

ITEM #	DESCRIPTION	QTY	PRICE
P7005	5 mL, Non-Sterile	300/cs	\$114.40
P7005-1S	5 mL, Sterile, Individ. Wrapped	100/cs	\$114.40
P7005-5S	5 mL, Sterile	5/bg, 20 bgs/cs	\$114.40
P8010	10 mL, Non-Sterile	300/cs	\$114.40
P8010-1S	10 mL, Sterile, Individ. Wrapped	100/cs	\$114.40
P8010-5S	10 mL, Sterile	5/bg, 20 bgs/cs	\$114.40
P8025	25 mL, Non-Sterile	100/cs	\$59.87
P8025-1S	25 mL, Sterile, Individ. Wrapped	100/cs	\$125.50
P8025-5S	25 mL, Sterile	5/bg, 40 bgs/cs	\$125.50
P8050	50 mL, Non-Sterile	100/cs	\$59.87
P8050-1S	50 mL, Sterile, Individ. Wrapped	100/cs	\$125.50
P8050-5S	50 mL, Sterile	5/bg, 40 bgs/cs	\$125.50
P8100	100 mL, Non-Sterile	100/cs	\$59.87
P8100-1S	100 mL, Sterile, Individ. Wrapped	100/cs	\$125.50
P8100-5S	100 mL, Sterile	5/bg, 40 bgs/cs	\$125.50

- Stable base limits tipping and flipping
- Pour corner spouts, to save on reagents
- Ideal for multi-channel pipetting



Same features as above + 2 chambers!



2 Chamber Reservoirs

ITEM #	DESCRIPTION	QTY	PRICE
P8051	2 Chamber, Total Vol 50 mL, Non-Sterile	100/cs	\$59.87
P8051-1S	2 Chamber, Total Vol 50 mL, Sterile, Individ. Wrapped	100/cs	\$125.50
P8051-5S	2 Chamber, Total Vol 50 mL, Sterile	5/bg, 40 bgs/cs	\$125.50

Green Reagent Reservoirs

ITEM #	DESCRIPTION	QTY	PRICE
RB25	25 mL Reservoir Base, Non-Sterile	2/pk	\$14.00
RR25	25 mL Green Reagent Reservoirs, Sterile	200/pk	\$172.00
RB100	100 mL Reservoir Base, Non-Sterile	10/pk	\$72.00
RR100	100 mL Green Reagent Reservoirs, Sterile	200/pk	\$172.00

Innovative Design Saves Space and Reduces Plastic Waste!



PETRI DISHES & SPREADERS

- Slideable design for automated filling and streaking systems
- Features side arrows and ISO mark target for automation
- Stackable design for more stability
- High quality, value priced product



PR1MA™ Petri Dishes

ITEM #	DESCRIPTION	QTY	PRICE
667665	60 x 15 mm, Sterile	25/bg, 20 bgs/cs	\$133.19
667694	100 x 15 mm, Sterile, Slideable	25/bg, 20 bgs/cs	\$120.22
667695	100 x 15 mm, Sterile, Stackable	25/bg, 20 bgs/cs	\$120.22
667682	100 x 15 mm, Sterile, I-Plate (2-Section)	25/bg, 20 bgs/cs	\$170.77
667683	100 x 15 mm, Sterile, Y-Plate (3-Section)	25/bg, 20 bgs/cs	\$173.41
667684	100 x 15 mm, Sterile, X-Plate (4-Section)	25/bg, 20 bgs/cs	\$176.26




Standard Petri Dishes



ITEM #	DESCRIPTION	QTY	PRICE
900	100 x 15 mm, Sterile, Stackable	25/pk, 500/cs	\$162.00
901	60 x 15 mm, Sterile, Stackable	20/pk, 500/cs	\$135.00
902	150 x 15 mm, Sterile, Slippable	10/pk, 100/cs	\$130.00
906	100 x 25 mm, Sterile, Extra Deep Dish, Stackable	17/pk, 340/cs	\$164.00
910	100 x 15 mm, Sterile, Slippable	25/pk, 500/cs	\$162.00
911	100 x 15 mm, Sterile, I-Plate (2-Section), Slippable	25/pk, 500/cs	\$190.00
912	100 x 15 mm, Sterile, Y-Plate (3 Section), Stackable	25/pk, 500/cs	\$190.00

Disposable Spreaders

ITEM #	DESCRIPTION	QTY	PRICE
LS-SIP	L-Shape, 38 x 156 mm, Sterile, Individ. Wrapped	1/pk, 500/cs	\$114.67
LS-SP	L-Shape, 38 x 156 mm, Sterile	10/pk, 500/cs	\$98.12
TS-SIP	T-Shape, 34 x 140 mm, Sterile, Individ. Wrapped	1/pk, 500/cs	\$117.04
TS-SP	T-Shape, 34 x 140 mm, Sterile	10/pk, 500/cs	\$108.82
TR-SB-S	Triangle Shape, 30 x 207 mm, Sterile	10/bg, 100/cs	\$85.12
TR-SB-L	Triangle Shape, 60 x 235 mm, Sterile	10/bg, 100/cs	\$85.12



ERLENMEYER FLASKS

Erlenmeyer Shaker Flasks

ITEM #	DESCRIPTION	QTY	PRICE
F4060-F	125 mL [PETG], Flat Base, Sterile, Disposable	24/pk	\$219.58
F4061-F	250 mL [PETG], Flat Base, Sterile, Disposable	12/pk	\$123.36
F4062-F	500 mL [PETG], Flat Base, Sterile, Disposable	12/pk	\$159.13
F4063-F	1,000 mL [PETG], Flat Base, Sterile, Disposable	6/pk	\$120.28
F4060-B	125 mL [PC], Baffled Base, Sterile, Autoclavable/ Reusable	24/pk	\$245.69
F4061-B	250 mL [PC], Baffled Base, Sterile, Autoclavable/ Reusable	12/pk	\$139.81
F4062-B	500 mL [PC], Baffled Base, Sterile, Autoclavable/ Reusable	12/pk	\$184.01
F4063-B	1,000 mL [PC], Baffled Base, Sterile, Autoclavable/ Reusable	6/pk	\$127.88
F4086	2.8L Jumbo Erlen/Fernbach, Styrene/Butadiene, Baffled Base, Sterile	6/pk	\$290.00
F4087	5.0L Jumbo Erlen/Fernbach, Styrene/Butadiene, Baffled Base, Sterile	4/pk	\$390.64
F4067	Caps w/ 0.2 µm Filter, for 125 & 250 mL Flasks	10/pk	\$29.81
F4068	Caps w/ 0.2 µm Filter, for 500 & 1,000 mL Flasks	10/pk	\$34.95
F4069	Caps w/ 0.2 µm Filter, for 2.8 L & 5.0 L Jumbo Flasks	10/pk	\$40.09



PR1MA™ Erlenmeyer Flasks

ITEM #	DESCRIPTION	QTY	PRICE
667801	125 mL [PETG], Clear, Vent Cap, Plain Bottom, Sterile	1/pk, 12 pks/cs	\$116.00
667805	250 mL [PETG], Clear, Vent Cap, Plain Bottom, Sterile	1/pk, 12 pks/cs	\$131.00
667809	500 mL [PETG], Clear, Vent Cap, Plain Bottom, Sterile	1/pk, 12 pks/cs	\$167.00
667813	1,000 mL [PETG], Clear, Vent Cap, Plain Bottom, Sterile	1/pk, 6 pks/cs	\$127.00

High Efficiency Erlenmeyer Flasks

ITEM #	DESCRIPTION	QTY	PRICE
MSI-785101	2 L [PC], HDPE Vent Cap, Indiv. Wrapped, Sterile	4/pk	\$163.39
MSI-786111	3 L [PC], HDPE Vent Cap, Indiv. Wrapped, Sterile	5/pk	\$282.75
MSI-786511	3 L [PC], Wide Mouth, HDPE Vent Cap, Indiv. Wrapped, Sterile	6/pk	\$282.75
MSI-787011	5 L, HDPE Vent Cap, Indiv. Wrapped, Sterile	7/pk	\$422.85



FILTRATION

- Polypropylene housing with nylon mesh filter
- Extended gripping tab on housing for easy handling

- The elongated lip on the housing extends beyond the edge of the centrifuge tube enabling aseptic handling with forceps



PR1MA™ Cell Strainers

ITEM #	DESCRIPTION	QTY	PRICE
40ICS	40 µm, Blue, Individ. Wrapped, Sterile	1/pk, 50 pks/cs	\$127.69
70ICS	70 µm, White, Individ. Wrapped, Sterile	1/pk, 50 pks/cs	\$127.69
100ICS	100 µm, Yellow, Individ. Wrapped, Sterile	1/pk, 50 pks/cs	\$127.69
SG-PEST	5.5 in. Pestle, Blue, Individ. Wrapped, Sterile	1/pk, 50 pks/cs	\$61.54

Same great products,
Just a new name!

- Advanced design - **faster flow through**
- Superior ventilation - **prevents clogging/spills**
- Designed for use with 50 mL tubes (adapter for 5 & 15 mL tubes included)



PR1MA™ Signature Cell Strainers

ITEM #	DESCRIPTION	QTY	PRICE
40CELL	40 µm, Blue, w/ Reducing adapter for 15 mL Tube, Sterile	50/pk	\$112.05
70CELL	70 µm, White, w/ Reducing adapter for 15 mL Tube, Sterile	50/pk	\$112.05
100CELL	100 µm, Yellow, w/ Reducing adapter for 15 mL Tube, Sterile	50/pk	\$112.05
CELLADAPT	Reducing Adapters for 5 mL and 15 mL Tubes, Sterile, Individ. Wrapped	25/pk	\$49.22



Same great product,
Just a new name!

PR1MA™ Signature Double Sided Pestle

ITEM #	DESCRIPTION	QTY	PRICE
PEST2	Double Sided Pestle for Cell Strainers and Tubes, Individ. Wrapped, Sterile	50/pk	\$52.43

TISSUE CELL CULTURE

- Unique design allows for more of what you want through and less of what you don't
- 0.22 µm PES membrane is the lowest protein membrane available
- 45 mm neck thread on the filter - Filter tops can fit bottles with 45 mm neck thread
- Gamma-sterilized, DNase/RNase-Free, and pyrogen-Free
- Individually wrapped in a peel-off plastic wrapper
- Screw cap & vacuum accessories included in each sterile package



TPP® Filtermax Filters

ITEM #	DESCRIPTION	QTY	PRICE
TP99150	150 mL, Filter Top & Reservoir, Sterile	1/pk, 18 pks/cs	\$406.04
TP99157	150 mL, Reservoir, Sterile	1/pk, 24 pks/cs	\$218.38
TP99250	250 mL, Filter Top & Reservoir, Sterile	1/pk, 12 pks/cs	\$307.13
TP99257	250 mL, Reservoir, Sterile	1/pk, 24 pks/cs	\$267.11
TP99500	500 mL, Filter Top & Reservoir, Sterile	1/pk, 10 pks/cs	\$323.13
TP99505	500 mL, Filter Top, Sterile	1/pk, 21 pks/cs	\$509.63
TP99507	500 mL, Reservoir, Sterile	1/pk, 36 pks/cs	\$677.58
TP99950	1 L, Filter Top & Reservoir, Sterile	1/pk, 9 pks/cs	\$568.85
TP99955	1 L, Filter Top, Sterile	1/pk, 12 pks/cs	\$539.97
TP99957	1 L, Reservoir, Sterile	1/pk, 20 pks/cs	\$575.75

- Design allows for faster filtration & less clogging
- Fits the rim of a 50 mL perfectly

- 0.22 µm or 0.45 µm PES membrane is the lowest protein membrane available



TPP® Syringe Filters

ITEM #	DESCRIPTION	QTY	PRICE
TP99722	0.22 µm, Sterile	1/pk, 40 pk/cs	\$213.81
TP99745	0.45 µm, Sterile	1/pk, 40 pk/cs	\$213.81

TISSUE CELL CULTURE



- Ideal for starting cultures & primary cells
- Culture & spin down - same tube
- Reduces plastics, cost, & lose fewer cells during transfer

TPP® Tiny Flasks

ITEM #	DESCRIPTION	QTY	PRICE
TP91106	20 cm ² (16 x 121 mm) Mini Roller Bottle [PS], Vent Cap [PE], Sterile	20/bg, 40 bgs/cs	\$1,073.44
TP91243	10 cm ² (30 x 105 mm) Tube Flask [PP], Filter Cap [PE] 0.22 µm PTFE Membrane, Sterile	4/pk, 54 pks/cs	\$523.14

- Culture & image on a (COP) slide
- Suitable for: Confocal fluorescence microscopy, individual coatings, staining, and immunofluorescence measurements
- Resistant to solvents (e.g. acetone, ethanol, xylene)
- Removable medium chamber with a "click" and without the need of a removing tool

TPP® Clipmax Chamber Slide System

ITEM #	DESCRIPTION	QTY	PRICE
TP70010	10 cm ² (93 x 33 x 32 mm) Tube Flask [PS], Filter Cap [PE] 0.22 µm PTFE Membrane, Sterile	1/pk, 5 pks/cs	\$112.84



Need other chamber slide options? Just ask!



- Scraper heads made from a flexible polyethylene material with sharp edges
- Single-wrapped in blister packaging
- **Swivel head technology:**
Allows access to entire flask!

TPP® Scraper & Spatula

ITEM #	DESCRIPTION	QTY	PRICE
TP99002	24 cm Scraper "S" [PE], Rotating, Sterile	1/pk, 150 pks/cs	\$649.35
TP99003	30 cm Scraper "M" [PE], Rotating, Sterile	1/pk, 150 pks/cs	\$734.76
TP99004	38 cm Scraper "L" [PE], Rotating, Sterile	1/pk, 100 pks/cs	\$606.98
TP99010	20 cm Spatula [PP], Fixed, Sterile	1/pk, 100 pks/cs	\$573.26

TISSUE CELL CULTURE

- Ideal for suspension cultures
- Scale up or scale down ability to do large in parallel studies
- Better viability & proliferation than flasks or bottles
- Better gas exchange control & cyclonic stirring than flasks or bottles



The standard in bioreacting labs across the world!

TPP® TubeSpin® Bioreactor

ITEM #	DESCRIPTION	QTY	PRICE
TP87015	15 mL, Conical, [PP], Sterile	40/bg, 20 bgs/cs	\$1,154.58
TP87017	15 mL, Round, [PP], Sterile	40/bg, 20 bgs/cs	\$1,730.05
TP87050	50 mL, Conical, [PP], Sterile	20/bg, 9 bgs/cs	\$503.07
TP87056	50 mL, Round, [PP], Sterile	20/bg, 9 bgs/cs	\$784.47
TP87450	400 mL, Conical, [PP], Sterile	1/bg, 32 bgs/cs	\$668.70
TP87600	600 mL, Conical, [PP], Sterile	1/bg, 26 bgs/cs	\$1,055.00



Need different quantities or samples?
Contact us for more options!

- Ideal for high throughput & robotics
- Re-closeable septum allows repeated access while keeping sterile
- Compatible for injection needles up to 2.1 mm diameter
- Ideal for PRP kits



TPP® TubeSpin® Bioreactor with Septum

ITEM #	DESCRIPTION	QTY	PRICE
TP86050	50 mL, Conical, with septum, [PP], Sterile	20/bg, 9 bgs/cs	\$2,892.55
TP186550	50 mL, Conical, with septum, no barcode, Individ. Wrapped, [PP], Sterile	1/bg, 220 bgs/cs	\$1,758.94



Barcoded for automated tracking

TISSUE CELL CULTURE

WHY TPP?

TPP® is Swiss made with excellent flatness and clarity

- **Flatness:** Reduces clumping, ensuring even growth
- **Vents:** Promotes rapid temperature equilibration for better cell health
- **Clarity:** Ideal for clear microscopy
- **Stacked Design:** Dishes, flasks, and plates interlock with vents between levels
- **Stability:** Interlocking feature prevents drops when handling stacks



- **Minimized edge effect in 96-well plates**
- **Easy-to-read yellow marking spot and black alphanumeric**
- **Made from optically clear polystyrene for excellent viewing**
- **Sterile, packaged individually or in convenient multi-packs**

Comfortably stack your plates!

TPP® Tissue Culture Plates

ITEM #	DESCRIPTION	QTY	PRICE
TP92006	6-Well TC Plate [PS], Flat Bottom, Sterile	1/pk, 126/cs	\$738.43
TP92012	12-Well TC Plate [PS], Flat Bottom, Sterile	1/pk, 126/cs	\$739.11
TP92024	24-Well TC Plate [PS], Flat Bottom, Sterile	1/pk, 126/cs	\$737.79
TP92048	48-Well TC Plate [PS], Flat Bottom, Sterile	1/pk, 126/cs	\$738.44
TP92096	96-Well TC Plate [PS], Flat Bottom, Sterile	1/pk, 162/cs	\$911.24
TP92097	96-Well TC Plate [PS], Round Bottom, Sterile	1/pk, 162/cs	\$906.44

Also available in multi-packs!
Contact your trusted MIDSCI™ Lab Consultant for more information

- **Easy-to-read painted yellow marking spot**
- **Ridged grip around base greatly improves handling. No more dropped dishes!**

TPP® Tissue Culture Dishes

ITEM #	DESCRIPTION	QTY	PRICE
TP93040	40 x 11 mm TC Dish [PS], Sterile	20/pk, 45 pk/cs	\$730.89
TP93060	60 x 16 mm TC Dish [PS], Sterile	14/pk, 60 pk/cs	\$752.35
TP93100	96 x 21 mm TC Dish [PS], Sterile	10/pk, 24 pk/cs	\$408.61
TP93150	146 x 21 mm TC Dish [PS], Sterile	5/pk, 20 pk/cs	\$485.12



TISSUE CELL CULTURE

TPP® Tissue Culture Flasks

- Wider neck allows for better access for 100% cell retrieval when scraping
- Short access to cell floor reduces chance of contamination



ITEM #	DESCRIPTION	QTY	PRICE
TP90025	25 cm ² (92 x 51 x 29 mm) TC Flask [PS], Solid Vent Cap [PE], Sterile	10/pk, 36 pks/cs	\$824.60
TP90026	25 cm ² (92 x 51 x 29 mm) TC Flask [PS], Filter Cap [PE], Sterile	10/pk, 36 pks/cs	\$900.79
TP90075	75 cm ² (155 x 87 x 40 mm) TC Flask [PS], Solid Vent Cap [PE], Sterile	5/pk, 20 pks/cs	\$413.95
TP90076	75 cm ² (155 x 87 x 40 mm) TC Flask [PS], Filter Cap [PE], Sterile	5/pk, 20 pks/cs	\$485.57
TP90150	150 cm ² (210 x 122 x 50 mm) TC Flask [PS], Solid Vent Cap [PE], Sterile	3/pk, 13 pks/cs	\$379.05
TP90151	150 cm ² (210 x 122 x 50 mm) TC Flask [PS], Filter Cap [PE], Sterile	3/pk, 13 pks/cs	\$403.28
TP90300	300 cm ² (275 x 170 x 50 mm) TC Flask [PS], Solid Vent Cap [PE], Sterile	3/pk, 6 pks/cs	\$404.46
TP90301	300 cm ² (275 x 170 x 50 mm) TC Flask [PS], Filter Cap [PE], Sterile	3/pk, 6 pks/cs	\$423.96

Filter caps come with a 0.22 µm PTFE Membrane

- Access of a dish & security of a flask

TPP® Re-Closable Flasks

- Ideal for long incubations, colony picking, organ growth, & so much more

ITEM #	DESCRIPTION	QTY	PRICE
TP90552	150 cm ² (210 x 122 x 50 mm) TC Flask [PS], Filter Cap [PE], Sterile	3/pk, 6 pks/cs	\$556.71
TP90652	115 cm ² (210 x 122 x 50 mm) TC Flask [PS] w/ Barrier, Filter Cap [PE], Sterile	3/pk, 6 pks/cs	\$558.69

Filter caps come with a 0.22 µm PTFE Membrane



TPP® Peel-Off Flasks



ITEM #	DESCRIPTION	QTY	PRICE
TP90028	25 cm ² (92 x 51 x 31 mm) TC Flask [PS], Filter Cap [PE], Sterile	5/pk, 14 pks/cs	\$975.94
TP90078	75 cm ² (155 x 87 x 45 mm) TC Flask [PS], Filter Cap [PE], Sterile	2/pk, 16 pks/cs	\$508.51
TP90653	115 cm ² (210 x 122 x 50 mm) TC Flask [PS] w/ Barrier, Filter Cap [PE], Sterile	3/pk, 6 pks/cs	\$470.22
TP90153	150 cm ² (210 x 122 x 50 mm) TC Flask [PS], Filter Cap [PE], Sterile	3/pk, 6 pks/cs	\$462.33
TP90303	300 cm ² (275 x 170 x 55 mm) TC Flask [PS], Filter Cap [PE], Sterile	2/pk, 4 pks/cs	\$208.18

MIXING



- Powerful vortexing for tubes up to 50 mL
- Pressure activated operation
- "Mini sized", less than 4 inches wide
- 2-year warranty

PR1MA™ Mini Vortexer

ITEM #	DESCRIPTION	QTY	PRICE
VM1-B	Mini Vortexer, Blue, 4 mm Orbit, Max Speed 2,800 rpm	1/ea	\$229.09
VM1-G	Mini Vortexer, Green, 4 mm Orbit, Max Speed 2,800 rpm	1/ea	\$229.09
VM1-GR	Mini Vortexer, Gray, 4 mm Orbit, Max Speed 2,800 rpm	1/ea	\$229.09
VM1-P	Mini Vortexer, Purple, 4 mm Orbit, Max Speed 2,800 rpm	1/ea	\$229.09
VM1-R	Mini Vortexer, Red, 4 mm Orbit, Max Speed 2,800 rpm	1/ea	\$229.09

- Powerful motor for INSTANT vortexing
- Continuous or "touch" operation
- Q-Drive™ dynamic balancing system
- Wide variety of accessories
- 2-year warranty

PR1MA™ Vortex Mixers

ITEM #	DESCRIPTION	QTY	PRICE
VM-3200	Adjustable Vortex Mixer, 3 mm Orbit, Max Speed 3,200 rpm	1/ea	\$296.62
VM-COMBO	Combo Head Attachment, 1.5 mL, 0.5 mL Microtubes & 1 Microplate	1/ea	\$80.27
VM-FLAT	Flat Head Attachment, 3 in.	1/ea	\$44.76
VM-H15	Horizontal Head, 12 x 1.5 mL Microtubes	1/ea	\$113.99
VM-H150	Horizontal Head, 4 x 15 mL Microtubes	1/ea	\$113.99
VM-H500	Horizontal Head, 2 x 50 mL Microtubes	1/ea	\$113.99



- Reversible top platform, cup, or flat-dimpled surface
- Powerful motor for INSTANT vortexing
- Continuous or touch operation
- Q Drive dynamic balancing system
- Variable speed from 200 to 3,200 rpm
- Wide variety of accessories, for tubes up to 50 mL
- 2-year warranty



PR1MA™ V2 Vortex Mixer

ITEM #	DESCRIPTION	QTY	PRICE
VM-1003	Vortex Mixer, 200 to 3200 rpm	1/ea	\$335.38

MIXING

- Accepts a wide range of microplates, including deep well, PCR, Elisa plates, etc.
- Includes tube accessories for vortexing of 1.5 / 2.0 mL, 15 mL, or 50 mL tubes
- Converts to a traditional "touch" mode vortexer with included cup head attachment
- Powerful vortexing at speeds up to 3,000 rpm
- 2-year warranty

NEW!



Benchmark™ BenchMixer™ MP Vortexer

ITEM #	DESCRIPTION	QTY	PRICE
BV1007	Vortexer, 3,000 rpm (115 V)	1/ea	\$499.00

- Blocks available from 96 x 0.2 mL to 6 x 50 mL, custom, or multi-purpose
- Speed 200 to 1,500 rpm
- Heating rate 5°C/min - Cooling rate 2.5°C/min
- 2-year warranty

Multi-Therm™ Incubated Vortexer



ITEM #	DESCRIPTION	QTY	PRICE
MTH-1500	Temp Controlled Vortexer, Ambient +5°C to 100°C	1/ea	\$2,449.92
MTHC-1500	Temp Controlled Vortexer, Ambient -20°C to 100°C	1/ea	\$3,218.86
MTB-500	Block, 50 mL x 6 Tubes*	1/ea	\$499.14
MTB-5MT	Block, 5 mL x 15 conical Tubes*	1/ea	\$499.14
MTB-02	Block, 0.2 mL x 96 Tubes or 1 PCR Plate*	1/ea	\$421.01
MTB-05	Block, 0.5 mL x 54 Tubes*	1/ea	\$354.33
MTB-15	Block, 1.5 mL x 35 Tubes*	1/ea	\$354.33

*Blocks sold separately



Additional Blocks Available - visit www.midsci.com for more options!

- Adjustable speed & tilt for a wide range of applications
- Large 14 x 12 in. work surface
- Dual platform accessory doubles work surface
- Cold room and incubator safe
- 2-year warranty

PR1MA™ Variable Speed 2D Rocker



ITEM #	DESCRIPTION	QTY	PRICE
R2D-30	2D Rocker, 14 x 12 in. w/ Flat Mat	1/ea	\$892.81
BR2000-STACK	14 x 12 in. Stacking Platform w/ Flat Mat	1/ea	\$155.17

MIXING

- Gentle orbital motion for gel applications



- Horizontally circular 19 mm orbit for aeration & mixing
- Large 14 x 12 inch work surface

- Variable speed control from 3 to 70 rpm
- Incubator and cold room safe
- 2-year warranty

Orbi-Blotter™

ITEM #	DESCRIPTION	QTY	PRICE
BT30	14 x 12 in. Orbital Shaker, 3 x 2 L or 9 lbs Capacity	1/ea	\$1,014.59
BR2000-STACK	14 x 12 in. Stacking Platform	1/ea	\$155.17
BR2000-FLAT	14 x 12 in. Flat Mat	1/ea	\$53.58

- Optimized speed and tilt for gel blotting, washing & staining
- Tri-directional motion for thorough, gentle mixing

- Optional stacking platforms available
- Cold room and incubator safe
- 2-year warranty



BlotBoy™ Low Speed 3D Rockers

ITEM #	DESCRIPTION	QTY	PRICE
R3D-710	7.5 x 10 in. Low Speed 3D Rocker, 12 rpm	1/ea	\$485.30
R3D-710ST	7.5 x 10 in. Stacking Platform	1/ea	\$142.67
R3D-1212	12 x 12 in. Low Speed 3D Rocker, 12 rpm	1/ea	\$587.40
R3D-1212ST	12 x 12 in. Stacking Platform	1/ea	\$142.67



- Adjustable speed & tilt for a wide range of applications
- Insta-TILT™ technology for instant adjustment to the tilt angle

- Dual platform accessory doubles work surface
- Cold room and incubator safe
- 2-year warranty

BenchWaver™ 3-D Rocker

ITEM #	DESCRIPTION	QTY	PRICE
R3D-5000	13.5 X 13 in. 3D Rocker, 5 to 105 rpm	1/ea	\$1,548.97
R3D-5000-ST	13 in. Stacking Platform, 5.4 in. Clearance	1/ea	\$162.32
R3D-5000-DMP	13 in. Dimpled Mat	1/ea	\$108.00

MIXING

- Adjustable speed & tilt for a wide range of applications
- Cold room and incubator safe
- Insta-TILT™ technology for instant adjustment to the tilt angle
- Dual platform accessory doubles work surface
- 2-year warranty



PR1MA™ 3D Nutating Shaker

ITEM #	DESCRIPTION	QTY	PRICE
R3D-30	12 x 12 in. 3D Nutating Shaker, 2 to 30 rpm	1/ea	\$1,120.98
BR1000-STACK	12 x 12 in. Stacking Platform	1/ea	\$142.67
BR1000-DIMPLED	12 x 12 in. Dimpled Mat	1/ea	\$49.69

- Ideal for all tube sizes 1.5 to 50 mL
- Includes grooved mat for tub
- Cold room and incubator safe
- No adapters or mounting hardware required
- 2-year warranty

Ideal for a variety of tube sizes



TubeRocker Mini

ITEM #	DESCRIPTION	QTY	PRICE
M2100	12 x 4.5 in. Rocker w/ Grooved Mat, Fixed 24 rpm, 115V	1/ea	\$427.07

- Variable speed has mix/pause/rock function
- Includes tube holders for 11 x 1.5/2.0 mL, 5 x 10-15 mL, 5 x 5-7 mL, and 6 x 50 mL
- Light weight (3 lbs)
- 2-year warranty

Need Temperature Control?
Contact us!



PR1MA™ Roto-Mini Series



ITEM #	DESCRIPTION	QTY	PRICE
RFP2020	10.8 x 6 x 6.5 in., 24 rpm, Fixed Speed	1/ea	\$441.50
RMP2024	10.8 x 6 x 6.5 in., 5 to 70 rpm, Variable Speed	1/ea	\$530.71

MIXING

- Digital speed control up to 1,500 rpm
- 3 mm orbit for thorough mixing in microplates
- Includes platform for up to 4 microplates
- Cold room and incubator safe
- Optional adapters for mixing microtubes
- 2-year warranty



PR1MA™ Orbital Microplate Shaker

ITEM #	DESCRIPTION	QTY	PRICE
BT1500	13.1 x 11.7 x 4.7 in. Shaker, Holds Up to 4 Microplates	1/ea	\$1,168.82
BT1500-A05	20 x 0.5 mL Optional Adapter	1/ea	\$97.50
BT1500-A1520	20 x 1.5/2.0 mL Optional Adapter	1/ea	\$97.50

- Gentle orbital motion, ideal for gel staining
- 13 x 12 in. non-slip platform included, clamps sold separately
- 30 to 300 rpm
- Capacity 5 x 1 L or 14 lbs
- 2-year warranty

PR1MA™ Orbital Shaker

ITEM #	DESCRIPTION	QTY	PRICE
OS-3000	13 x 12 in. Shaker, w/Touch Screen Display 3 x 2 L Flasks or 14 lbs Capacity	1/ea	\$1,901.20
OS-PLAT	14 x 12 in. MAGic Clamp™ Universal Platform	1/ea	\$327.68
H1000-MR-CMB	MAGic Clamp™ Combo, 125 mL, 250 mL, 500 mL Flasks	1/ea	\$54.21



- High-capacity benchtop shaker with 18 x 18 in. platform
- Circular 19 mm orbit
- Recommended to use with MAGic Clamp™ platform
- 2-year warranty



Orbi-Shaker™ XL

ITEM #	DESCRIPTION	QTY	PRICE
OSXL-1818	18 x 18 in. Shaker w/ Touch Screen Display & Rubber Mat Platform	1/ea	\$3,780.10
IS-1010-PLAT	18 x 18 in. MAGic Clamp™ Universal Platform	1/ea	\$421.92
H1000-MR-CMB	MAGic Clamp™ Combo, 125 mL, 250 mL, 500 mL Flasks	1/ea	\$54.21

MIXING

- Redesigned to accept two 1 L Flasks
- Increased temperature control up to 70°C
- 2-year warranty



PR1MA™ Mini Incubated Shaker

ITEM #	DESCRIPTION	QTY	PRICE
IS-300	9.5 x 11.5 in. Mini Incubated Shaker, Ambient + 5°C to 70°C, Max Speed 300 rpm	1/ea	\$3,634.96
IS-PLAT	9.5 x 11.5 in. MAGic Clamp™ Platform For PR1MA™ Mini Incubated Shaker	1/ea	\$302.05
H1000-MR-CMB	MAGic Clamp™ Combo, 125 mL, 250 mL, 500 mL Flasks	1/ea	\$54.21


Platforms, clamps and racks available for 1.5 mL to 2 L and microplates

- Exceptional temperature uniformity
- Touch screen control over time, speed and temperature
- Horizontally circular 19 mm orbit for aeration & mixing
- 2-year warranty

10 L Incubated Shaker

ITEM #	DESCRIPTION	QTY	PRICE
IS-1010	18 x 18 in. Incubated Shaker, Ambient +5°C to 70°C, Max Speed 300 rpm	1/ea	\$7,570.26
IS-1010R	18 x 18 in. Refrigerated Incubated Shaker, Ambient -15°C to 60°C, Max Speed 300 rpm,	1/ea	\$10,701.98
IS-1010-PLAT	18 x 18 in. MAGic Clamp™ Universal Platform	1/ea	\$421.92
H1000-MR-CMB	MAGic Clamp™ Combo, 125 mL, 250 mL, 500 mL Flasks	1/ea	\$54.21

Platforms, clamps and racks available for 1.5 mL to 2 L and microplates



MAGic Clamp™ system allows for fast, customizable sharing options!



Never screw in a clamp again!

- Large white ceramic top plate
- Digital speed and temperature control
- Safety LED indicates hot surface
- Available with stir capacities up to 15 L
- 2-year warranty



Digital Hotplates, Stirrers and Hotplate Stirrers

ITEM #	DESCRIPTION	QTY	PRICE
H3710-H	10 x 10 in. Digital Hotplate	1/ea	\$567.32
H3710-S	10 x 10 in. Magnetic Stirrer	1/ea	\$567.32
H3710-HS	10 x 10 in. Digital Hotplate Magnetic Stirrer	1/ea	\$789.85
H3770-H	7 x 7 in. Digital Hotplate	1/ea	\$326.05
H3770-S	7 x 7 in. Magnetic Stirrer	1/ea	\$326.05
H3770-HS	7 x 7 in. Digital Hotplate Magnetic Stirrer	1/ea	\$495.43

CENTRIFUGES

- Eight position microtube rotor
- A second rotor for 0.2 mL PCR tubes and strips is also included
- Storage compartment for PCR rotor
- 2-year warranty

PR1MA™ 8 Mini Centrifuge



Storage compartment for PCR strip rotor

ITEM #	DESCRIPTION	QTY	PRICE
MF8-R	8 Mini, Blue, 8 x 1.5/2.0 mL or 32 x 0.2 mL Tubes or 4 x PCR, Max Speed 6,000 rpm	1/ea	\$346.90
MF8-B	8 Mini, Clear, 8 x 1.5/2.0 mL or 32 x 0.2 mL Tubes or 4 x PCR, Max Speed 6,000 rpm	1/ea	\$346.90
MF8-G	8 Mini, Green, 8 x 1.5/2.0 mL or 32 x 0.2 mL Tubes or 4 x PCR, Max Speed 6,000 rpm	1/ea	\$346.90
MF8-P	8 Mini, Purple, 8 x 1.5/2.0 mL or 32 x 0.2 mL Tubes or 4 x PCR, Max Speed 6,000 rpm	1/ea	\$346.90
MF8-C	8 Mini, Red, 8 x 1.5/2.0 mL or 32 x 0.2 mL Tubes or 4 x PCR, Max Speed 6,000 rpm	1/ea	\$346.90

Includes both microtube and PCR strip rotors



COMBI-Rotor™ accepts tubes and strips!

- Twice the capacity of traditional mini centrifuges
- Starts & stops with closing/opening of the lid
- Unique COMBI-Rotor™ is all that is required for running 12 microtubes and 4 PCR strips simultaneously
- 2-year warranty

PR1MA™ 12 Mini Microcentrifuge

ITEM #	DESCRIPTION	QTY	PRICE
MF12	12 Mini, Capacity 12 x 1.5 mL Tubes, Max Speed 5,500 rpm	1/ea	\$462.04

- Includes rotor for all common centrifuge tubes up to 5 mL
- Digital control up to 5,500 rpm / 1,690 xg
- Conserves valuable bench space
- Includes adapters for common blood tubes (12 x 75 mm)
- 2-year warranty

MyFuge™ 5D

NEW!

ITEM #	DESCRIPTION	QTY	PRICE
C2205	5 Micro w/ Combo Rotor, Capacity 4 x 5 mL & 4 x 1.5/2.0 mL, Max Speed 5,500 rpm	1/ea	\$595.00



Spins 5 mL tubes!
See tubes on page 12!

CENTRIFUGES



- Quickly spins down droplets in PCR plates (less than 20 seconds)
- Accepts skirted, non-skirted, standard, and most microtiter PCR plates
- 50% smaller than traditional centrifuges
- Unique rotor design prevents sample spillage
- 2-year warranty

PRIMA™ Microplate Microcentrifuge



ITEM #	DESCRIPTION	QTY	PRICE
PCRCENT	Microplate Centrifuge w/ Rotor & Plate Carriers, Capacity 2 x 2 PCR Plates or 24 x 0.2 mL PCR strips, Max Speed 2,550 rpm	1/ea	\$894.54
PCRCENT-A02	Adapter Plate for 96 x 0.2 mL Tubes or 12 PCR Strips	2/pk	\$68.60
PCRCENT-DWMP	Rotor for Deep Well Plates (Up to 35 mm Tall)	1/ea	\$235.96

- Unique rotor for 12-position PCR strips
- Starts and stops with opening/closing the lid
- Tiny benchtop footprint and affordability
- Place a StripSpin™ at the bench, thermal cycler & water bath for convenience
- 2-year warranty

StripSpin™ 12D

ITEM #	DESCRIPTION	QTY	PRICE
C2248	Mini 12-Strip PCR Tubes Centrifuge, Capacity 4 x 12 PCR Strips, Max Speed 5,500 rpm	1/ea	\$595.00



- Full color touch screen
- COMBI-Rotor™ accepts tubes and strips
- Compact, only 9 x 12 x 7.75 inches
- Variable speed control, set in RPM or RCF
- 2-year warranty



PRIMA™ MC24 Touch Microcentrifuge

ITEM #	DESCRIPTION	QTY	PRICE
MC24	High Speed Microcentrifuge w/ COMBI-Rotor™, Capacity Varies, Max Speed 13,500 rpm	1/ea	\$2,693.23



COMBI-Rotor™ accepts tubes and strips!

CENTRIFUGES



- 12 tube capacity / up to 16,100 xg
- Less than 9 inches wide
- Quick deceleration
- Whisper quiet operation
- 2-year warranty

PR1MA™ High Speed Microcentrifuge

ITEM #	DESCRIPTION	QTY	PRICE
MF1612	High Speed 15,500 rpm, 12 x 1.5 mL Rotor Microcentrifuge	1/ea	\$1,715.69

- Digital speed control up to 5,000 rpm
- Safety lid lock, prevents opening lid until rotor stops
- Unique airflow design keeps samples cool
- 8 tube capacity, for tubes up to 15 mL
- 2-year warranty

LC-8 & LC-8 Plus Compact Centrifuges

ITEM #	DESCRIPTION	QTY	PRICE
MSLD8	LC-8 w/ Rotor, Capacity 8 x 15 mL, Max Speed 3,500 rpm	1/ea	\$560.26
MSLD8P	LC-8 Plus w/ Rotor, Capacity 8 x 15 mL, Max Speed 5,000 rpm	1/ea	\$809.81



- COMBI-Rotor™ for tubes and strips
- Digital control up to 7,000 rpm/2,739 xg!
- Near silent operation
- Timer up to 99 mins.
- Twice the capacity of traditional mini centrifuges
- 2-year warranty



NEW!

MyFuge™ 12 Plus

ITEM #	DESCRIPTION	QTY	PRICE
C2212	Mini Microcentrifuge, 12 x 1.5/2.0 mL, 32 x 0.2 mL, 4 x PCR Strips	1/ea	\$649.00
C1008-A2-6	Adapters for 0.2 mL Tubes	6/pk	\$27.53
C1008-A5-6	Adapters for 0.5 mL Tubes	6/pk	\$27.53

CENTRIFUGES

- Powerful refrigeration system, Quick-COOL to 4°C in 8 minutes
- Includes easy access rotor (24 x 1.5 mL) up to 17,135 x g

- Exceptionally quiet and compact
- Optional StripSpin™ adapter for 0.2 mL tubes and strips
- 2-year warranty



Best Seller!

Prism™ R Refrigerated Microcentrifuge

ITEM #	DESCRIPTION	QTY	PRICE
C2500-R	Refrigerated Microcentrifuge, Capacity 24 x 1.5/2.0 mL Tubes, Max Speed 13,500 rpm	1/ea	\$5,580.24

- Included COMBI-Rotor™ accepts tubes and strips (lid included)

- Powerful refrigeration, as low as -20°C
- Up to 16,800 xg!

- 2-year warranty
- Capacity:
24 x 1.5/2.0 mL
4 x PCR strips (32 x 0.2 mL)
24 x spin column



NEW!

MC-24R Refrigerated Microcentrifuge

ITEM #	DESCRIPTION	QTY	PRICE
C2417-R	Refrigerated H/S Microfuge w/ COMBI-rotor, Max Speed 13,500 rpm	1/ea	\$5,780.43
C2417-2420H	Hermetically Sealed Rotor, 24 x 1.5/2.0 mL Tubes	1/ea	\$858.23
C2417-3202	PCR Strip Rotor, 4 x 32 x 0.2 mL Strips	1/ea	\$810.90
C2417-RCMB	COMBI-Rotor™, 24 x 1.5/2.0 mL &; 2 x PCR Strips (16 x 0.2 mL)	1/ea	\$576.30

- Temp Range: Room temp. to 4°C in <10 min

- 5-year warranty

Z216-MK Refrigerated Microcentrifuge

ITEM #	DESCRIPTION	QTY	PRICE
Z216-K-CMB	High Speed w/ Combi 24 x 1.5 & 4 PCR Strip Rotor, Capacity 44 x 1.5/2.0 mL, Max Speed 15,000 rpm	1/ea	\$8,677.01
Z216-K-BND44	High Speed w/ 44 x 1.5 mL Rotor, Capacity 44 x 1.5/2.0 mL, Max Speed 15,000 rpm	1/ea	\$7,981.55

** For more rotors and accessories visit midsci.com **

Looking for other centrifuge options? Let us know!



24 x 1.5/2.0 mL Rotor
Z216-4420



COMBI-Rotor™
Z216-2420S



HOODS



- Energy Efficient: DC ECM blower saves 70% energy; standby mode reduces power by 60%
- Quiet Operation: Stable airflow with minimal noise, even during voltage fluctuations or filter loading
- Superior Filtration: ULPA filter (10x more efficient than HEPA) provides ISO Class 3 protection
- Easy Maintenance & Safety: One-piece interior, removable tray for cleaning; airflow sensor and UV timer for safety
- 5-Year Warranty

Esco® Airstream® Biosafety Cabinet

ITEM #	DESCRIPTION	QTY	PRICE
AR2-3S9-G4	3 ft. Class II, Type A2, 8-inch Opening, 110-130 VAC, 50/60 Hz	1/ea	Contact for Quote
AR2-4S9-G4	4 ft. Class II, Type A2, 8-inch Opening, 110-130 VAC, 50/60 Hz	1/ea	Contact for Quote
AR2-5S9-G4	5 ft. Class II, Type A2, 8-inch Opening, 110-130 VAC, 50/60 Hz	1/ea	Contact for Quote
AR2-6S9-G4	6 ft. Class II, Type A2, 8-inch Opening, 110-130 VAC, 50/60 Hz	1/ea	Contact for Quote

- User-Friendly Touchscreen: Intuitive controls with clear safety info and quickstart mode
- Energy Efficient: Saves up to 70% energy with a DC ECM motor and low power consumption
- Adjustable LED Lighting: Optimal brightness (1200 lux) with dimmable settings for comfort
- Ergonomic Design: Easy sash opening, tempered glass, and stainless steel armrest for comfort
- Superior Filtration: ULPA filter (99.999% efficiency) for ISO Class 3 air quality
- Durable & Easy to Maintain: Stainless steel components and simple drain pan cleaning
- Antimicrobial Coating: Isocide™ coating eliminates 99.9% of bacteria within 24 hours
- 5-Year Warranty

NEW!

Airstream® NS G4 Class II Type A2 Bio Safety Cabinet- AC2 Models

ITEM #	DESCRIPTION	QTY	PRICE
AC2-3S9-NS-G4	3 ft. Type A2, Class II, 8-inch Opening, 110-130 VAC, 50/60 Hz	1/ea	Contact for Quote
AC2-4S9-NS-G4	4 ft. Type A2, Class II, 8-inch Opening, 110-130 VAC, 50/60 Hz	1/ea	Contact for Quote
AC2-5S9-NS-G4	5 ft. Type A2, Class II, 8-inch Opening, 110-130 VAC, 50/60 Hz	1/ea	Contact for Quote
AC2-6S9-NS-G4	6 ft. Type A2, Class II, 8-inch Opening, 110-130 VAC, 50/60 Hz	1/ea	Contact for Quote



Looking for other options?
Let us know!

SPECTROPHOTOMETERS



- Standalone or remote operated
- Accurate, consistent results
- Easy sample application
- No maintenance and recalibration required ever
- Effortless integration into any LIMS system
- GMP, GLP, GxP Ready 21 CFR Part 11

Implen NanoPhotometer® Spectrometers

ITEM #	DESCRIPTION	QTY	PRICE
N50-TOUCH	NanoPhotometer NP50 Touch UV/Vis, 200-650 nm	1/ea	\$8,900.00
N60-TOUCH	NanoPhotometer NP60 Touch UV/Vis, 200-900 nm	1/ea	\$11,200.00
NP80-TOUCH	NanoPhotometer NP80 Touch UV/Vis, 200-900 nm	1/ea	\$13,250.00

- Quick and accurate results from a sample droplet
- Touch screen control panel, no PC or software required
- Built-in printer and cuvette port for OD600 readings
- XF model includes fluorometer for sensitive DNA quantification

SmartDrop™ Nano SpectroPhotometers

ITEM #	DESCRIPTION	QTY	PRICE
NS1000	SmartDrop™ L Spectrophotometer	1/ea	\$6,150.90
NS1010	SmartDrop™ X UV-Vis Spectrophotometer	1/ea	\$11,539.50
NS1020	SmartDrop™ XF UV-Vis Spectrophotometer, Fluorometer	1/ea	\$13,639.50



- Higher concentration range, accuracy and reproducibility!
- User-friendly, multilingual interface
- High precision versatile measurements, linear up to 4 A
- Stand-alone operation
- Extremely lightweight, compact and portable
- No maintenance or recalibration needed ever!



New & Improved!

Implen OD600®

ITEM #	DESCRIPTION	QTY	PRICE
OD600-10-US	Photometric, OD600: 0 - 4 A, McFarland: 0 - 16 MFU, k 600 nm	1/ea	\$1,890.00
DC10	Disposable Cuvettes, 1/10 Dilutions	96/pk	\$130.00

INCUBATORS



- Digital control of temperature and CO₂
- Temp range ambient -10°C to 55°C
- Accepts bottles and flasks up to 2 L
- Includes 2 stainless steel shelves
- Dual beam, infrared sensor
- 2-year warranty

PR1MA™ Mini CO₂ Incubator

ITEM #	DESCRIPTION	QTY	PRICE
INCO2-H	0.8 CuFt Incubator, 13 x 14.5 in. Footprint (20 L)	1/ea	\$4,586.60
H2300-REG	Optional CO ₂ Gas Regulator	1/ea	\$391.31

- Digital temperature control
- Personal-sized, 13 x 14.5 inch footprint (20 L)
- Accepts bottles and flasks up to 2 L
- Internal outlet, for installation of shaker, rockers, etc.
- 2-year warranty

PR1MA™ Mini Digital Incubators

ITEM #	DESCRIPTION	QTY	PRICE
IN20-H	Heat Only (Ambient +5 to 60°C), Incl. 2 Shelves	1/ea	\$605.50
IN20-HC	Mini Digital Incubator - Heat and Cool (Ambient +5 to 60°C), Incl. 2 Shelves	1/ea	\$921.04
H2200-SH	Extra Shelf, Plastic	1/ea	\$48.40



- Innovative chamber design features smooth rounded internal corners for easy cleaning
- Infra-Red(IR) CO₂ sensor for precise and reliable CO₂ detection
- Stackable design up to 3 units
- 25 mm and 50 mm orbit models available
- 2-year warranty

AEOLUS Compact CO₂ Incubator Shakers

ITEM #	DESCRIPTION	QTY	PRICE
IS-AMA	970 x 780 x 530 mm Incubator Shaker, 7 in. High-Res Touch Screen, Air Jacket, 25 mm (1.02 in.) Orbit	1/ea	\$17,550.00
IS-A5MA	970 x 780 x 530 mm Incubator Shaker, 7 in. High-Res Touch Screen, Air Jacket, 50 mm (2 in.) Orbit	1/ea	\$19,500.00
IS-ACMA	970 x 780 x 530 mm CO ₂ Incubator Shaker, 7 in. High-Res Touch Screen, Air Jacket, 25 mm (1.02 in.) Orbit	1/ea	\$22,725.00
IS-A5CMA	970 x 780 x 530 mm CO ₂ Incubator Shaker, 7 in. High-Res Touch Screen, Air Jacket, 50 mm (2 in.) Orbit	1/ea	\$23,730.00

INCUBATORS

- Contamination control with long life ULPA filters and 90°C moist heat decontamination cycle
- Direct heat & air jacket design with gentle yet effective forced convection
- Intuitive interfaces and intelligent data and event logger
- 2-year warranty



CelCulture® CO₂ Incubator

ITEM #	DESCRIPTION	QTY	PRICE
CCL-170B-9-SP	170 L CelCulture® CO ₂ Incubator, IR Sensor, CO ₂ Control, ULPA, High Temp Decontamination	1/ea	Contact for Quote

- 6-sided heating, ambient +5°C to 60°C
- Quick recovery of temperature and CO₂
- Internal microscope option
- Monitor cells on the display or remotely on a smart device
- 2-year warranty

SureTherm™ CO₂ Incubators



ITEM #	DESCRIPTION	QTY	PRICE
H3565-45	45 L CO ₂ Incubator w/ Two Shelves	1/ea	\$6,277.51
H3565-180	180 L CO ₂ Incubator w/ Three Shelves	1/ea	\$7,582.86
H3565-180HD	180 L CO ₂ Incubator w/ High Heat Decontamination	1/ea	\$9,480.04
H3565-180HDO2	180 L CO ₂ Incubator w/ High Heat Decontamination, Split Window Door and O ₂ Control	1/ea	\$13,101.51
H3565-45-LCI	Integrated Microscope Shelf, for 45 L Incubator	1/ea	\$4,719.50
H3565-180-LCI	Integrated Microscope Shelf, for 180 L Incubator	1/ea	\$4,916.81

- 7 in. High-Res color touch screen with intuitive graphic interface
- Capacity is 6 x 5 L
- Infra-Red(IR) CO₂ sensor for precise and reliable CO₂ detection
- Optional WiFi and RS485 monitoring and control
- Two sets of UV lights for easy sterilization
- Stackable design up to 3 units
- 2-year warranty



AEOLUS Stackable 302 L CO₂ Incubator Shaker

ITEM #	DESCRIPTION	QTY	PRICE
IS-18A	302 L Stackable Incu Shaker, 25 mm Orbit	1/ea	\$23,010.00
IS-18CA	302 L Stackable CO ₂ Incu Shaker, 25 mm Orbit	1/ea	\$33,240.00
IS-18B5A	302 L Stackable Incu Shaker, 50 mm Orbit	1/ea	\$35,623.14
IS-18C5A	302 L Stackable CO ₂ Incu Shaker, 50 mm Orbit	1/ea	\$34,240.00

PCR EQUIPMENT



- Gradient capabilities- ideal for protocol optimization
- Intuitive touchscreen control- works even when wearing gloves!
- 2-year warranty
- No need for interchangeable blocks - multi-format accepts:
 - 0.2 mL tubes or strips (96)
 - 0.5 mL tubes (39)
 - 96-well PCR plate

PRIMA™ TC9639 Thermal Cycler

ITEM #	DESCRIPTION	QTY	PRICE
BSTC-96	10.2 x 14.2 in. Cycler w/ Combination Block, 96 x 0.2 mL	1/ea	\$5,995.00
BSTC-384	10.2 x 14.2 in. Cycler w/ 384-Well Block	1/ea	\$6,979.04
BSTC-AIS	Combo Adapter (4 slides)	1/ea	\$587.20

- Unique flexible programming with FlexTemp™ technology
- Fast ramp rates: 5°C heating, 3.5°C cooling
- No condensation after overnight cooling at 4°C
- Capacity:
 - 1 x 96 well plate
 - 12 - 8 x 0.2 mL strip tubes
 - 96 x 0.2 mL tubes
- 3-year warranty

Best Seller!



MultiGene™ Optimax Thermal Cycler

ITEM #	DESCRIPTION	QTY	PRICE
TC9610	96 Sample Thermal Cycler, 115V Model	1/ea	\$4,328.89

The Azure Cielo is a 96-well real-time PCR machine that features innovative and high-performance optical technology with broad application compatibility.



- Unparalleled data uniformity
- Shorter run times
- Data emailed after run completion
- Multiplex qPCR of up to six-targets

Real Time qPCR Azure Systems

ITEM #	DESCRIPTION	QTY	PRICE
AIQ030	Cielo 3 qPCR System, 0.1 mL x 96	1/ea	\$27,850.00
AIQ060	Cielo 6 qPCR System, 0.1 mL x 96	1/ea	\$34,850.00
AIQ1000	Cielo 3 Software Plate Setup, PCR Compatible	1/ea	\$1,499.00

SWABS AND SALIVA KITS

- Easy-to-use with anti-spill funnel and pictorial user instructions

- Total DNA yield from 2 mL of saliva of 100 µg stable at room temperature for over 5 years
- RNA, microbiome and stool collection also available



Isohelix™ GeneFix DNA/RNA Collection and Stabilization

ITEM #	DESCRIPTION	QTY	PRICE
GFX-02-12	2 mL DNA Saliva Collector	12 Samples	\$314.37
GFX-02-50	2 mL DNA Saliva Collector	50 Samples	\$1,129.38

- Larger textured swab reliably collects over 10 µg of DNA

- RapiDri™ Kits integrate barcodes into collection and offer 3 months of DNA Stability at room temperature

- Dri-Capsules work with 5 mL tube for over 3 years of room temperature DNA Stability

Isohelix™ DNA/RNA Buccal Swabs and Dri-Capsules

ITEM #	DESCRIPTION	QTY	PRICE
SK-1S	9 mm Swab w/ 5 mL Tube & Cap, Indiv. Wrap	100/pk	\$220.08
SK-2S	9 mm Swab w/ 2 mL Tube & Special Release Cap, Indiv. Wrap	100/pk	\$297.87
SK-3S	9 mm Swab, Indiv. Wrap	250/pk	\$484.05
SK-4S	9 mm Swabs, 2/pk, Indiv. Wrap	250/cs	\$880.07
RD-01	9 mm Swab w/ RapiDri™ Pouch & Barcode, Indiv. Wrap	200/pk	\$772.00
SGC-50	Dri-Capsules, Individually Foil Wrapped	50/pk	\$194.98



All swabs are Ethylene Oxide treated

- BuccalFix is a liquid lysis and stability buffer for up to 2 years of DNA stability at room temperature



Isohelix™ BuccalFix DNA Stabilization and Lysis Kits

ITEM #	DESCRIPTION	QTY	PRICE
BFX-S1-05-50	2 mL Screw Cap Tubes Pre-Filled w/ 0.5 mL Stabilization Buffer	50/pk	\$346.67
BFX-25	25 mL BuccalFix Stabilization Buffer Bulk Pack	50/pk	\$194.98

ELECTROPHORESIS

New!

Azure Electrophoresis Equipment


ITEM #	DESCRIPTION	QTY	PRICE
AC4200	Vertical Gel Running System, 10 x 8 cm	1/ea	\$994.00
AC4201	Aqua Transfer Cell, 7.5 x 10 cm	1/ea	\$639.00
AC4203	Aqua Power Supply, 600V, 1.2 A, 500W	1/ea	\$2,345.00



- Molded, durable construction
- Rubber gates for easy, tapeless gel casting
- Adjustable height combs
- Loading guides aid in well visualization

myGel InstaView System

ITEM #	DESCRIPTION	QTY	PRICE
E1201	Electrophoresis System w/Transilluminator, 10.5 x 10 cm or 0.5 x 6 cm	1/ea	\$1,016.73
E1200	Gel Tank w/ Transilluminator, 10.5 x 10 cm or 0.5 x 6 cm	1/ea	\$778.51



- Molded, durable construction
- Rubber gates for easy, tapeless gel casting
- Adjustable height combs
- Loading guides aid in well visualization



Enduro™ Horizontal Electrophoresis System

ITEM #	DESCRIPTION	QTY	PRICE
E1007-7	7 x 7 cm (Gel Tank), Buffer Vol 225 mL, Max Sample Cap 32	1/ea	\$453.17
E1007-10	7 x 10 cm (Gel Tank), Buffer Vol 225 mL, Max Sample Cap 64	1/ea	\$471.96
E1010-10	10 x 10 cm (Gel Tank), Buffer Vol 300 mL, Max Sample Cap 100	1/ea	\$497.50
E1015-10	15 x 10 cm (Gel Tank), Buffer Vol 500 mL, Max Sample Cap 140	1/ea	\$576.85
E1015-15	15 x 15 cm (Gel Tank), Buffer Vol 500 mL, Max Sample Cap 210	1/ea	\$746.03
E1020-20	20 x 20 cm (Gel Tank), Buffer Vol 1,200 mL, Max Sample Cap 450	1/ea	\$779.89

Visit midsci.com for more options!

ELECTROPHORESIS



Galileo gel boxes and accessories are fully compatible with Owl systems!

Galileo Electrophoresis Systems

ITEM #	DESCRIPTION	QTY	PRICE
81-2025	20 x 25 cm Horizontal Electrophoresis System	1/ea	\$812.67
81-2325	23.5 x 25 cm Extra Wide Format Gel System	1/ea	\$849.75
81-2340-BEP	23.5 x 40 cm Extra Wide Large Format Gel System	1/ea	\$915.67

Contact us for more options!

- Complete system, includes gel box, power supply and casting sets
- Quick and easy gel casting, no tape, clamps or spacers required
- Compact design, footprint: 7 x 5.1 in.

PR1MA™ Mini Gel Electrophoresis System

ITEM #	DESCRIPTION	QTY	PRICE
EMINIGEL	Horizontal System, One 10.5 x 6 cm or two 5 x 6 cm Gels, Variable Teeth	1/ea	\$645.92
E1101-COMB1	22/12 Reversible Comb	2/pk	\$9.84
E1101-COMB2	18/10 Reversible Comb	2/pk	\$12.32
E1101-CS1	10.5 x 6 cm Gel Casting Stand	1/ea	\$49.43
E1101-CS2	5 x 6 cm Gel Casting Stand	1/ea	\$79.93



- Touch screen control
- Dual orientation, auto rotating display
- Bright color display
- 2-year warranty




PR1MA™ Touch Power Supply

ITEM #	DESCRIPTION	QTY	PRICE
PS2300	Compact Electrophoresis Power Supply, 100 W	1/ea	\$679.59

TRANSILLUMINATORS

- Space-saving design with three models to fit various laboratory needs and gel sizes
- Hinged UV-blocking safety cover for added protection
- 302 nm UV wavelength, scratch-resistant UV filter glass, and optimal for EtBr, SmartGlow™, and SYBR® Stains

Accuris™ UV Transilluminators



ITEM #	DESCRIPTION	QTY	PRICE
BE3000	16 x 19 cm, Peak at 302 nm, Low/High Intensity, 115V	1/ea	\$1,287.09
E3100	16.5 x 13.5 cm, Peak at 302 nm, 115V	1/ea	\$808.47
E3200	21 x 26 cm, Peak at 302 nm, Variable Intensity, 115V	1/ea	\$1,250.00



Accuris™ SmartBlue™ Transilluminators

ITEM #	DESCRIPTION	QTY	PRICE
SB4500	Bench Top Transilluminator, 465 nm, Includes Amber Cover 115V	1/ea	\$675.22
E4101	Mini Transilluminator, 465 nm, Includes Amber Cover 115V	1/ea	\$496.56



- Excellent lot-to-lot consistency
- Exceptional clarity with extremely low backgrounds
- Classic push-button cassettes
- Constructed from lightweight aluminum & nylon
- Special curved door gives exceptional screen-film contact

Autoradiography Film and Cassettes



ITEM #	DESCRIPTION	QTY	PRICE
BX57	5 x 7 in. PR1MA™ Autoradiography & Chemiluminescence Film	100/pk	\$72.10
BX810	8 x 10 in. PR1MA™ Autoradiography & Chemiluminescence Film	100/pk	\$94.76
CPB57	5 x 7 in. Aluminum Autoradiography Cassette	1/ea	\$93.38
CPB810	8 x 10 in. Aluminum Autoradiography Cassette	1/ea	\$104.42
CPB1012	10 x 12 in. Aluminum Autoradiography Cassette	1/ea	\$120.52

IMAGERS

SmartDoc™ Gel Imaging

ITEM #	DESCRIPTION	QTY	PRICE
E5001-SDB	SmartDoc™ Complete System	1/ea	\$823.89
E5001-SD	SmartDoc™ Enclosure	1/ea	\$471.32
E4000-VG1	SmartBlue Plus Viewing Glasses	1/ea	\$53.48



ENDURO™ GDS II Gel Documentation System

ITEM #	DESCRIPTION	QTY	PRICE
GDST2-1302	Enduro™ GDS Touch II with 302 nm	1/ea	\$10,176.23
GDST2-1365	Enduro™ GDS Touch II with 365 nm	1/ea	\$10,176.23

- Iris Camera: 9.2mp, -60°C cooling, 16 bit camera, Lens f/0.75
- Chromascan technology for up to 9 channels
- Includes 15.6 in. touch screen
- Plug and play upgradeable
- Superior imaging especially for IR and NIR fluorescence
- Optional LED UV transilluminator
- Unlimited license acquisition, editing and analysis software included



New! UVITEC Alliance IRIS

ITEM #	DESCRIPTION	QTY	PRICE
1541-1000-1	Chemiluminescence & Fluorescence Western Blot Imaging System, Motorized Lens f/0.75	1/ea	Contact for Quote



- Fully automated system
- Real on-board computer
- Giant 12.2 in. HD touch-screen
- 6-megapixel native resolution
- Voice Control® option
- No learning curve
- 65,535 gray levels
- 26 x 21 cm imaging area

UVIDOC-HD6 Plus 26MX

ITEM #	DESCRIPTION	QTY	PRICE
164126031	Stand-Alone Gel Doc for DNA Gel Imaging	1/ea	Contact for Quote

HOMOGENIZERS

Bullet Blender 5E Pro

- Homogenizes 12 samples in 5 mL tubes (adapter for smaller tubes)
- Powerful motor, touch screen, air cooling™



Bullet Blender Storm Pro

- Homogenizes 24 samples in 1.5 mL tubes
- Powerful motor, touch screen, air cooling™

Bullet Blender Lite

- Homogenizes 12 samples in 1.5 mL GATOR™ tubes
- Touch screen, air cooling™

Bullet Blender® Storm Homogenizers

ITEM #	DESCRIPTION	QTY	PRICE
BT5E	5E Pro, Air Cooling, 12 Samples, 1.5 mL, 2 mL, 5 mL	1/ea	\$5,225.00
BT24M	Storm Pro, Air Cooling, 24 Samples, 1.5 mL	1/ea	\$5,125.00
BT12LT	Lite, Air Cooling, 12 Samples, 1.5 mL	1/ea	\$3,995.00

2-year warranty on all units

- Ideal for tough tissues (e.g., mouse femur, tumor, skin)
- Compatible with dry ice to prevent sample heating

- Consistent, high-yield homogenization with minimal variation
- No probes or blades, eliminating cross-contamination
- No cleanup required between runs

New!

Bullet Blender® Gold+ Homogenizers

ITEM #	DESCRIPTION	QTY	PRICE
BB50AUP	4°C to 40°C, Dry Ice and Liquid Nitrogen Cooling, 8 Samples, 50 mL	1/ea	\$9,885.00
BB5EAUP	4°C to 40°C, Dry Ice and Liquid Nitrogen Cooling, 12 Samples, 1.5 mL, 2 mL, 5 mL	1/ea	\$9,985.00
BB24AUP	4°C to 40°C, Dry Ice and Liquid Nitrogen Cooling, 24 Samples, 1.5 mL	1/ea	\$9,685.00

2-year warranty on all units



- Extra-powerful, efficient homogenization
- Processes up to 4 tough or 8 soft samples

- Handles up to 3.5 g tissue or 3.5 mL cells with 20 mL buffer
- Compatible with HIPPO™ 50 mL tubes
- Air cooling™ prevents overheating

New!

Bullet Blender® 50DX+

ITEM #	DESCRIPTION	QTY	PRICE
BB50DXP	Air Cooling, 4-8 Samples, 50 mL	1/ea	\$7,085.00

2-year warranty



BATHS



- Temp range ambient +5°C to 100°C
- 2-year warranty
- Blocks available from 0.2 mL to 50 mL to plates to custom and more!

PR1MA™ Mini Dry Baths

ITEM #	DESCRIPTION	QTY	PRICE
MB100	Mini Digital Dry Bath	1/ea	\$320.01
MB100-HL	Mini Digital Dry Bath, Heated Lid	1/ea	\$530.36
MB300	Mini Digital Dry Bath, Heat and Cool	1/ea	\$854.76

Blocks sold separately

- Temp range ambient +5°C to 105°C
- 2-year warranty

PR1MA™ Digital Dry Baths

ITEM #	DESCRIPTION	QTY	PRICE
LDMB1	Dry Bath, 1 Block Capacity	1/ea	\$554.82
LDMB2	Dual Dry Bath, 2 Block Capacity	1/ea	\$690.75

Blocks Not Included



Contact us for other block options!



- Works with thermal beads or water
- Quick-Cal™ simple in lab calibration function
- Runs constant - No burn out!
- Temp Control:
 - Ambient +5°C to 80°C (water mode)
 - Ambient +5°C to 70°C (bead mode)
- 2-year warranty

BeadBath Duo™

ITEM #	DESCRIPTION	QTY	PRICE
B2402	2 L Combination Water & Bead Bath	1/ea	\$847.96
B2410	10 L Combination Water & Bead Bath	1/ea	\$1,446.07

Need thermal beads? We've got 'em. Give us a call!

Lab Armor® Bath Beads

ITEM #	DESCRIPTION	QTY	PRICE
42370-002	2 Liter Beads, (Beads Only)	1 x 2L Bag	\$233.06
42370-004	4 Liter Beads, (Beads Only)	1 x 4L Bag	\$466.13
42370-005	5 Liter Beads, (Beads Only)	1 x 5L Bag	\$582.66
42370-008	8 Liter Beads, (Beads Only)	2 x 4L Bags	\$932.27
42370-140	140 Liter Beads, (Beads Only)	2 X 70L Barrels	\$16,314.20



Wave farewell to water bath mess, Beads are here to save the day!

COLD STORAGE

- Precise temperature control with better performance than mechanical thermostats
- One interior duplex outlet and side all port hole for easier setup

- Fan-driven air circulation for uniform temperature and quick recovery
- Uses natural HC refrigerants, reducing GWP and meeting EPA/SNAP standards
- Cycle defrost, audible/visual alarms, remote contacts, and two temperature probes

Glass Door Chromatography Refrigerators

ITEM #	DESCRIPTION	QTY	PRICE
ABT-HC-CS-16	16 CuFt, 1°C to 10°C, Swing Door, 115 V/60 HZ	1/ea	\$4,836.00
ABT-HC-CS-23	23 CuFt, 1°C to 10°C, Swing Door, 115 V/60 HZ	1/ea	\$5,534.00
ABT-HC-CS-26	26 CuFt, 1°C to 10°C, Swing Door, 115 V/60 HZ	1/ea	\$5,844.00
ABT-HC-CS-33	33 CuFt, 1°C to 10°C, Two Sliding Doors, 115 V/60 HZ	1/ea	\$7,386.00
ABT-HC-CS-36	36 CuFt, 1°C to 10°C, Two Swing Doors, 115 V/60 HZ	1/ea	\$7,242.00
ABT-HC-CS-47	47 CuFt, 1°C to 10°C, Two Sliding Doors, 115 V/60 HZ	1/ea	\$7,982.00
ABT-HC-CS-49	49 CuFt, 1°C to 10°C, Two Swing Doors, 115 V/60 HZ	1/ea	\$7,894.00
ABT-HC-CS-69	69 CuFt, 1°C to 10°C, Three Sliding Doors, 115 V/60 HZ	1/ea	\$9,544.00
ABT-HC-CS-72	72 CuFt, 1°C to 10°C, Three Swing Doors, 115 V/60 HZ	1/ea	\$9,486.00



- Precise temperature control, digital microprocessor for accurate management
- Superior airflow, forced air system for uniform temperature and quick recovery

- Eco-Friendly, uses natural HC refrigerants, meeting EPA/SNAP standards
- Safety features, audible/visual alarms, two probes, and remote contacts for monitoring

Standard Glass Door Laboratory Refrigerator

ITEM #	DESCRIPTION	QTY	PRICE
ABT-HC-LS-23	23 CuFt, 1°C to 10°C, Swing Door, 115 V/60 HZ	1/ea	\$5,032.00
ABT-HC-LS-26	26 CuFt, 1°C to 10°C, Swing Door, 115 V/60 HZ	1/ea	\$5,092.00
ABT-HC-LS-33	33 CuFt, 1°C to 10°C, Two Sliding Doors, 115 V/60 HZ	1/ea	\$6,600.00
ABT-HC-LS-36	36 CuFt, 1°C to 10°C, Two Swing Doors, 115 V/60 HZ	1/ea	\$6,478.00
ABT-HC-LS-47	47 CuFt, 1°C to 10°C, Two Sliding Doors, 115 V/60 HZ	1/ea	\$7,128.00

Explore even more options!
Contact us or visit our website for additional details.

COLD STORAGE

- Features a microprocessor controller, digital display, and dual probes for accurate management of critical samples

- Uses natural hydrocarbon refrigerants (R290) to reduce GWP, meeting EPA/SNAP and UL/ASHRAE standards

- UL/C-UL or ETL Listed, Energy Star certified, with a 2-year warranty (3-year for compressor) and alarms for high/low temperatures



Standard Manual Defrost Laboratory Freezer

ITEM #	DESCRIPTION	QTY	PRICE
ABT-HC-MFP-14	14 CuFt, -15°C to -25°C, Swing Door, 115 V/60 HZ	1/ea	\$4,100.00
ABT-HC-MFP-17	17 CuFt, -15°C to -25°C, Swing Door, 115 V/60 HZ	1/ea	\$4,242.00
ABT-HC-MFP-20	20 CuFt, -15°C to -25°C, Swing Door, 115 V/60 HZ	1/ea	\$4,416.00



- Manual defrost design, built for general-purpose use, with periodic manual defrost to prevent ice build-up
- Eco-Friendly refrigerants, uses natural hydrocarbon refrigerants (R290 or R600a), reducing GWP and meeting EPA/SNAP standards

- Features a digital display and microprocessor controller, offering improved performance over past solid door models
- 2-year warranty parts & labor, 3-year for compressor

Corepoint -20°C Manual Defrost Freezer

ITEM #	DESCRIPTION	QTY	PRICE
NSPF14	14 CuFt, -15°C to -25°C, Swing Door, 115 V/60 HZ	1/ea	\$3,451.00
NSPF20	20 CuFt, -15°C to -25°C, Swing Door, 115 V/60 HZ	1/ea	\$3,691.00
LF30	30 CuFt, -15°C to -25°C, Swing Door, 115 V/60 HZ	1/ea	\$5,305.00



COLD STORAGE

- Precise Temperature Control: -40°C to -86°C, with 1°C increments and clear LED display
- Advanced Alarms & Safety: Multiple alarms, remote monitoring, and UL/CE certified
- Reliable & Durable: High-quality components, frost-reducing inner doors
- Energy-Efficient: VIP insulation, low energy use, Energy Star rated
- User-Friendly: Adjustable shelves, USB data download, quiet operation



Haier -86°C ULT Freezers

ITEM #	DESCRIPTION	QTY	PRICE
DW-86L338J	11.9 CuFt, -40°C to -86°C, Upright, 115 V/60 HZ	1/ea	Contact for Quote
DW-86L578J	20.4 CuFt, -40°C to -86°C, Upright, 115 V/60 HZ	1/ea	Contact for Quote
DW-86L728J	25.7 CuFt, -40°C to -86°C, Upright, 115 V/60 HZ	1/ea	Contact for Quote
DW-86L828JA	29.2 CuFt, -40°C to -86°C, Upright, 220 V/60 HZ	1/ea	Contact for Quote

3-year warranty parts & labor, 5-year warranty compressor



- 10 in. Touchscreen: Intuitive UI for easy sample management
- IoT System: Real-time monitoring, fault diagnosis, and barcode scanning
- Fast Retrieval: One-key, one-code operation for quick access
- Security: Multiple lock options for added sample safety
- Energy-Efficient: <10kWh/day, noise reduction, and advanced insulation

Haier -80°C Smart Frequency ULT Freezers

ITEM #	DESCRIPTION	QTY	PRICE
DW-86L579BPT	20.5 CuFt, -86°C, Upright, 110-230 V/60 Hz	1/ea	Contact for Quote
DW-86L729BPT	26 CuFt, -86°C, Upright, 110-230 V/60 Hz	1/ea	Contact for Quote
DW-86L829BPT	29 CuFt, -86°C, Upright, 208-230 V/60 Hz	1/ea	Contact for Quote
DW-86L959BPT	34 CuFt, -86°C, Upright, 110-230 V/60 Hz	1/ea	Contact for Quote

5-year warranty parts & labor, 10-year warranty compressor

- Energy Efficient: High-efficiency cooling for long-term savings
- Dual Systems: Two independent refrigerants ensure reliable -80°C performance
- Multiple Alarms: Alerts for temperature, power, and other potential issues
- Eco-Friendly: Uses HC refrigerants with zero GWP and low-GWP foam
- Stable Temps: Frequency conversion system ensures temperature consistency

Haier -86°C TwinCool Frequency Conversion ULT Freezers

ITEM #	DESCRIPTION	QTY	PRICE
DW-86L578BPST	20.5 CuFt, -40°C to -86°C, 110-230V/60Hz	1/ea	Contact for Quote
DW-86L728BPST	25.7 CuFt, -40°C to -86°C, 110-230V/60Hz	1/ea	Contact for Quote
DW-86L828BPST	29.2 CuFt, -40°C to -86°C, 110-230V/60Hz	1/ea	Contact for Quote

5-year warranty parts & labor, 10-year warranty compressor



COLD STORAGE

- Temperature Control: Designed for 2 to 8°C storage of medicines, vaccines, and more
- Energy Efficient: Eco-friendly refrigerants (R290, R600a) with low power consumption
- Reliable Cooling: Hermetically sealed compressors and forced-air circulation with auto defrost
- Security & Accessibility: Key lock, adjustable shelves, and flexible storage
- Quiet Operation: Low noise (35-43 dB(A)) and precise digital controls for easy monitoring



HYC-85GD

HYC-85FD



HYC-461GD



HYC-1031GD



HYC-1031FD

Haier Pharmacy Refrigerators

Looking for other options? Contact us!

ITEM #	DESCRIPTION	QTY	PRICE
HYC-85GD	3 CuFt, Upright, Glass Door, N, 16~32°C	1/ea	Contact for Quote
HYC-461GD	16.3 CuFt, Upright, Glass Door, N, 16~32°C	1/ea	Contact for Quote
HYC-639	22.2 CuFt, Upright, Glass Door, ST, 10°C~43°C	1/ea	Contact for Quote
HYC-1031GD	36.4 CuFt, Upright, Glass Door, N, 16~32°C	1/ea	Contact for Quote
HYC-85FD	3 CuFt, Upright, Solid Door, N, 16~32°C	1/ea	Contact for Quote
HYC-461FD	16.3 CuFt, Upright, Solid Door, N, 16~32°C	1/ea	Contact for Quote
HYC-1031FD	22.2 CuFt, Upright, Solid Door, N, 16~32°C	1/ea	Contact for Quote

3-year warranty parts & labor, 5-year warranty compressor

Product Advantages

- Easy to access
- Brand-new intelligent monitoring system
- Reduces frosting and freezing

Key Design Features

- High degree of safety
- Ergonomic design

NEW!

Haier Biobank Series LN₂ Container



ITEM #	DESCRIPTION	QTY	PRICE
CRYOBIO-11Z	200 L Volume, Holds 11,000 or 13,000 Vials	1/ea	Contact for Quote
CRYOBIO-20Z	340 L Volume, Holds 20,350 or 24,050 Vials	1/ea	Contact for Quote
CRYOBIO-34Z	550 L Volume, Holds 33,550 or 39,650 Vials	1/ea	Contact for Quote

2-year warranty parts & labor, 5-year warranty vacuum

COLD STORAGE



- Rugged and durable aluminum construction
- Vapor phase cryogenic storage
- Straw and cryovial storage options
- Designed to ensure no LN₂ spillage
- 3-year vacuum warranty

Haier LN₂ Dryshipper Series for Transportation

ITEM #	DESCRIPTION	PRICE
YDH-3	3 L, 132 Straws (0.5 mL), Incl. Lockable Cap, Protective Jacket	Contact for Quote
YDH-6-80	6 L, 374 Straws (0.5 mL), Incl. Lockable Cap, Protective Jacket	Contact for Quote
YDH-10-125-F	10 L, 100 Vials (2 mL), Incl. Rack, Cryoboxes, Cap, Jacket	Contact for Quote
YDH-15-216-F	15 L, 300 Vials (2 mL), Incl. Rack, Cryoboxes, Cap	Contact for Quote
YDH-25-216-F	25 L, 500 Vials (2 mL), Incl. Rack, Cryoboxes, Cap	Contact for Quote

- Durable aluminum construction
- Temperature monitoring available
- Liquid or vapor phase option available
- High thermal efficiency
- 5-year vacuum warranty

Haier LN₂ Medium Sized Storage Series



ITEM #	DESCRIPTION	PRICE
YDS-65-216	65 L, 2,400 Vials (2 mL), Incl. Lockable Cap	Contact for Quote
YDS-95-216	95 L, 3,000 Vials (2 mL), Incl. Lockable Cap	Contact for Quote
YDS-115-216	115 L, 3,600 Vials (2 mL), Incl. Lockable Cap	Contact for Quote
YDS-145-216	145 L, 4,800 Vials (2 mL), Incl. Lockable Cap	Contact for Quote
YDS-175-216	175 L, 6,000 Vials (2 mL), Incl. Lockable Cap	Contact for Quote

- LN₂ splash proof
- Data logging
- Hot gas bypass
- 5-year vacuum warranty



Haier Smart Type LN₂ Biobank Series

ITEM #	DESCRIPTION	PRICE
YDD-350-326	350 L, 13,000 Vials, w/ Vacuum Insulated Hose	Contact for Quote
YDD-550-445	587 L, 27,000 Vials, w/ Vacuum Insulated Hose	Contact for Quote
YDD-750-445	783 L, 37,800 Vials, w/ Vacuum Insulated Hose	Contact for Quote
YDD-850-465	890 L, 42,900 Vials, w/ Vacuum Insulated Hose	Contact for Quote
YDD-1300-635	1340 L, 58,500 Vials, w/ Vacuum Insulated Hose	Contact for Quote
YDD-1800-635	1880 L, 94,875 Vials, w/ Vacuum Insulated Hose	Contact for Quote

Don't see what you need?
Contact us for other cryo storage solutions!

WATER SYSTEMS



Genie PURIST ultrapure water system produces Type I ultrapure water from RO, distilled or deionized water.
Ultrapure water for critical laboratory applications.

Rephile PURIST® UV Ultrapure

ITEM #	MAX VOL/DAY	FEED WATER SOURCE	WATER QUALITY	QTY	PRICE
RS2200QUVPK	2.0 L/min	RO/DI/Distilled	Type 1 (18.2MΩ)	1/ea	Contact for Quote

Set Includes Dispenser

NuZar Q water system set is a compact water purification system, producing ultrapure water directly from pure water (EDI, RO, distilled, deionized water).

Rephile NuZar Q



ITEM #	MAX VOL/DAY	FEED WATER SOURCE	WATER QUALITY	QTY	PRICE
RNOQ00000K	Up to 2 L/min	EDI, RO, Distilled, DI	Type 1 (18.2MΩ)	1/ea	Contact for Quote

NuZar U24 water system set is a two-in-one integrated water purification system, producing ultrapure water and reverse osmosis water directly from tap water.

NEW!

Rephile NuZar U24



ITEM #	MAX VOL/DAY	FEED WATER SOURCE	WATER QUALITY	QTY	PRICE
RNOU020T1K	Up to 2 L/min	Tap Water	Ultrapure/RO Type 1 (18.2MΩ)	1/ea	Contact for Quote

Rephile Genie PURIST®

ITEM #	MAX VOL/DAY	FEED WATER SOURCE	WATER QUALITY	QTY	PRICE
RG0S00000K	Up to 2 L/min	EDI, RO, Distilled, DI	Type 1 Ultrapure (18.2MΩ)	1/ea	Contact for Quote



Visit our website for more water systems details!

**Not sure what you need?
Contact us, and our Water Systems Specialist will assist you!**

BALANCES

- External calibration
- Ideal for bulk chemicals, liquids and animals
- Easy in-lab calibration feature
- RS232 port for data output to a computer or printer
- Readability from 0.001 g to 0.01 g
- 2-year warranty

Accuris™ Precision Balances

ITEM #	DESCRIPTION	QTY	PRICE
W3200-120	120 g Max Cap., 0.001 g Readability	1/ea	\$1,138.52
W3200-320	320 g Max Cap., 0.001 g Readability	1/ea	\$1,220.23
W3200-500	500 g Max Cap., 0.001 g Readability	1/ea	\$1,412.90
W3200-1200	1,200 g Max Cap., 0.01 g Readability	1/ea	\$1,144.34
W3200-3200	3,200 g Max Cap., 0.01 g Readability	1/ea	\$1,220.23
W3200-5K	5,000 g Max Cap., 0.01 g Readability	1/ea	\$1,681.47



- Readability from 0.1 to 0.01 g
- Ideal for liquids, animals, general weighing and transportability
- Easy in-lab calibration feature
- RS232 port for data output to a computer or printer
- 2-year warranty



Accuris™ Compact Balances

ITEM #	DESCRIPTION	QTY	PRICE
W3300-120	120 g Max Cap., 0.01 g Readability	1/ea	\$330.49
W3300-300	300 g Max Cap., 0.01 g Readability	1/ea	\$385.50
W3300-500	500 g Max Cap., 0.01 g Readability	1/ea	\$385.50
W3300-1200	1,200 g Max Cap., 0.1 g Readability	1/ea	\$428.36
W3300-5K	3,200 g Max Cap., 0.1 g Readability	1/ea	\$428.36
W3300-10K	5,000 g Max Cap., 0.1 g Readability	1/ea	\$575.26

- Ideal for reagents and basic analytical balance needs
- Multiple mode options (weighting, count mode, percent)
- RS232 port for data output to a computer or printer
- 2-year warranty

Accuris™ Analytical Balances

ITEM #	DESCRIPTION	QTY	PRICE
W3100-120	120 g Max Cap., 0.0001 g Readability, External Calibration	1/ea	\$2,164.45
W3100-210	220 g Max Cap., 0.0001 g Readability, External Calibration	1/ea	\$2,513.19



MICROSCOPES

Research Inverted Live Cell Microscopes



i600 Series

Advanced Optical Performance

NIS60 optical design for Nomarski DIC, phase, Fluorescent

Manual Cell Culture

Users control the microscope for cell culture imaging and simple cell observation\documentation

Manual Fluorescence

User controlled manual components for up to 6 fluorescence wavelengths

Automated and Fast

Selectable High Speed motorized components for fast and software coordinated image acquisition

Multi-Wavelength and Time-Lapse

Integration of Multi Line-High Speed LED sources and Motorized stages

Long Term Experiments

Optional Incubation systems for Temperature and CO₂ control

Multiple Cameras

Acquisition of Multi-Parameter experiments via multiple camera ports

Advanced Imaging

Integration of Confocals, TIRF illuminators or Photo Activation lasers

Software Control

Seamless integration in to 3i SlideBook for acquisition, deconvolution and analysis

Inverted Microscopes

ITEM #	DESCRIPTION	QTY	PRICE
I410PH2XY	4x, 10x, 20x Apodized Plan Phase Objectives: Mech	1/ea	\$5,495.00
I410PH3XY	4x, 10x, 20x, 40x Apodized Plan Phase Objectives	1/ea	\$6,150.00



Stereo Microscopes

ITEM #	DESCRIPTION	QTY	PRICE
Z645	0.67x to 4.5x Continuous Zoom, Binocular, Platform Stand	1/ea	\$1,795.00
Z730	0.7x to 3.5x Continuous Zoom, Binocular, Platform Stand	1/ea	\$1,471.00
Z850	0.8x to 5.0x Continuous Zoom, Binocular, Platform Stand	1/ea	\$2,205.00

Contact your trusted MIDSCI™ Lab Consultant for more information, other options and customization!

GENERAL LAB SUPPLIES



Aurelia™ Amazing®

ITEM #	SIZE	QTY	PRICE
92885	X-Small	300/bx	\$25.05
92886	Small	300/bx	\$25.05
92887	Medium	300/bx	\$25.05
92888	Large	300/bx	\$25.05
92889	X-Large	270/bx	\$25.05

Non-Sterile, Powder-Free, Ambidextrous, Finger Textured, Transblue Color, 2 mil Palm Thickness, 100% Latex-Free



Aurelia™ Quest®

ITEM #	SIZE	QTY	PRICE
92895	X-Small	100/bx	\$8.57
92896	Small	100/bx	\$8.57
92897	Medium	100/bx	\$8.57
92898	Large	100/bx	\$8.57
92899	X-Large	100/bx	\$8.57

Non-Sterile, Powder-Free, Ambidextrous, Finger Textured, White Color, 3.5 mil Palm Thickness, 100% Latex-Free



Aurelia™ Blush®

ITEM #	SIZE	QTY	PRICE
78885	X-Small	200/bx	\$17.14
78886	Small	200/bx	\$17.14
78887	Medium	200/bx	\$17.14
78888	Large	200/bx	\$17.14
78889	X-Large	200/bx	\$17.14

Non-Sterile, Powder-Free, Ambidextrous, Finger Textured, Pink Color, 2.5 mil Palm Thickness, 100% Latex-Free



Aurelia™ Protégé®

ITEM #	SIZE	QTY	PRICE
93995	X-Small	100/bx	\$8.57
93996	Small	100/bx	\$8.57
93997	Medium	100/bx	\$8.57
93998	Large	100/bx	\$8.57
93999	X-Large	100/bx	\$8.57

Non-Sterile, Powder-Free, Ambidextrous, Finger Textured, Sky Blue Color, 4.0 mil Palm Thickness, 100% Latex-Free



Aurelia™ Perform®

ITEM #	SIZE	QTY	PRICE
92395	X-Small	200/bx	\$17.14
92396	Small	200/bx	\$17.14
92397	Medium	200/bx	\$17.14
92398	Large	200/bx	\$17.14
92399	X-Large	200/bx	\$17.14

Non-Sterile, Powder-Free, Ambidextrous, Finger Textured, Teal Color, 2.5 mil Palm Thickness, 100% Latex-Free

*Palm Thickness.



Aurelia™ Bold®

ITEM #	SIZE	QTY	PRICE
73995	X-Small	100/bx	\$9.01
73996	Small	100/bx	\$9.01
73997	Medium	100/bx	\$9.01
73998	Large	100/bx	\$9.01
73999	X-Large	100/bx	\$9.01
73994	XX-Large	100/bx	\$9.01

Non-Sterile, Powder-Free, Ambidextrous, Finger Textured, Black Color, 5 mil Palm Thickness, 100% Latex-Free

GENERAL LAB SUPPLIES

- Safely dispose of tips, tubes, gloves, etc.
- Easy push open/close spring lid
- Secures and covers hazardous waste
- Compliant with OSHA standard 29 CFR 1910.1030

BenchBin™ Benchtop Biohazard Bins

ITEM #	DESCRIPTION	QTY	PRICE
A8401B	Benchtop Biohazard Bin, 1.5L Capacity	1/ea	\$40.91
New! A8403B	Benchtop Biohazard Bin, 3L Capacity	1/ea	\$48.83
A9004R	1.5 L Red Biohazard Bags w/ Black Printed markings, 8.5 x 11 in.	100 bgs/pk	\$31.87
A9006R	3 L Red Biohazard Bags w/ Black Printed markings, 10 x 12.5 in.	100 bgs/pk	\$40.09



Bowtie™ Biohazard Bins

ITEM #	DESCRIPTION	QTY	PRICE
A8000B	8 gallon capacity (30L) Bin, w/ Hands-Free Foot Pedal, 10 x 14 x 19 in. [PP], Attached Lid	1/ea	\$80.73
New! A8006B	13 gallon capacity (50L) Bin, w/ Hands-Free Foot Pedal, 16 x 16 x 24 in. [PP], Attached Lid	1/ea	\$109.00
A8002B	26 gallon capacity (100L) XL Bin, w/ Foot Pedal Operated Lid, 19 x 19 x 33 in.	1/ea	\$194.29
A8001R	27 x 35 in. Biohazard, Autoclave Bag w/ Drawstring, [PE], Marking Area, Sterilization Indicator	100/pk	\$163.89
A8003R	33 x 41 in. Biohazard, Autoclave Bag w/ Drawstring, [PE], Marking Area, Sterilization Indicator	100/pk	\$192.75

Laboratory Chairs



Laboratory Chairs

Our lab chairs feature the latest advancements in laboratory ergonomics, such as contoured waterfall cushions and 6-way armrests.



Lab Stools With Backs

Each lab stool features fully adjustable ergonomics for optimal comfort.

Laboratory Stools

Lab stools feature a 360° all-around pneumatic height adjustment for easy adjustability, making them perfect for many environments.

Speciality Seating

Laboratory speciality seating provides the same comfortable ergonomics as our conventional laboratory seating product line.



We have so many options, it's hard to list them all!
Please contact us, we'd be happy to assist you in selecting the perfect chair for your needs.



MIDSCI Freebies with Purchase



Visit midsci.com
for more items!

ITEM #	DESCRIPTION	QTY	QUALIFYING PURCHASE
MS-SHARK-R-F	MIDSCI 23 in. Inflatable Shark, Red	1/ea	\$300
MS-SHARK-B-F	MIDSCI 23 in. Inflatable Shark, Blue	1/ea	\$300
MS-LABLNYRDS-F	MIDSCI Lab Scientist Dog 3/4 in. Lanyard	1/ea	\$300
MS-RBRSCIDUCK-F	MIDSCI Rubber Scientist Duck	1/ea	\$600
MS-SPOONSCALE-F	MIDSCI Electronic Measuring Spoon, Blue	1/ea	\$600
MSTIME-F	MIDSCI Magnetic Timer, 0 to 99 Minutes	1/ea	\$900
MS-WATRBG-BL-F	MIDSCI Cocona Water Bottle Bag*	1/ea	\$900
MS-ADVBAG-R-F	MIDSCI Adventure Duffle Bag	1/ea	\$1,200
MS-TOOL26PC-S-F	MIDSCI Precision 26Pc Tool Kit, Silver	1/ea	\$1,200
MS-TUMBLEBLUE-F	MIDSCI 40 oz. Insulated Tumbler, Blue	1/ea	\$1,500
MS-TUMBLERED-F	MIDSCI 40 oz. Insulated Tumbler, Red	1/ea	\$1,500
LCTD-S-F	Tie-Dye Lab Coat, S, 42 in. Length	1/ea	\$3,000
LCTD-M-F	Tie-Dye Lab Coat, M, 42 in. Length	1/ea	\$3,000
LCTD-L-F	Tie-Dye Lab Coat, L, 42 in. Length	1/ea	\$3,000
LCTD-XL-F	Tie-Dye Lab Coat, XL, 42 in. Length	1/ea	\$3,000
LCTD-XXL-F	Tie-Dye Lab Coat, XXL, 42 in. Length	1/ea	\$3,000

Disclaimer: Free gifts are available with qualifying purchases only, no exceptions.

Limit of 3 free gifts per order, based on the total amount spent. Offer cannot be combined with any other promotions or discounts. All terms and conditions apply.

* Bag only, does not include cup or phone

Images are not to scale



Family owned and trusted since 1983



Stocked products ship in 24 hours



Responses within 24 hours



REAL people answering your questions and placing your orders

We're here to help.

Place your order with us today!

ORDER TODAY! midsci.com • custserv@midsci.com • 800.227.9997 • f: 636.225.9998