

MIDSCI™



C. Elegans



Microbiology



Cold Storage



Cell Culture

Drosophila



Large Enough to Deliver, Small Enough to Care

ORDER TODAY! midsci.com • custserv@midsci.com • 800.227.9997 • f: 636.225.9998

TABLE OF CONTENTS

3 - 4	•••••	C. Elegans Essentials
5 - 15	•••••	Cell Culture Essentials
16 - 19	•••••	Cold Storage Essentials
20 - 22	•••••	Drosophila Essentials
23 - 28	•••••	Microbiology Essentials
29 - 36	•••••	Lab Essentials

Greetings to our valued customers,

Over the last 40 years MIDSCI™ has established a reputation of being a “go-to” for scientific supplies and equipment. We realize people are our most important asset. We make certain to take care of our MIDSCI™ family, so we can ensure we all work together to take care of you.

We believe in providing the best value to our customers by offering quality products and exceptional service. Our knowledge and expertise have allowed us to earn the trust of our customers daily, and we take pride in maintaining that trust.

In addition, we still have real people answering our phones and providing you with real answers. We not only value our relationship with you but take the responsibility of being the supplier you can count on.

To our long-time customers, we would like to thank you for your incredible support and your belief and trust in us.

To our new customers, we look forward to earning your trust and business well into the future.

Sincerely,
Mike Degenhart
President



If you can't find it here, give us a call at **800-227-9997** or email us at tech@midsci.com



MIDSCI assumes no liability for printing errors, mistakes, or product changes. All prices subject to change without notice.

***All promotions valid through 09/30/24, unless otherwise noted. Promos cannot be combined with any other promotions or discounts. All terms and conditions apply. Visit our website for more information: <https://midsci.com/terms>*

This catalog's cover page features a Visual Artist Concept; created by using AI-generated images, included for illustration and decorative purposes only. Accuracy, quality, and appropriateness may vary.

V.07.2024

C. ELEGANS ESSENTIALS

PR1MA™ Nematode Growth Medium

ITEM #	DESCRIPTION	QTY	PRICE
KCN81800-500	NGM Pre-Mixed Growth Media, 500 g	1/ea	\$86.25
KCN81800-1000	NGM Pre-Mixed Growth Media, 1 Kg	1/ea	\$159.23
KCN81800-5000	NGM Pre-Mixed Growth Media, 5 Kg	1/ea	\$749.42



PR1MA™ NGM Media Component Essentials

ITEM #	DESCRIPTION	QTY	PRICE
KCA20030-500	Agar, Bacteriological Grade, Powder, 500 g	1/ea	\$195.29
KCA20030-1000	Agar, Bacteriological Grade, Powder, 1 Kg	1/ea	\$323.08
KCA20030-5000	Agar, Bacteriological Grade, Powder, 5 Kg	1/ea	\$1,053.46
KCA20030-10000	Agar, Bacteriological Grade, Powder, 10 Kg	1/ea	\$1,606.00
KCS23020-500	Sodium Chloride, 500 g	1/ea	\$21.35
KCS23020-1000	Sodium Chloride, 1 Kg	1/ea	\$38.65
KCS23020-5000	Sodium Chloride, 5 Kg	1/ea	\$102.12
KCP20240-500	Peptone Powder, 500 g	1/ea	\$67.21
KCP20240-1000	Peptone Powder, 1 Kg	1/ea	\$124.33
KCP20240-5000	Peptone Powder, 5 Kg	1/ea	\$543.17

Standard Petri Dishes

ITEM #	DESCRIPTION	QTY	PRICE
706011	NEST® Dish 35 mm x 12 mm, 1 Compartment, Sterile	20/pk, 25 pks/cs	\$115.32
667665	PR1MA™ 60 mm x 15 mm, 1 Compartment, Sterile	25/pk, 20 pks/cs	\$133.19



C. ELEGANS EQUIPMENT ESSENTIALS

- Digital temperature control
- Personal-sized, 13 x 14.5 inch footprint (20 L)
- Accepts bottles and flasks up to 2 L
- Internal outlet, for installation of shaker, rockers, etc.

PRIMA[™] Mini Digital Incubators

ITEM #	DESCRIPTION	QTY	PRICE
IN20-H	Heat Only (Ambient +5 to 60°C), Incl. 2 Shelves	1/ea	\$596.55
IN20-HC	Mini Digital Incubator - Heat and Cool, Incl. 2 Shelves	1/ea	\$907.43
H2200-SH	Extra Shelf, Plastic	1/ea	\$47.69



- Accurate temperature control from 0° to 60°C
- Powerful compressor maintains 0°C
- Stainless-steel internal chamber
- Internal power outlet for shaker, mixers, etc.
- Sized for floor standing use or on benchtops less than 32" tall

MyTemp[™] 65 L Incubators

ITEM #	DESCRIPTION	QTY	PRICE
H2265-HC	Digital Refrigerated Incubator, 65 L	1/ea	\$3,728.40

Contact us for a quote!

Stereo Microscopes

ITEM #	DESCRIPTION	QTY	PRICE
Z850	0.8X-5.0X Continuous Zoom, Binocular, Course/Fine Focus	1/ea	\$2,060.00
Z850LED	0.8X-5.0X Continuous Zoom, Binocular w/ Epi/Dia LED, Course/Fine Focus	1/ea	\$2,286.00
Z850T	0.8X-5.0X Continuous Zoom, Trinocular, Course/Fine Focus	1/ea	\$2,281.00
Z850TLED	0.8X-5.0X Continuous Zoom, Trinocular w/ Epi/Dia LED, Course/Fine Focus	1/ea	\$2,507.00



ADHERENT CELL CULTURE

- Wider neck allows for better access for 100% cell retrieval when scraping
- Short access to cell floor reduces chance of contamination



TPP® Tissue Culture Flasks

ITEM #	DESCRIPTION	QTY	PRICE
TP90025	25 cm ² (92 x 51 x 29 mm) TC Flask [PS], Solid Vent Cap [PE], Sterile	10/pk, 36 pks/cs	\$794.00
TP90026	25 cm ² (92 x 51 x 29 mm) TC Flask [PS], Filter Cap [PE], Sterile	10/pk, 36 pks/cs	\$867.03
TP90075	75 cm ² (155 x 87 x 40 mm) TC Flask [PS], Solid Vent Cap [PE], Sterile	5/pk, 20 pks/cs	\$451.70
TP90076	75 cm ² (155 x 87 x 40 mm) TC Flask [PS], Filter Cap [PE], Sterile	5/pk, 20 pks/cs	\$485.57
TP90150	150 cm ² (210 x 122 x 50 mm) TC Flask [PS], Solid Vent Cap [PE], Sterile	3/pk, 13 pks/cs	\$364.88
TP90151	150 cm ² (210 x 122 x 50 mm) TC Flask [PS], Filter Cap [PE], Sterile	3/pk, 13 pks/cs	\$388.30
TP90300	300 cm ² (275 x 170 x 50 mm) TC Flask [PS], Solid Vent Cap [PE], Sterile	3/pk, 6 pks/cs	\$389.52
TP90301	300 cm ² (275 x 170 x 50 mm) TC Flask [PS], Filter Cap [PE], Sterile	3/pk, 6 pks/cs	\$408.34

Filter caps come with a 0.22 µm PTFE Membrane

- Access of a dish & security of a flask
- Ideal for long incubations, colony picking, organ growth, & so much more

TPP® Re-Closable Flasks

ITEM #	DESCRIPTION	QTY	PRICE
TP90552	150 cm ² (210 x 122 x 50 mm) TC Flask [PS], Filter Cap [PE], Sterile	3/pk, 6 pks/cs	\$536.00
TP90652	115 cm ² (210 x 122 x 50 mm) TC Flask [PS] w/ Barrier, Filter Cap [PE], Sterile	3/pk, 6 pks/cs	\$537.90

Filter caps come with a 0.22 µm PTFE Membrane



- Access of a dish & security of a flask




TPP® Peel-Off Flasks

ITEM #	DESCRIPTION	QTY	PRICE
TP90028	25 cm ² (92 x 51 x 31 mm) TC Flask [PS], Filter Cap [PE], Sterile	5/pk, 14 pks/cs	\$939.70
TP90078	75 cm ² (155 x 87 x 45 mm) TC Flask [PS], Filter Cap [PE], Sterile	2/pk, 16 pks/cs	\$489.69
TP90653	115 cm ² (210 x 122 x 50 mm) TC Flask [PS] w/ Barrier, Filter Cap [PE], Sterile	3/pk, 6 pks/cs	\$452.88
TP90153	150 cm ² (210 x 122 x 50 mm) TC Flask [PS], Filter Cap [PE], Sterile	3/pk, 6 pks/cs	\$445.28
TP90303	300 cm ² (275 x 170 x 55 mm) TC Flask [PS], Filter Cap [PE], Sterile	2/pk, 4 pks/cs	\$200.51

ADHERENT CELL CULTURE

You can always trust MIDSCI™ for quality TPP® products that help you save money where it matters most, cost per experiment.

- TPP® is Swiss made with excellent flatness and clarity
- Flatness = Less clumping and even growth
- Clarity = Great for microscopy
- Vents = Rapid temperature equilibration for better cell consistency and health
- Stacked dishes, flasks, and plates interlock and create vents between levels
- Interlock makes handling stacks much less prone to dropping



TPP® Tissue Culture Plates

- Minimized edge effect in 96-well plates
- Easy-to-read yellow marking spot and black alphanumeric

ITEM #	DESCRIPTION	QTY	PRICE
TP92006	6-Well TC Plate [PS], Flat Bottom, Sterile	1/pk, 126 pks/cs	\$711.31
TP92012	12-Well TC Plate [PS], Flat Bottom, Sterile	1/pk, 126 pks/cs	\$711.96
TP92024	24-Well TC Plate [PS], Flat Bottom, Sterile	1/pk, 126 pks/cs	\$710.69
TP92048	48-Well TC Plate [PS], Flat Bottom, Sterile	1/pk, 126 pks/cs	\$711.32
TP92096	96-Well TC Plate [PS], Flat Bottom, Sterile	1/pk, 162 pks/cs	\$877.09
TP92097	96-Well TC Plate [PS], Round Bottom, Sterile	1/pk, 162 pks/cs	\$872.47

Also available in multipacks!

Contact your trusted MIDSCI™ Lab Consultant for more information

- Easy to read painted yellow marking spot
- Ridged grip around base greatly improves handling. No more dropped dishes



TPP® Tissue Culture Dishes

ITEM #	DESCRIPTION	QTY	PRICE
TP93040	40 x 11 mm TC Dish [PS], Sterile	20/pk, 45 pks/cs	\$701.57
TP93060	60 x 16 mm TC Dish [PS], Sterile	14/pk, 60 pks/cs	\$725.52
TP93100	96 x 21 mm TC Dish [PS], Sterile	10/pk, 24 pks/cs	\$393.88
TP93150	146 x 21 mm TC Dish [PS], Sterile	5/pk, 20 pks/cs	\$466.96

ADHERENT CELL CULTURE



- Ideal for starting cultures & primary cells
- Culture & spin down - same tube
- Reduces plastics, cost, & lose fewer cells during transfer

TPP® Tiny Flasks

ITEM #	DESCRIPTION	QTY	PRICE
TP91106	20 cm ² (16 x 121 mm) Mini Roller Bottle [PS], Vent Cap [PE], Sterile	20/pk, 40 pks/cs	\$1,031.94
TP91243	10 cm ² (30 x 105 mm) Tube Flask [PP], Filter Cap [PE] 0.22 µm PTFE Membrane, Sterile	4/pk, 54 pks/cs	\$503.44

- Culture & image on a (COP) slide
- Suitable for: confocal fluorescence microscopy, individual coatings, staining and especially for immunofluorescence measurements
- Resistant to solvents (e.g. acetone, ethanol, xylene)
- Removable medium chamber with a "click" and without the need of a removing tool

TPP® Clipmax Chamber Slide System

ITEM #	DESCRIPTION	QTY	PRICE
TP70010	10 cm ² (93 x 33 x 32 mm) Tube Flask [PS], Filter Cap [PE] 0.22 µm PTFE Membrane, Sterile	1/pk, 5 pks/cs	\$108.63



Need other chamber slide options?
Just ask!



- Scraper heads made from a flexible polyethylene material with sharp edges
- Small raised knobs on handle for secure grip
- Single-wrapped in blister packaging

TPP® Scraper & Spatula

ITEM #	DESCRIPTION	QTY	PRICE
TP99002	24 cm Scraper "S" [PE], Rotating, Sterile	1/pk, 150 pks/cs	\$625.53
TP99003	30 cm Scraper "M" [PE], Rotating, Sterile	1/pk, 150 pks/cs	\$707.45
TP99004	38 cm Scraper "L" [PE], Rotating, Sterile	1/pk, 100 pks/cs	\$584.72
TP99010	20 cm Spatula [PP], Fixed, Sterile	1/pk, 100 pks/cs	\$551.53

ADHERENT CELL CULTURE



- Made of high clarity, 100% virgin polystyrene
- Sterilized by E-beam, SAL=10⁻⁶
- Non-Pyrogenic, DNase/RNase free
- Individually packaged in double-layer sterile bags

NEST® BioFactory™

ITEM #	DESCRIPTION	QTY	PRICE
771001	1 Layer, 2 Wide Plug Seal Caps, Plus 16 Sterile Wide Vent Caps	8/cs	\$552.25
771101	2 Layers, 2 Wide Plug Seal Caps, Plus 16 Sterile Wide Vent Caps	8/cs	\$851.00
771204	5 Layers, 2 Wide Plug Seal Caps, Plus 8 Sterile Wide Vent Caps	4/cs	\$622.41
771302	10 Layers, 2 Wide Plug Seal Caps, Plus 12 Sterile Wide Vent Caps	6/cs	\$1,051.59
771403	40 Layers, 2 Wide Plug Seal Caps, Plus 4 Sterile Wide Vent Caps	2/cs	\$1,001.51
772001	1 Layer, 2 Narrow Plug Seal Caps, Plus 16 Sterile Wide Vent Caps	8/cs	\$552.25
772101	2 Layers, 2 Narrow Plug Seal Caps, Plus 16 Sterile Wide Vent Caps	8/cs	\$851.00

Vent caps packaged separately

A BioFactory™ with pre-installed set of transfer tubes and filters, which can directly transfer culture media and cells to or out of BioFactory™ in a sterile environment by connecting tubes, without change caps or connect tubes, greatly reduce the risk of exogenous pollution in the drug development and production process, simultaneously omit the operations of tube design, assemble and sterilization, improve the production efficiency.



NEST® Biofactory™ 10 Chamber with Pre-installed with Transfer Cap System

ITEM #	DESCRIPTION	QTY	PRICE
774301	0.22 µm, 62 mm Diameter, 80 mm Hose (1/4" ID, 3/8" OD)	2/cs	\$336.00

Solid overcaps can be used directly on the adaptor caps to block dust and bacteria.



NEST® BioFactory™ with 2 Solid Overcaps

ITEM #	DESCRIPTION	QTY	PRICE
771322	10 Chamber, Plus 12 Sterile Vented Overcaps	6/cs	\$1,001.51
771422	40 Chamber, Plus 4 Sterile Vented Overcaps	2/cs	\$1,001.51

Overcaps packaged separately

SUSPENSION CELL CULTURE

- Ideal for suspension cultures
- Scale up or scale down ability to do large in parallel studies
- Better viability & proliferation than flasks or bottles
- Better gas exchange control & cyclonic stirring than flasks or bottles

TPP® TubeSpin® Bioreactor

ITEM #	DESCRIPTION	QTY	PRICE
TP87015	15 mL Conical [PP], Sterile	40/pk, 20 pks/cs	\$1,111.48
TP87017	15 mL Round [PP], Sterile	40/pk, 20 pks/cs	\$1,665.50
TP87017-200	15 mL Round [PP], Sterile	40/pk, 5 pks/cs	\$414.91
TP87017-400	15 mL Round [PP], Sterile	40/pk, 10 pks/cs	\$834.22
TP87050	50 mL Conical [PP], Sterile	20/pk, 9 pks/cs	\$484.26
TP87050-60	50 mL Conical [PP], Sterile	20/pk, 3 pks/cs	\$177.91
TP87056	50 mL Round [PP], Sterile	20/pk, 9 pks/cs	\$755.42
TP87450	450 mL Conical [PP], Sterile	1/pk, 32 pks/cs	\$643.86
TP87600	600 mL Conical [PP], Sterile	1/pk, 26 pks/cs	\$1,015.84
TP87600-13	600 mL Conical [PP], Sterile	1/pk, 13 pks/cs	\$506.13
TP87600-3	600 mL Conical [PP], Sterile	1/pk, 3 pks/cs	\$117.33



- Ideal for high throughput & robotics
- Re-closeable septum allows repeated access while keeping sterile
- Compatible for injection needles up to 2.1 mm diameter
- Ideal for PRP kits

TPP® TubeSpin® Bioreactor w/ Septum

ITEM #	DESCRIPTION	QTY	PRICE
TP86050	50 mL Conical [PP] w/ Septum, Sterile	20/pk, 9 pks/cs	\$2,785.71



Barcoded for automated tracking



PR1MA™ Erlenmeyer Flasks

ITEM #	DESCRIPTION	QTY	PRICE
667801	125 mL, Clear, Vent Cap, Plain Bottom, PETG, Sterile	1/pk, 12 pks/cs	\$113.00
667805	250 mL, Clear, Vent Cap, Plain Bottom, PETG, Sterile	1/pk, 12 pks/cs	\$128.00
667809	500 mL, Clear, Vent Cap, Plain Bottom, PETG, Sterile	1/pk, 12 pks/cs	\$163.00
667813	1,000 mL, Clear, Vent Cap, Plain Bottom, PETG, Sterile	1/pk, 6 pks/cs	\$124.00

SUSPENSION CELL CULTURE

- Available in 6-, 12-, 24-, 48-, and 96-well format
- Polystyrene, NOT tissue culture treated (suspension culture)
- Sterile (gamma irradiated)

PR1MA™ Multi-Well Plates

ITEM #	DESCRIPTION	QTY	PRICE
667506	6-Well, [PS] Non-Treated, w/ Lid, Idiv. Wrapped, Sterile	1/pk, 100 pks/cs	\$205.05
667512	12-Well, [PS] Non-Treated, w/ Lid, Idiv. Wrapped, Sterile	1/pk, 100 pks/cs	\$256.00
667524	24-Well, [PS] Non-Treated, w/ Lid, Idiv. Wrapped, Sterile	1/pk, 100 pks/cs	\$236.48
667548	48-Well, [PS] Non-Treated, w/ Lid, Idiv. Wrapped, Sterile	1/pk, 100 pks/cs	\$250.77
667596	96-Well, [PS] Non-Treated, w/ Lid, Idiv. Wrapped, Sterile	1/pk, 100 pks/cs	\$285.00



Suspension Cell Culture Flasks

ITEM #	DESCRIPTION	QTY	PRICE
229500	25 mL, [PS] Non-Treated, HDPE Vent Cap, Sterile	10/bg, 20 bgs/cs	\$271.00
229510	50 mL, [PS] Non-Treated, HDPE Vent Cap, Sterile	10/bg, 20 bgs/cs	\$284.00
229520	250 mL, [PS] Non-Treated, HDPE Vent Cap, Sterile	5/bg, 20 bgs/cs	\$268.00
229530	600 mL, [PS] Non-Treated, HDPE Vent Cap, Sterile	5/bg, 8 bgs/cs	\$241.00
229540	850 mL, [PS] Non-Treated, HDPE Vent Cap, Sterile	3/bg, 6 bgs/cs	\$251.00



- Polycarbonate (PC) body, meets USP Class 6 Standard
- Withstands high temperature, up to 121°C
- No DNA or RNA enzymes, no pyrogens
- Sterilized by E-Beam

NEST High Efficiency Erlenmeyer Flasks

ITEM #	DESCRIPTION	QTY	PRICE
785101	2 L, HDPE Vent Cap, Individ. Wrapped, Sterile	4/pk	\$149.79
786111	3 L, HDPE Vent Cap, Individ. Wrapped, Sterile	5/pk	\$259.19
786511	3 L, Wide Mouth, HDPE Vent Cap, Individ. Wrapped, Sterile	6/pk	\$259.19
787011	5 L, HDPE Vent Cap, Individ. Wrapped, Sterile	7/pk	\$387.61



CELL CULTURE MEDIA ESSENTIALS

- Unique design allows for more of what you want through and less of what you don't
- Less tipping & faster filtration
- 0.22 µm PES membrane is the lowest protein membrane available
- 45 mm neck

TPP® Filtermax Filters

ITEM #	DESCRIPTION	QTY	PRICE
TP99150	150 mL, Filter Top & Reservoir, Sterile	1/pk, 18 pks/cs	\$391.07
TP99157	150 mL, Reservoir, Sterile	1/pk, 24 pks/cs	\$210.34
TP99250	250 mL, Filter Top & Reservoir, Sterile	1/pk, 12 pks/cs	\$295.67
TP99257	250 mL, Reservoir, Sterile	1/pk, 24 pks/cs	\$267.11
TP99500	500 mL, Filter Top & Reservoir, Sterile	1/pk, 10 pks/cs	\$310.97
TP99505	500 mL, Filter Top, Sterile	1/pk, 21 pks/cs	\$490.79
TP99507	500 mL, Reservoir, Sterile	1/pk, 36 pks/cs	\$652.45
TP99950	1 L, Filter Top & Reservoir, Sterile	1/pk, 9 pks/cs	\$547.67
TP99955	1 L, Filter Top, Sterile	1/pk, 12 pks/cs	\$519.83
TP99957	1 L, Reservoir, Sterile	1/pk, 20 pks/cs	\$554.36



- Design allows for faster filtration & less clogging
- Fits the rim of a 50 mL perfectly
- 0.22 µm or 0.45 µm PES membrane is the lowest protein membrane available



TPP® Syringe Filters

ITEM #	DESCRIPTION	QTY	PRICE
TP99722	0.22 µm, Sterile	1/pk, 40 pk/cs	\$205.90
TP99745	0.45 µm, Sterile	1/pk, 40 pk/cs	\$205.90

Harvest Fetal Bovine Serum (FBS)

ITEM #	DESCRIPTION	QTY	PRICE
USFBS	US Origin, Fetal Bovine Serum (FBS)	500 mL	\$555.00
USDAFBS	USDA Certified, Fetal Bovine Serum (FBS)	500 mL	\$524.48
USFBSHI	US Origin, Fetal Bovine Serum (FBS), Heat Inactivated	500 mL	\$560.00
USDAFBSHI	USDA Certified, Fetal Bovine Serum (FBS), Heat Inactivated	500 mL	\$530.00



FILTRATION

- Polypropylene housing with nylon mesh filter
- Extended gripping tab on housing for easy handling
- The elongated lip on the housing extends beyond the edge of the centrifuge tube enabling aseptic handling with forceps



PRIMA™ Economy Cell Strainer

ITEM #	DESCRIPTION	QTY	PRICE
40ICS	40 µm, Blue, Individ. Wrapped, Sterile	1/pk, 50 pks/cs	\$124.84
70ICS	70 µm, White, Individ. Wrapped, Sterile	1/pk, 50 pks/cs	\$124.84
100ICS	100 µm, Yellow, Individ. Wrapped, Sterile	1/pk, 50 pks/cs	\$124.84
SG-PEST	5.5 in. Pestle, Blue, Individ. Wrapped, Sterile	1/pk, 50 pks/cs	\$60.00

- Advanced design - faster flow through
- Superior ventilation - prevents clogging/spills
- Designed for use with 50 mL tubes (adaptor for 15 mL tubes included)



Harvest Cell Strainer

ITEM #	DESCRIPTION	QTY	PRICE
40CELL	40 µm, Blue, w/ Reducing adapter for 15 mL Tube, Sterile	50/pk	\$109.00
70CELL	70 µm, White, w/ Reducing adapter for 15 mL Tube, Sterile	50/pk	\$109.00
100CELL	100 µm, Yellow, w/ Reducing adapter for 15 mL Tube, Sterile	50/pk	\$109.00
CELLADAPT	Reducing Adaptors for 5 mL and 15 mL Tubes, Sterile, Individ. Wrapped	25/pk	\$47.88



Harvest Double Sided Pestle

ITEM #	DESCRIPTION	QTY	PRICE
PEST2	Harvest Double Sided Pestle for Cell Strainers and Tubes, Individ. Wrapped, Sterile	50/pk	\$51.00

CELL CULTURE EQUIPMENT

- Centurion Touchscreen Controller
- Dimmable LED
- Ergonomic Features
- Equipped with energy-efficient DC ECM
- ULPA Filter Isocide™ Antimicrobial Powder Coating



Airstream® NS G4 Class II Type A2 Bio Safety Cabinet- AC2 Models

ITEM #	DESCRIPTION	QTY	PRICE
AC2-3S9-NS-G4	Gen 4 Biosafety Cabinet, 3 ft., 8 in. Sash	1/ea	\$7,800.00
AC2-4S9-NS-G4	Gen 4 Biosafety Cabinet, 4 ft., 8 in. Sash	1/ea	\$8,117.00
AC2-5S9-NS-G4	Gen 4 Biosafety Cabinet, 5 ft., 8 in. Sash	1/ea	\$8,810.00
AC2-6S9-NS-G4	Gen 4 Biosafety Cabinet, 6 ft., 8 in. Sash	1/ea	\$9,748.00
STA-3A0	Support Stand for 3 ft. G4 cabinets	1/ea	\$460.00
STA-4A0	Support Stand for 4 ft. G4 cabinets	1/ea	\$482.00
STA-5A0	Support Stand for 5 ft. G4 cabinets	1/ea	\$526.00
STA-6A0	Support Stand for 6 ft. G4 cabinets	1/ea	\$548.00

- Contamination control with long life ULPA filters and 90°C moist heat decontamination cycle.
- Direct heat & air jacket design with gentle yet effective forced convection.
- Intuitive interfaces and intelligent data and event logger.

CelCulture® CO₂ Incubator

ITEM #	DESCRIPTION	QTY	PRICE
CCL-170B-9-SP	170L CelCulture® CO ₂ Incubator, IR Sensor, CO ₂ Control, ULPA, High Temp Decontamination	1/ea	\$6,363.00



- Smooth rounded internal corners for easy cleaning
- Infra-Red(IR) CO₂ sensor for precise and reliable CO₂ detection
- Stackable design up to 3 units
- 25 mm and 50 mm orbit models available



AEOLUS™ Incubator Shakers

ITEM #	DESCRIPTION	QTY	PRICE
IS-AMA	Incubator Shaker, 7" High-Res Touch Screen, Air Jacket, 25 mm (1.02 in.) Orbit, Dims 970 x 780 x 530 mm	1/ea	\$17,550.00
IS-A5MA	Incubator Shaker, 7" High-Res Touch Screen, Air Jacket, 50 mm (2 in.) Orbit, Dims 970 x 780 x 530 mm	1/ea	\$19,500.00
IS-ACMA	CO ₂ Incubator Shaker, 7" High-Res Touch Screen, Air Jacket, 25 mm (1.02 in.) Orbit, 970 x 780 x 530 mm	1/ea	\$22,725.00
IS-A5CMA	CO ₂ Incubator Shaker, 7" High-Res Touch Screen, Air Jacket, 50 mm (2 in.) Orbit, 970 x 780 x 530 mm	1/ea	\$23,730.00

CELL CULTURE EQUIPMENT

Inverted Microscopes

ITEM #	DESCRIPTION	QTY	PRICE
I410PH2XY	4x, 10x, 20x Apodized Plan Phase Objectives: Mech	1/ea	\$5,475.00
I410PH3XY	4x, 10x, 20x, 40x Apodized Plan Phase Objectives	1/ea	\$6,101.00



**Need florescence?
Contact your trusted MIDSCI™ Lab
Consultant for more information!**

Research Inverted Live Cell Microscopes



i600 Series
Research Inverted Live
Cell Microscope

Advanced Optical Performance

NIS60 optical design for Nomarski DIC, phase, Fluorescent

Manual Cell Culture

Users control the microscope for cell culture imaging and simple cell observation\documentation

Manual Fluorescence

User controlled manual components for up to 6 fluorescence wavelengths

Automated and Fast

Selectable High Speed motorized components for fast and software coordinated image acquisition

Multi-Wavelength and Time-Lapse

Integration of Multi Line-High Speed LED sources and Motorized stages

Long Term Experiments

Optional Incubation systems for Temperature and CO₂ control

Multiple Cameras

Acquisition of Multi-Parameter experiments via multiple camera ports

Advanced Imaging

Integration of Confocals, TIRF illuminators or Photo Activation lasers

Software Control

Seamless integration in to 3i SlideBook for acquisition, deconvolution and analysis

**Need a Live Cell Imaging Scope?
Contact us!**

LIVE CELL IMAGING



- Real-time live cell monitoring of bottom layer for multi-layer flasks up to 10 levels.
- Checks cell quality and detects confluency rate while culturing
- Improved reproducibility for mass cell culture and production
- No need to take the cells out of the incubator thus eradicating contamination worries
- Multi-position time-lapse imaging & recording

Curiosis® Celloger Stack

ITEM #	DESCRIPTION	QTY	PRICE
CRCLG-SB02	Real-Time Cell Monitoring System	1/ea	\$40,000.00

- Compact size that easily fits into standard CO₂ incubators
- Fully motorized multi-position imaging up to 96 well plate
- Compatible with various cell and tissue culture vessel types
- Multiple focal planes can be captured through Z-stack imaging



Curiosis® Celloger Mini Plus

ITEM #	DESCRIPTION	QTY	PRICE
CRCLG-MPB04	Live Cell Imaging (Bright Field, 4X)	1/ea	\$25,000.00
CRCLG-MPB10	Live Cell Imaging (Bright Field, 10X)	1/ea	\$26,500.00
CRCLG-MPBG04	Live Cell Imaging (Bright Field+Green Fluore, 4X)	1/ea	\$33,500.00
CRCLG-MPBG10	Live Cell Imaging (Bright Field+Green Fluore, 10X)	1/ea	\$35,000.00
CRCLG-MPBR04	Live Cell Imaging (Bright Field+Red Fluore, 4X)	1/ea	\$33,500.00
CRCLG-MPBR10	Live Cell Imaging (Bright Field+Red Fluore, 10X)	1/ea	\$35,000.00

- One color fluorescence (GFP or RFP) and brightfield imaging technology
- Compact size that easily fits into standard CO₂ incubators
- High quality system yet cost-effective
- Compatible with various cell and tissue culture vessel types
- Intuitive UI/UX and easy to acquire confluency data



Curiosis® Celloger Nano

ITEM #	DESCRIPTION	QTY	PRICE
CRCLG-NB04	Live Cell Imaging (Bright Field, 4X)	1/ea	\$5,250.00
CRCLG-NB10	Live Cell Imaging (Bright Field, 10X)	1/ea	\$7,500.00
CRCLG-NBG04	Live Cell Imaging (Bright Field+Green Fluore, 4X)	1/ea	\$8,100.00
CRCLG-NBG10	Live Cell Imaging (Bright Field+Green Fluore, 10X)	1/ea	\$8,850.00
CRCLG-NBR04	Live Cell Imaging (Bright Field+Red Fluore, 4X)	1/ea	\$8,100.00
CRCLG-NBR10	Live Cell Imaging (Bright Field+Red Fluore, 10X)	1/ea	\$8,850.00

COLD STORAGE ESSENTIALS



- Unique cap & vial design eliminates silicone gasket pops
- Self standing
- Externally threaded for reduced contamination
- Color coded cap inserts available

TPP® CryoTubes

ITEM #	DESCRIPTION	QTY	PRICE
TP89012-100	1.2 mL [PP], External Thread, Self Standing, Natural, Sterile	100/pk	\$82.49
TP89020-100	2 mL [PP], External Thread, Self Standing, Natural, Sterile	100/pk	\$84.77
TP89050-100	4.5 mL [PP], External Thread, Self Standing, Natural, Sterile	100/pk	\$105.14
TP99020	Cap Inserts, Six Assorted Colors	600/cs	\$91.21

- Manufactured from high-clarity, ultra-pure polypropylene (meets FDA 21 CFR 177.1520 & USP Class VI)
- Tubes comply with the In Vitro Devices Directive (IVD)
- Contrasting black 100 µL graduations are bold and easy to read
- Large white writing area on all tubes
- Resealable Packaging - Tubes are packaged in resealable stand-up pouches that minimize contamination and provide easy access

New!

PR1MA™ Cryogenic Tubes



ITEM #	DESCRIPTION	QTY	PRICE	PROMO PRICE
PR-CVET10-SB	1.0 mL, Cryo Tubes, Ext. Thread, Star Base, Sterile, Clear	100/pk	\$26.57	\$24.57
PR-CVET18-RB	1.8 mL, Cryo Tubes, Ext. Thread, Round Base, Sterile, Clear	100/pk	\$27.16	\$25.16
PR-CVET18-SB	1.8 mL, Cryo Tubes, Ext. Thread, Star Base, Sterile, Clear	100/pk	\$27.16	\$25.16
PR-CVET45-SB	4.5 mL, Cryo Tubes, Ext. Thread, Star Base, Sterile, Clear	100/pk	\$35.14	\$33.14
PR-CVIT10-SB	1.0 mL, Cryo Tubes, Int. Thread, Star Base, Sterile, Clear	100/pk	\$30.20	\$28.20
PR-CVIT18-RB	1.8 mL, Cryo Tubes, Int. Thread, Round Base, Sterile, Clear	100/pk	\$30.20	\$28.20
PR-CVIT18-SB	1.8 mL, Cryo Tubes, Int. Thread, Star Base, Sterile, Clear	100/pk	\$30.20	\$28.20
PR-CVIT36-RB	3.6 mL, Cryo Tubes, Int. Thread, Round Base, Sterile, Clear	100/pk	\$39.62	\$37.62
PR-CVIT36-SB	3.6 mL, Cryo Tubes, Int. Thread, Star Base, Sterile, Clear	100/pk	\$39.62	\$37.62
PR-CVIT45-RB	4.5 mL, Cryo Tubes, Int. Thread, Round Base, Sterile, Clear	100/pk	\$40.52	\$38.52
PR-CVIT45-SB	4.5 mL, Cryo Tubes, Int. Thread, Star Base, Sterile, Clear	100/pk	\$40.52	\$38.52
PR-CVCI	Cap Inserts Cryo Tubes, Assorted	500/pk	\$31.19	\$29.19



Note: not suitable for immersion in liquid nitrogen

Promo valid 07/01/24 - 09/30/24

COLD STORAGE ESSENTIALS

Same great Product,
with a New Name!



- Used for upright freezer or chest freezer racks
- Temperature range: -196°C to +121°C
- Drain holes on bottom and writing area on lid
- Standard 5.25 x 5.25 x 2 in.
- Available cell dividers: 144, 100, 81, 64, 49, 36, 25, or 16

PR1MA[™] Cardboard Freezer Boxes

ITEM #	DESCRIPTION	QTY	PRICE	PROMO PRICE
SB2CC-81	Standard White Cardboard 5.25 x 5.25 x 2 in. Box w/ 81 Dividers	1/ea	\$5.00	\$2.99
SB2CC-100	Standard White Cardboard 5.25 x 5.25 x 2 in. Box w/ 100 Dividers	1/ea	\$5.00	\$2.99

Promo valid 07/01/24 - 09/30/24

Need other options?
Just let us know what you need.

- Stores tubes up to 2.0 mL
- Suitable for temperatures as low as -90°C
- Rugged polypropylene
- High-contrast alphanumeric indexing on lid



Polypropylene Freezer Boxes

ITEM #	DESCRIPTION	QTY	PRICE
FB-81	81-Well 5.25 x 5.25 x 2 in. [PP], w/ Lift-Off Lid, Clear	5/pk	\$36.91
FB-81C	81-Well 5.25 x 5.25 x 2 in. [PP], w/ Lift-Off Lid, Assrt	5/pk	\$40.05

- Fit all the most popular cryo and cold storage vials
- Fit into standard freezer racks
- Temperature range: -196°C to 121°C
- Drain holes for venting and safety



Polycarbonate Freezer Box

ITEM #	DESCRIPTION	QTY	PRICE
MPB2-25	25-Well, 3 x 3 x 2.0625 in. [PC], White	1/ea	\$8.35
PB2-81	81-Well, 5.25 x 5.25 x 2 in. [PC], White	1/ea	\$10.99
PB4-81	81-Well, 5.25 x 5.25 x 3.25 in. [PC], White	1/ea	\$13.29
PB2-100	100-Well, 5.25 x 5.25 x 2 in. [PC], White	1/ea	\$10.99

Have you seen our
Freezer rack Configurator?
Scan the code to check it out!



COLD STORAGE ESSENTIALS

Haier -86°C ULT Freezers

- Microprocessor-controlled system designed for controllable range of -40°C to -86°C for cabinet space with 1°C increment

ITEM #	DESCRIPTION	PRICE
DW-86L338J	11.9 CuFt, -40°C to -86°C, Upright, 115 V/60 Hz	\$11,996.00
DW-86L578J	20.4 CuFt, -40°C to -86°C, Upright, 115 V/60 Hz	\$14,995.00
DW-86L728J	25.7 CuFt, -40°C to -86°C, Upright, 115 V/60 Hz	\$16,995.00
DW-86L828JA	29.2 CuFt, -40°C to -86°C, Upright, 220 V/60 Hz	\$17,995.00

Energy Star® Certified



For a limited time, Purchase any Haier ULT series with MIDSCI™
(TwinCool, Smart Frequency, Energy Savings Salvum J Series)

And get: 1 FREE SHELF OF STORAGE RACKS, DOCK TO DOCK SHIPPING, & RECEIVE \$500 OFF OF 2nd ULT FREEZER!

Visit www.midsci.com for more details. Promo valid 07/01-12/31/2024. All terms apply.

Freezer Inventory Racks



UF-442



UF-452



UF-452

ITEM #	TYPE	FREEZER STYLE	D x H x W	BOX HEIGHT	RACK CAPACITY (DxH)	PRICE
UFE-332	Eco	Upright	16.31 x 6.63 x 5.5 in.	2 in.	3 x 3 = 9	\$68.39
UFE-342	Eco	Upright	16.31 x 8.88 x 5.5 in.	2 in.	3 x 4 = 12	\$78.42
UFE-352	Eco	Upright	16.31 x 11.00 x 5.5 in.	2 in.	3 x 5 = 15	\$92.48
UFE-362	Eco	Upright	16.31 x 12.63 x 5.5 in.	2 in.	3 x 6 = 18	\$119.74
UFE-432	Eco	Upright	21.31 x 6.63 x 5.5 in.	2 in.	4 x 3 = 12	\$83.02
UFE-442	Eco	Upright	21.31 x 8.88 x 5.5 in.	2 in.	4 x 4 = 16	\$122.88
UFE-452	Eco	Upright	21.63 x 11.00 x 5.5 in.	2 in.	4 x 5 = 20	\$139.60
UFE-462	Eco	Upright	21.63 x 12.63 x 5.5 in.	2 in.	4 x 6 = 24	\$153.62
UFE-532	Eco	Upright	26.75 x 6.63 x 5.5 in.	2 in.	5 x 3 = 15	\$107.81
UFE-542	Eco	Upright	26.75 x 8.88 x 5.5 in.	2 in.	5 x 4 = 20	\$170.64
UFE-552	Eco	Upright	26.75 x 11.00 x 5.5 in.	2 in.	5 x 5 = 25	\$175.07
UF-332	Std	Upright	16.31 x 6.63 x 5.5 in.	2 in.	3 x 3 = 9	\$75.97
UF-342	Std	Upright	16.31 x 6.63 x 5.5 in.	2 in.	3 x 4 = 12	\$97.52
UF-352	Std	Upright	16.31 x 11.00 x 5.5 in.	2 in.	3 x 5 = 15	\$115.56
* UF-442	Std	Upright	21.63 x 8.88 x 5.5 in.	2 in.	4 x 4 = 16	\$108.97
* UF-452	Std	Upright	21.63 x 11.00 x 5.5 in.	2 in.	4 x 5 = 20	\$141.09
* UF-462	Std	Upright	21.63 x 12.63 x 5.5 in.	2 in.	4 x 6 = 24	\$192.04
UF-542	Std	Upright	26.75 x 8.88 x 5.5 in.	3 in.	5 x 4 = 20	\$155.04
UF-433	Std	Upright	21.63 x 9.38 x 5.5 in.	2 in.	4 x 3 = 12	\$96.27
CF-10-2	Std	Chest	22.00 x 5.63 x 5.5 in.	2 in.	1 x 10 = 10	\$76.98
CF-12-2	Std	Chest	26.38 x 5.63 x 5.5 in.	2 in.	1 x 12 = 12	\$92.24

*Most Popular!

Not sure what rack will maximize storage in your freezer? – No Problem.
Visit our Free Freezer Rack Configurator



COLD STORAGE ESSENTIALS

- Microprocessor temperature controller
- Touch screen digital temperature display and alarm module
- Cycle defrost
- Password-protected control panel



Templog Premier Glass Door Lab Refrigerator

ITEM #	DESCRIPTION	PRICE
ABT-HC-LP-12-TS	12 CuFt, 1°C to 10°C, 29.75 x 25.00 x 62.25 in.	\$7,692.00
ABT-HC-LP-16-TS	16 CuFt, 1°C to 10°C, 29.75 x 25.00 x 79.00 in.	\$7,834.00
ABT-HC-LP-23-TS	23 CuFt, 1°C to 10°C, 34.75 x 26.88 x 81.75 in.	\$8,086.00
ABT-HC-LP-26-TS	26 CuFt, 1°C to 10°C, 36.50 x 28.38 x 78.25 in.	\$8,484.00
ABT-HC-LP-33-TS	33 CuFt, 1°C to 10°C, 33.25 x 39.63 x 81.75 in.	\$11,876.00
ABT-HC-LP-36-TS	36 CuFt, 1°C to 10°C, 34.75 x 39.63 x 81.75 in.	\$11,456.00
ABT-HC-LP-47-TS	47 CuFt, 1°C to 10°C, 33.25 x 54.13 x 81.75 in.	\$11,746.00



- 2 Probes; 1 in air and 1 in glycol bottle
- Keyed door lock
- Six fixed shelves
- Magnetic door gasket for positive seal
- Includes monitoring system

-20°C Lab Freezers Standard Manual Defrost

ITEM #	DESCRIPTION	PRICE
ABT-HC-MFP-14	14 CuFt, -15°C to -25°C, w/ Natural Refrigerants	\$3,980.00
ABT-HC-MFP-17	17 CuFt, -15°C to -25°C, w/ Natural Refrigerants	\$4,118.00
ABT-HC-MFP-20	20 CuFt, -15°C to -25°C, w/ Natural Refrigerants	\$4,288.00

**Don't see what you need?
Call us, we have many more options available for you!**

- Dual independent refrigeration systems for maximum sample safety
- High speed refrigeration system for faster pull down and temperature recovery after door opening
- Smart full-size touchscreen and wireless monitoring connectivity

Haier -86°C TwinCool ULT Freezers

ITEM #	DESCRIPTION	PRICE
DW-86L578SAT	20.4 CuFt, -40°C to -86°C, Upright, 220 V/60 Hz	\$21,995.00
DW-86L578ST	20.4 CuFt, -40°C to -86°C, Upright, 115 V/60 Hz	\$21,995.00
DW-86L728SAT	25.7 CuFt, -40°C to -86°C, Upright, 220 V/60 Hz	\$23,895.00
DW-86L728ST	25.7 CuFt, -40°C to -86°C, Upright, 115 V/60 Hz	\$23,895.00

Energy Star® Certified



For a limited time, Purchase any Haier ULT series with MIDSCI™
(TwinCool, Smart Frequency, Energy Savings Salvum J Series)

And get: 1 FREE SHELF OF STORAGE RACKS, DOCK TO DOCK SHIPPING, & RECEIVE \$500 OFF OF 2nd ULT FREEZER!

Visit www.midsci.com for more details. Promo valid 07/01-12/31/2024. All terms apply.

DROSOPHILA ESSENTIALS

- Polystyrene: highly transparent and conducive to sunlight exposure, easy to observe and store, cost-effective
- Disposable vials are safe and cost-effective alternatives to glass vials
- Available in bulk or tray packing

Drosophila Vials - Polystyrene



ITEM #	DESCRIPTION	QTY	PRICE
DRK-51-0500	Narrow Vials, 25 x 95 mm, Bulk	25/bg, 20 bgs/cs	\$87.03
DRK-51-0501	Narrow Vials Combo, 25 x 95 mm, Tray	100/pk, 5 pks/cs	\$104.84
DRK-51-0800	Wide Vials, 28.5 x 95 mm, Bulk	25/bg, 20 bgs/cs	\$90.99
DRK-51-0801	Wide Vials Combo, 28.5 x 95 mm, Tray	100/pk, 5 pks/cs	\$110.77

- Durable and suitable for drosophila growth and development in most situations
- Reusable
- Available in bulk or tray packing
- Polypropylene: firm, clear, flexible, corrosion resistant, and fully autoclavable (121°C, 15 psi, 15 min)



Drosophila Vials - Polypropylene

ITEM #	DESCRIPTION	QTY	PRICE
DRK-51-0510	Narrow Vials, 25 x 95 mm, Bulk	25/bg, 20 bgs/cs	\$71.21
DRK-51-0511	Narrow Vials Combo, 25 x 95 mm, Tray	100/pk, 5 pks/cs	\$92.97
DRK-51-0810	Wide Vials, 28.5 x 95 mm, Bulk	25/bg, 20 bgs/cs	\$81.10
DRK-51-0811	Wide Vials Combo, 28.5 x 95 mm, Tray	100/pk, 5 pks/cs	\$89.01

- K-Resin: High transparency, good toughness, and can be sterilized by ethylene oxide gas or ultraviolet rays
- Glass-like clarity for sample observation
- Available in bulk or tray packing

Drosophila Vials - K-Resin



ITEM #	DESCRIPTION	QTY	PRICE
DRK-51-0520	Narrow Vials, 25 x 95 mm, Bulk	25/bg, 20 bgs/cs	\$96.92
DRK-51-0521	Narrow Vials Combo, 25 x 95 mm, Tray	100/pk, 5 pks/cs	\$106.81
DRK-51-0820	Wide Vials, 28.5 x 95 mm, Bulk	25/bg, 20 bgs/cs	\$100.88
DRK-51-0821	Wide Vials Combo, 28.5 x 95 mm, Tray	100/pk, 5 pks/cs	\$120.66

Vials, bottles and closures sold separately.

DROSOPHILA ESSENTIALS

- 6 oz square bottom
- Polypropylene and Polyethylene
- Specially designed neck to prevent the plug from falling in
- Universal closure compatibility



Drosophila Bottles - Square Bottom

ITEM #	DESCRIPTION	QTY	PRICE
DRK-51-17720	6 oz/177 mL [PP], 57 x 57 x 103 mm, Bulk	200/bg, 1 bg/cs	\$81.10
DRK-51-17721	6 oz/177 mL [PP], 57 x 57 x 103 mm, Tray	25/pk, 8 pks/cs	\$108.79
DRK-51-17730	6 oz/177 mL [PE], 57 x 57 x 103 mm, Bulk	200/bg, 1 bg/cs	\$85.05
DRK-51-17731	6 oz/177 mL [PE], 57 x 57 x 103 mm, Tray	25/pk, 8 pks/cs	\$114.73




- Polypropylene and Polyethylene
- Specially designed neck to prevent the plug from falling in
- Universal closure compatibility

Drosophila Bottles - Round Bottom

ITEM #	DESCRIPTION	QTY	PRICE
DRK-51-22720	8 oz/227 mL [PP], 60 x 130 mm, Bulk	100/bg, 1 bg/cs	\$61.32
DRK-51-22721	8 oz/227 mL [PP], 60 x 130 mm, Tray	25/pk, 4 pks/cs	\$79.12
DRK-51-22730	8 oz/227 mL [PE], 60 x 130 mm, Bulk	100/bg, 1 bg/cs	\$63.30
DRK-51-22731	8 oz/227 mL [PE], 60 x 130 mm, Tray	25/pk, 4 pks/cs	\$83.08

Drosophila Closures

ITEM #	DESCRIPTION	QTY	PRICE
DRK-51-0505	Narrow Vial Closures, Densest Foam, 30 x 45 mm, Bulk	100/bg, 10 bgs/cs	\$239.34
DRK-51-0515	Narrow Vial Closures, Densest Sponge, 25 x 25 mm, Bulk	100/bg, 10 bgs/cs	\$278.90
DRK-51-0525	Narrow Vial Closures, Degreasing Fiber, 24 x 30 mm, Bulk	100/bg, 10 bgs/cs	\$104.84
DRK-51-0805	Wide Vial Closures, Densest Foam, 35 x 45 mm, Bulk	100/bg, 10 bgs/cs	\$259.12
DRK-51-0815	Wide Vial Closures, Densest Sponge, 28 x 30 mm, Bulk	100/bg, 10 bgs/cs	\$358.02
DRK-51-0825	Wide Vial Closures, Degreasing Fiber, 27 x 30 mm, Bulk	100/bg, 10 bgs/cs	\$116.70
DRK-51-17725	Bottle Closures, Densest Foam, 45 x 45 mm Bulk	100/bg, 10 bgs/cs	\$239.34
DRK-51-17735	Bottle Closures, Densest Sponge, 40 x 30 mm, Bulk	100/bg, 10 bgs/cs	\$292.75



Vials, bottles and closures sold separately.

DROSOPHILA ESSENTIALS

Media Components

ITEM #	DESCRIPTION	QTY	PRICE
KCA20030-500	Agar, Bacteriological Grade, Powder	500 g	\$195.29
KCY20020-500	Yeast Extract, Powder	500 g	\$86.25
KC248700	Water Molecular Biology Grade	1 gal	\$64.90



Available in multiple volumes.
Contact your trusted MIDSCI Lab Consultant for more options!

- Our scopes are Greenough Optical designed
- Provide the enhanced depth of focus, color reproduction resolution, and image quality possible
- Building block designs to allow for selectable configurations and customization to meet specific requirements and budgets
- Optional LED illumination available



Stereo Microscopes

ITEM #	DESCRIPTION	QTY	PRICE
Z645	0.67X - 4.5X Continuous Zoom, Binocular, Platform Stand	1/ea	\$1,677.00
Z730	0.7x - 3.5x Continuous Zoom, Binocular, Platform Stand	1/ea	\$1,374.00
Z850	0.8X - 5.0X Continuous Zoom, Binocular, Platform Stand	1/ea	\$2,060.00

Contact your trusted MIDSCI™ Lab Consultant for a quote!

- Precise Temperature Control - Superior Uniformity
- Independent Over Temperature Thermostat
- Door Ajar Alarm
- Temperature Range 15°C to 40°C
- 2-year warranty

Shel-Lab Drosophila Incubator

ITEM #	DESCRIPTION	QTY	PRICE
SRI20PF	19.3 cu.ft., Self Standing, Programmable Day/Night Temp Control	1/ea	\$10,218.34



MICROBIOLOGY SUPPLIES

- Stackable design for more stability
- Slideable design for automated filling and streaking systems
- Features side arrows and ISO mark target for automation
- High quality, value priced product




PR1MA[™] Petri Dishes

ITEM #	DESCRIPTION	QTY	PRICE
667665	60 x 15 mm, [PS] Sterile	25/pk, 20 pks/cs	\$133.19
667694	100 x 15 mm, [PS] Sterile, Slideable	25/pk, 20 pks/cs	\$120.22
667695	100 x 15 mm, [PS] Sterile, Stackable	25/pk, 20 pks/cs	\$120.22
667697	100 x 15 mm, [PS] Sterile, Slippable	25/pk, 20 pks/cs	\$123.74
667682	100 x 15 mm, [PS] Sterile, 2 Compartment	25/pk, 20 pks/cs	\$166.81
667683	100 x 15 mm, [PS] Sterile, 3 Compartment	25/pk, 20 pks/cs	\$169.45
667684	100 x 15 mm, [PS] Sterile, 4 Compartment	25/pk, 20 pks/cs	\$172.31



DuoClick[™] Two Position Screw-Cap Culture Tubes

ITEM #	DESCRIPTION	QTY	PRICE
T8730	5 mL, w/ Screw Cap, Clear [PP], Sterile, Bagged	25/bg, 20 bgs/cs	\$144.45
T8732	5 mL, w/ Screw Cap, Clear [PP], Sterile, Racked	50/rk, 10 rk/cs	\$170.13
T8734	5 mL, w/ Screw Cap, Clear [PP], Sterile, Individ. Wrapped	1/pk, 500 pks/cs	\$207.58
T8750	5 mL, w/ Screw Cap, Clear [PP], Non-Sterile, Bulk	1,000/cs	\$130.80
T8830	14 mL, w/ Screw Cap, Clear [PP], Sterile, Bagged	25/bg, 20 bgs/cs	\$207.58
T8832	14 mL, w/ Screw Cap, Clear [PP], Sterile, Racked	50/rk, 10 rk/cs	\$255.73
T8834	14 mL, w/ Screw Cap, Clear [PP], Sterile, Individ. Wrapped	1/pk, 500 pks/cs	\$277.13
T8850	14 mL, w/ Screw Cap, Clear [PP], Non-Sterile, Bulk	1,000/cs	\$223.63



MICROBIOLOGY SUPPLIES

- Polycarbonate (PC) body, meets USP Class 6 Standard
- The wide-mouth neck design is convenient for sampling and liquid filling
- Sterile, RNase, DNase, pyrogen free
- Packaged in sterile bags, Sterilized by E-Beam
- Seal cap and vent cap options available




NEST® Conical Erlenmeyer Flasks

ITEM #	DESCRIPTION	QTY	PRICE
785001	2 L Conical, HDPE Seal Cap, Sterile, Indiv. Wrapped	6/pk	\$224.65
786001	3 L Conical, HDPE Seal Cap, Sterile, Indiv. Wrapped	4/pk	\$263.15
785011	2 L Conical, HDPE Vent Cap, Sterile, Indiv. Wrapped	6/pk	\$230.95
786011	3 L Conical, HDPE Vent Cap, Sterile, Indiv. Wrapped	4/pk	\$266.30


- Spread cultures and liquids evenly
- Smooth-edges – minimize "gouging"

- 3 configurations: "L", "T", or Triangle
- Ideal for bacterial CFU Counts
- Pre-sterilized – no flaming required



Disposable Spreaders

ITEM #	DESCRIPTION	QTY	PRICE
LS-SIP	L-Shape, 38 x 156 mm, Sterile, Indiv. Wrapped	1/pk, 500/cs	\$111.55
LS-SP	L-Shape, 38 x 156 mm, Sterile	10/pk, 500/cs	\$95.45
TS-SIP	T-Shape, 34 x 140 mm, Sterile, Indiv. Wrapped	1/pk, 500/cs	\$113.85
TS-SP	T-Shape, 34 x 140 mm, Sterile	10/pk, 500/cs	\$105.86
TR-SB-S	Triangle Shape, 30 x 207 mm, Sterile	10/bg, 100/cs	\$82.80
TR-SB-L	Triangle Shape, 60 x 235 mm, Sterile	10/bg, 100/cs	\$82.80



Sterile Disposable Inoculating Loops & Needles

ITEM #	DESCRIPTION	QTY	PRICE
M4015	Loop size: 1uL, Length: 200 mm, 20 Per Peel, 50 Per Pack	1,000/cs	\$75.90
M4025	Loop size: 10uL, Length: 200 mm, 20 Per Peel, 50 Per Pack	1,000/cs	\$75.90
M4030	200 mm, 20 Per Peel, 50 Per Pack	1,000/cs	\$75.90

MICROBIOLOGY SUPPLIES

PR1MA™ Growth

Available in Capsules, Granulated, and Powder Form

Agar, High Gel Strength

ITEM #	QTY	PRICE
Powder		
KCA20020-500	500 g	\$195.29
KCA20020-1000	1 kg	\$323.08
KCA20020-5000	5 kg	\$1,053.46
Granulated		
KCA20250-500	500 g	\$231.63
KCA20250-1000	1 kg	\$409.90
KCA20250-5000	5 kg	\$1,127.75

Agar, Noble

ITEM #	QTY	PRICE
KCA20460-100	100 g	\$128.37
KCA20460-500	500 g	\$459.52
KCA20460-1000	1 kg	\$874.04

Agar, Pytho

ITEM #	QTY	PRICE
KCA20300-1000	1 kg	\$293.94

Brain Heart Infusion Broth

Brain Heart Infusion:	17.5 g/l
Pancreatic Digest, Gelatin:	10 g/l
Dextrose:	2.0 g/l
Sodium Chloride:	5.0 g/l
Disodium Phosphate:	2.5 g/l

ITEM #	QTY	PRICE
KCB11000-500	500 g	\$140.19
KCB11000-1000	1 kg	\$231.63
KCB11000-5000	5 kg	\$1,005.87

Casamino Acid CAS 77-92-9

ITEM #	QTY	PRICE
KCC45000-500	500 g	\$128.37
KCC45000-1000	1 kg	\$251.83
KCC45000-2000	2 kg	\$492.40

Gelrite™ CAS 71010-52-1

ITEM #	QTY	PRICE
KCG35020-100	100 g	\$78.75
KCG35020-250	250 g	\$133.27
KCG35020-500	500 g	\$223.85

Glycerol

CAS 56-81-5 | C₃H₈O₃ | MW 92.1
Purity: >99%

ITEM #	QTY	PRICE
KCG22020-0.5	500 mL	\$38.65
KCG22020-1.0	1 L	\$51.35
KCG22020-4.0	4 L	\$166.15

Lactobacillus MRS Broth

Peptic Digest of Animal Tissue:	10 g/l
Beef Extract:	10 g/l
Yeast Extract:	5 g/l
Ammonium Citrate:	2 g/l
Dextrose:	20 g/l
Polysorbate 80:	1 g/l
Sodium acetate:	5 g/l
Magnesium Sulfate:	0.1 g/l
Manganese Sulfate:	0.05 g/l
Potassium Phosphate, Dibasic:	2 g/l

ITEM #	QTY	PRICE
KCL11000-500	500 g	\$109.90
KCL11000-1000	1 kg	\$186.35

LB Agar [Lennox L Agar]

Casein Digest Peptone:	10 g/l
Yeast Extract:	5 g/l
Sodium Chloride:	5 g/l
Agar:	15 g/l

ITEM #	QTY	PRICE
Powder		
KCL24030-500	500 g	\$93.46
KCL24030-1000	1 kg	\$178.27
KCL24030-2000	2 kg	\$331.44
Granulated		
KCL24033-100	100 g	\$47.31
KCL24033-500	500 g	\$131.83
KCL24033-1000	1 kg	\$231.63
KCL24033-2000	2 kg	\$386.54

LB Broth [Lennox L Broth]

Casein Digest Peptone:	10 g/l
Yeast Extract:	5 g/l
Sodium Chloride:	5 g/l

ITEM #	QTY	PRICE
Powder		
KCL24060-500	500 g	\$67.21
KCL24060-1000	1 kg	\$124.33
KCL24060-5000	5 kg	\$543.17
Granulated		
KCL24066-500	500 g	\$86.25
KCL24066-1000	1 kg	\$147.69
KCL24066-5000	5 kg	\$696.92

LB Broth [Lennox L Broth]

1L Capsules		
Tryptone:	10 g/l	
Yeast Extract:	5 g/l	
Sodium Chloride:	5 g/l	

ITEM #	QTY	PRICE
KCL24065-500	500 g	\$111.06
KCL24065-1000	1 kg	\$202.50
KCL24065-5000	5 kg	\$906.92

LB Broth [Lennox L Broth]

1g Capsules		
Tryptone:	10 g/l	
Yeast Extract:	5 g/l	
Sodium Chloride:	5 g/l	

ITEM #	QTY	PRICE
KCL24061-500	500 g	\$140.19
KCL24061-1000	1 kg	\$267.69

Luria Agar [Miller's LB Agar]

Capsules		
Tryptone:	10 g/l	
Yeast Extract:	5 g/l	
Sodium Chloride:	10 g/l	
Agar:	15 g/l	

ITEM #	QTY	PRICE
1g Capsules		
KCL24021-500	500 g	\$195.29
KCL24021-1000	1 kg	\$360.29

500mL Capsules		
KCL24025-500	500 g	\$178.27
KCL24025-1000	1 kg	\$277.79

Luria Broth [Miller's LB Broth]

Casein Digest Peptone:	10 g/l
Yeast Extract:	5 g/l
Sodium Chloride:	10 g/l

ITEM #	QTY	PRICE
Powder		
KCL24040-500	500 g	\$67.21
KCL24040-1000	1 kg	\$124.33
KCL24040-5000	5 kg	\$543.17

Granulated		
KCL24400-500	500 g	\$86.25
KCL24400-1000	1 kg	\$147.69
KCL24400-5000	5 kg	\$696.92

Luria Broth [Miller's LB Broth]

Capsules		
Tryptone:	10 g/l	
Yeast Extract:	5 g/l	
Sodium Chloride:	10 g/l	

ITEM #	QTY	PRICE
1g Capsules		
KCL24041-500	500 g	\$140.19
KCL24041-1000	1 kg	\$267.69

1L Capsules		
KCL24045-500	500 g	\$111.06
KCL24045-1000	1 kg	\$202.50
KCL24045-5000	5 kg	\$897.98

ORDER TODAY!

midsci.com • custserv@midsci.com • 800.227.9997 • f: 636.225.9998 25

MICROBIOLOGY SUPPLIES

PR1MA™ Growth

Available in Capsules, Granulated, and Powder Form

Nutrient Broth

Beef Extract:	3 g/l
Yeast Extract:	5 g/l
Sodium Chloride:	10 g/l

ITEM #	QTY	PRICE
KCN15100-500	500 g	\$93.46
KCN15100-1000	1 kg	\$159.23
KCN15100-5000	5 kg	\$609.52

NZCYM Broth

Casein Digest Protein:	10 g/l
Casamio Acid:	1 g/l
Sodium Chloride:	5 g/l
Yeast Extract:	5 g/l
Magnesium Chloride:	0.98 g/l

ITEM #	QTY	PRICE
KCN72000-500	500 g	\$77.88
KCN72000-1000	1 kg	\$135.58

Malt Extract

Maltose:	12.75 g/l
Dextrin:	2.75 g/l
Glycerol:	2.35 g/l
Peptone:	0.78 g/l
Agar 1:	5.0 g/l
Final pH(25°C):	4.6±0.2 g/l

ITEM #	QTY	PRICE
KCM21500-500	500 g	\$142.79
KCM21500-1000	1 kg	\$233.94
KCM21500-2000	2 kg	\$443.65

Peptone CAS 73049-73-79

ITEM #	QTY	PRICE
Powder		
KCP20240-500	500 g	\$67.21
KCP20240-1000	1 kg	\$124.33
KCP20240-5000	5 kg	\$543.17
Granulated		
KCP20250-500	500 g	\$86.25
KCP20250-1000	1 kg	\$147.69
KCP20250-5000	5 kg	\$696.92

Peptone, Soy CAS 73049-73-79

ITEM #	QTY	PRICE
Powder		
KCP50300-500	500 g	\$60.00

SOB Broth

Casein Enzymic Hydrolysate:	20 g/l
Yeast Extract:	5 g/l
Sodium Chloride:	0.5 g/l
Magnesium Sulfate:	2.4 g/l
Potassium Chloride:	0.186 g/l

ITEM #	QTY	PRICE
Powder		
KCS25000-100	100 g	\$42.98
KCS25000-500	500 g	\$109.90

500mL Capsules

KCS25040-500	500 g	\$140.19
KCS25040-1000	1 kg	\$231.63

Sodium Chloride

CAS 7647-14-5 | NaCl | MW 58.4
Purity: >99.9%

ITEM #	QTY	PRICE
KCS23020-1000	1 kg	\$38.65
KCS23020-5000	5 kg	\$102.12
KCS23020-10000	10 kg	\$194.13

Terrific Broth Modified

Casein Digest Peptone:	12 g/l
Yeast Extract:	24 g/l
Dipotassium phosphate:	9.5 g/l
Monopotassium phosphate:	2.2 g/l

ITEM #	QTY	PRICE
Powder		
KCT15000-500	500 g	\$67.21
KCT15000-1000	1 kg	\$124.33
KCT15000-5000	5 kg	\$543.17
Granulated		
KCT15100-500	500 g	\$86.25
KCT15100-1000	1 kg	\$147.69
KCT15100-5000	5 kg	\$696.92

Terrific Broth Modified

Capsules

Casein Digest Peptone:	12 g/l
Yeast Extract:	24 g/l
Dipotassium phosphate:	9.5 g/l
Monopotassium phosphate:	2.2 g/l

ITEM #	QTY	PRICE
500mL Capsules		
KCT15050-600	600 g	\$140.19
KCT15050-1200	1.2 kg	\$231.63
KCT15050-5000	5 kg	\$897.98

Tryptone CAS 91079-40-2

Total Nitrogen (TN):	12.7%
Amino Nitrogen (AN):	3.7%
AN/TN Ratio:	0.26
Sodium Chloride:	0.4%
Tryptophan:	1.0

ITEM #	QTY	PRICE
Powder		
KCT60060-500	500 g	\$86.25
KCT60060-1000	1 kg	\$140.20
KCT60060-5000	5 kg	\$514.33

Granulated

KCT60065-500	500 g	\$102.12
KCT60065-1000	1 kg	\$166.15
KCT60065-5000	5 kg	\$686.83

2XYT Broth

ITEM #	QTY	PRICE
Powder		
KCX15600-500	500 g	\$67.21
KCX15600-1000	1 kg	\$124.33
KCX15600-5000	5 kg	\$543.17

Granulated

KCX15680-500	500 g	\$93.46
KCX15680-1000	1 kg	\$157.79
KCX15680-5000	5 kg	\$751.73

500mL Capsules

KCX15640-500	500 g	\$124.33
KCX15640-1000	1 kg	\$222.69
KCX15640-5000	5 kg	\$1,053.46

Yeast Extract CAS 8013-01-2

ITEM #	QTY	PRICE
Powder		
KCY20020-500	500 g	\$86.25
KCY20020-1000	1 kg	\$144.52
KCY20020-5000	5 kg	\$600.00

Granulated

KCY20025-500	500 g	\$111.06
KCY20025-1000	1 kg	\$178.27
KCY20025-5000	5 kg	\$743.37

Yeast Nitrogen Base without Amino Acid

ITEM #	QTY	PRICE
KCY20040-100	100 g	\$86.25
KCY20040-500	500 g	\$304.04
KCY20040-1000	1 kg	\$514.33

Free Shipping on
PR1MA™ Chemicals!

Find More Growth
Media Chemicals Online!

MICROBIOLOGY SUPPLIES

- Digital control of temperature and CO₂
- Temp range ambient -10°C to 55°C
- Accepts bottles and flasks up to 2 L
- Includes 2 stainless steel shelves
- Dual beam, infrared sensor
- 2-year warranty

PR1MA™ Mini CO₂ Incubator

ITEM #	DESCRIPTION	QTY	PRICE
INCO2-H	13 x 14.5 in. (20 L) Incubator	1/ea	\$4,518.82
H2300-REG	Optional CO ₂ Gas Regulator	1/ea	\$385.53



- Redesigned to accept two 1 L Flasks
- Increased temperature control up to 70°C
- 2-year warranty



PR1MA™ Mini Incubated Shaker

ITEM #	DESCRIPTION	QTY	PRICE
IS-300	Mini Incubated Shaker, Ambient + 5°C to 70°C, Max Speed 300 rpm, 9.5 x 11.5 in. Platform	1/ea	\$3,598.97
IS-PLAT	MAGic Clamp™ Platform For PR1MA™ Mini Incubated Shaker	1/ea	\$296.13
H1000-MR-CMB	MAGic Clamp™ Combo, 125 mL, 250 mL, 500 mL Flasks	1/ea	\$52.63

Platforms, clamps and racks available for 1.5 mL to 2 L and microplates

- Exceptional temperature uniformity
- Touch screen control over time, speed and temperature
- Horizontally circular 19 mm orbit for aeration & mixing
- 2-year warranty

10 L Incubated Shaker

ITEM #	DESCRIPTION	QTY	PRICE
IS-1010	Incubated Shaker, Ambient +5°C to 70°C, Max Speed 300 rpm, 18 x 18 in. Platform	1/ea	\$7,495.31
IS-1010R	Refrigerated Incubated Shaker, Ambient -15°C to 60°C, Max Speed 300 rpm, 18 x 18 in. Platform	1/ea	\$10,596.02
IS-1010-PLAT	MAGic Clamp™ Platform For PR1MA™ Incubated Shaker	1/ea	\$413.65
H1000-MR-CMB	MAGic Clamp™ Combo, 125 mL, 250 mL, 500 mL Flasks	1/ea	\$52.63

Platforms, clamps and racks available for 1.5 mL to 2 L and microplates



MAGic Clamp™ system allows for fast, customizable sharing options!



Never screw in a clamp again!

MICROBIOLOGY SUPPLIES

- Precision temperature control, up to 75°C
- Dual-Convection technology, excellent uniformity and recovery
- Includes SureCheck™ temperature logger and software
- Stainless steel chamber with internal power outlet

SureTemp™ Incubators



ITEM #	DESCRIPTION	QTY	PRICE
H2505-40	40 L, Dual Convection Incubator	1/ea	\$2,019.58
H2505-70	70 L, Dual Convection Incubator	1/ea	\$2,603.43
H2505-130	130 L, Dual Convection Incubator	1/ea	\$3,471.62

- Designed specifically for the sterilization of laboratory research tools
- Fully automatic operation – just press start
- Internal memory for retrieving cycle data
- Built in water tank - no external connection needed
- 8 L, 16 L and 28 L models - all models small enough for benchtop use



BioClave™ Research Autoclaves

ITEM #	DESCRIPTION	QTY	PRICE
RA-8L	8 L Steam Sterilizer, 115V	1/ea	\$5,040.61
B4000-18	18 L Steam Sterilizer, 115V	1/ea	\$6,973.82
B4000-28-120	28 L Steam Sterilizer, 120V	1/ea	\$12,842.86

- Digital temperature control
- Personal-sized, 13 x 14.5 inch footprint (20 L)
- Accepts bottles and flasks up to 2 L
- Internal outlet, for installation of shaker, rockers, etc.

PR1MA™ Mini Digital Incubators

ITEM #	DESCRIPTION	QTY	PRICE
IN20-H	Heat Only (Ambient +5 to 60°C), Incl. 2 Shelves	1/ea	\$596.55
IN20-HC	Mini Digital Incubator - Heat and Cool, Incl. 2 Shelves	1/ea	\$907.43
H2200-SH	Extra Shelf, Plastic	1/ea	\$47.69



LAB ESSENTIALS



PR1MA™ Premium Microcentrifuge Tubes

ITEM #	DESCRIPTION	QTY	PRICE
PR-MCT6	0.5 mL, Snap Cap, Boil Proof, Clear, Non-Sterile	1,000/pk	\$32.59
PR-MCT6-LR	0.5 mL, Low-Retention, Snap Cap, Boil Proof, Clear, Non-Sterile	500/pk	\$24.13
PR-MCT6-S	0.5 mL, Snap Cap, Boil Proof, Clear, Sterile	500/pk	\$25.11
PR-MCT15	1.5 mL, Snap Cap, Boil Proof, Clear, Non-Sterile	500/pk	\$19.32
PR-MCT15-LR	1.5 mL, Low-Retention, Snap Cap, Boil Proof, Clear, Non-Sterile	250/pk	\$13.78
PR-MCT15-S	1.5 mL, Snap Cap, Boil Proof, Clear, Sterile	250/pk	\$14.51
PR-MCT17	1.5 mL, Large Snap Cap, Boil Resist, Clear, Non-Sterile	500/pk	\$19.32
PR-MCT2	2 mL, Snap Cap, Boil Proof, Clear, Non-Sterile	500/pk	\$21.84
PR-MCT2-LR	2 mL, Low-Retention, Snap Cap, Boil Proof, Clear, Non-Sterile	250/pk	\$17.48
PR-MCT2-S	2 mL, Snap Cap, Boil Proof, Clear, Sterile	250/pk	\$14.71

For colors add:

-B (blue), -G (green), -O (orange), -R (red), -P (purple), -Y (yellow), -RB (assorted), or -RBF (assorted neon) to the Item #
Sterile available for clear tubes only. Low-retention available for clear, non-sterile tubes only.

LARGE CAP: BOIL-RESISTANT, EASY TO OPEN AND CLOSE, AND A LARGE LABELING SURFACE

SMALL CAP: BOIL-PROOF, IDEAL FOR RIGOROUS HIGH-TEMP APPLICATIONS

OPTICAL CLARITY ENSURES MAXIMUM SAMPLE VISIBILITY

RATED TO 20,000 XG

GRADUATIONS EVERY 100 µL

AUTOCCLAVABLE

LARGE WRITING SURFACE ON CAP TOP



Lots more to choose from on our website!

*Certified free of detectable RNase, DNase, DNA, PCR inhibitors and tested pyrogen-free.
 Sterile product is certified pyrogen-free.*

LAB ESSENTIALS



- Outstanding for high centrifugation or ultra low freezing
- Liquid & temp resistant label on screw cap
- Leak proof
- White graduations down to 0.1 mL & 0.5 mL respectively
- Specs: -190°C to 100°C & 15,500 xg [PP]
-30°C to +70°C & 1,700 xg [PS], Sterile

TPP® Centrifuge Tubes



ITEM #	DESCRIPTION	QTY	PRICE
TP91016	13 mL, Round (17 x 105 mm) [PP], Sterile	40/pk, 20 pks/cs	\$748.22
TP91019	13 mL, Flat (17 x 100 mm) [PP], Sterile	40/pk, 20 pks/cs	\$748.87
TP91017	15 mL, Round (17 x 120 mm) [PP], Sterile	40/pk, 20 pks/cs	\$779.16
TP91015	15 mL, Conical, [PP], Sterile	40/pk, 20 pks/cs	\$653.83
TP91115	15 mL, Conical, [PS], Sterile	40/pk, 20 pks/cs	\$827.31
TP91050	50 mL, Conical, [PP], Sterile	20/pk, 18 pks/cs	\$351.23
TP91051	50 mL, Skirted Conical, [PP], Sterile	20/pk, 16 pks/cs	\$339.49
TP91056	50 mL, Round, [PP], Sterile	20/pk, 18 pks/cs	\$229.36

- DNase, RNase and pyrogen-FREE
- White shaded writing area

- Black graduations
- Temp Range: -80°C - 120°C

PR1MA™ Conical Centrifuge Tubes

ITEM #	DESCRIPTION	QTY	PRICE
C15B	15 mL [PP], Screw Cap, Bagged, Sterile	500/cs	\$155.82
C15R	15 mL [PP], Screw Cap, Racked, Sterile	500/cs	\$142.62
C50B	50 mL [PP], Screw Cap, Bagged, Sterile	500/cs	\$183.52
C50R	50 mL [PP], Screw Cap, Racked, Sterile	500/cs	\$176.27

**Tornado is changing to
PR1MA™!**
Same great
product, with new
name!



The 15 mL tubes can withstand up to 12,000 (RCF) g, and the 50 mL can withstand up to 16,000 (RCF) g.



- Fits most 50 mL rotors - Rim eliminates need for adapter
- Printed graduations - From 5 mL to 2 mL
- Compact low-profile design - Minimizes shipping, storage and waste
- Gamma sterilized - RNase, DNase, & Pyrogen free

25 mL Centrifuge Tubes

ITEM #	DESCRIPTION	QTY	PRICE
C2625-B	25 mL (29 x 77 mm) [PP], Screw Cap, Sterile, Bagged	25/pk, 8 pks/cs	\$101.20
C2625-R	25 mL (29 x 77 mm) [PP], Screw Cap, Sterile, Racked	25/pk, 8 pks/cs	\$104.65

LAB ESSENTIALS



PR1MA™ PCR Strips with Attached Caps

ITEM #	DESCRIPTION	QTY	PRICE
PR-PCR28ACD	0.2 mL x 8 Strip PCR Tubes, Clear w/ Indiv. Attached Dome Cap	120/cs	\$68.86
PR-PCR28ACF	0.2 mL x 8 Strip PCR Tubes, Clear w/ Indiv. Attached Flat Cap	120/cs	\$68.86
PR-PCR28ACFF	0.2 mL x 8 Strip "Easy Split" PCR Tubes, Clear w/ Indiv. Attached Flat Cap	120/cs	\$68.86
PR-PCR28ACZD	0.2 mL x 8 Strip PCR Tubes, Clear w/ Hinge Attached Dome Cap	125/cs	\$81.73
PR-PCR28ACZF	0.2 mL x 8 Strip PCR Tubes, Clear w/ Hinge Attached Flat Cap	125/cs	\$71.60

Tubes and caps are certified free of detectable RNase, DNase, PCR inhibitors, & pyrogen free.

PR1MA™ Individual PCR Tubes

ITEM #	DESCRIPTION	QTY	PRICE
PR-PCR2F	0.2 mL Individual PCR Tube w/Clear Flat Cap	1,000/cs	\$38.80
PR-PCR2FF	0.2 mL Individual PCR Tube w/Frosted Flat Cap	1,000/cs	\$38.80
PR-PCR2D	0.2 mL Individual PCR Tube w/Domed Cap	1,000/cs	\$45.33

Tubes and caps are certified free of detectable RNase, DNase, PCR inhibitors, & pyrogen free.



Not sure what will fit your equipment? **We have a solution for that!**
Check out our new **PR1MA™** Plates and Tubes Selection Tool



PR1MA™ Serological Pipettes

ITEM #	DESCRIPTION	QTY	PRICE
PR-SERO-1	1 mL, Indiv. Wrapped, Sterile	800/cs	\$120.44
PR-SERO-2	2 mL, Indiv. Wrapped, Sterile	600/cs	\$136.46
PR-SERO-5	5 mL, Indiv. Wrapped, Sterile	200/cs	\$74.19
PR-SERO-10	10 mL, Indiv. Wrapped, Sterile	200/cs	\$80.93
PR-SERO-25	25 mL, Indiv. Wrapped, Sterile	200/cs	\$99.12
PR-SERO-50	50 mL, Indiv. Wrapped, Sterile	90/cs	\$99.78

LAB ESSENTIALS

Supreme accuracy with the low-retention feature!

- Barrier tips have extra working volume, so liquid doesn't touch the filter
- Great stackable racks
- 10 μ L XL tips have 1 μ L graduations
- Clear & readable graduations
- Sterile, RNase free and DNase free



Improve consistency and reproducibility with low-retention tips!

PRIMA[™] Low-Retention Barrier Pipette Tips

ITEM #	DESCRIPTION	QTY	PRICE
PR-10RK-FL	10 μ L, Sterile	96/rk, 10 rks/cs	\$67.16
PR-10XLRK-FL	10 μ L XL, Sterile	96/rk, 10 rks/cs	\$66.96
PR-20RK-FL	20 μ L, Sterile	96/rk, 10 rks/cs	\$66.96
PR-100RK-FL	100 μ L, Sterile	96/rk, 10 rks/cs	\$66.96
PR-200RK-FL	200 μ L, Sterile	96/rk, 10 rks/cs	\$66.96
PR-200XLRK-FL	200 μ L XL 5 in. Long, Sterile	96/rk, 8 rks/cs	\$84.52
PR-300RK-FL	300 μ L, Sterile	96/rk, 10 rks/cs	\$66.96
PR-1000RK-FL	1000 μ L, Sterile	96/rk, 8 rks/cs	\$58.53
PR-1250RK-FL	1250 μ L, Sterile	96/rk, 8 rks/cs	\$55.41

1250 μ L

200 μ L XL



1 μ L graduation on 10 μ L XL tip

Graduations:

10 μ L: 2 & 10 μ L	200 μ L: 10, 50, & 100 μ L
10 μ L XL: 1, 2, & 10 μ L	300 μ L XL: 10, 50, 100, & 200 μ L
20 μ L: 10 μ L	1000 μ L: 100, 500, 1000 μ L
100 μ L: 10, 50, & 100 μ L	1250 μ L: 100, 500, & 1000 μ L

LAB ESSENTIALS



PR1MA™ Pipette Tip Reloads

ITEM #	DESCRIPTION	QTY	PRICE
PR-10RFL	10 µL, Non-Sterile	96/rk, 10 rks/cs	\$27.85
PR-10XLRFL	10 µL XL, Non-Sterile	96/rk, 10 rks/cs	\$27.85
PR-200RFL	200 µL, Non-Sterile	96/rk, 10 rks/cs	\$26.59
PR-300RFL	300 µL, Non-Sterile	96/rk, 10 rks/cs	\$29.58
PR-1250RFL	1250 µL, Non-Sterile	96/rk, 8 rks/cs	\$36.61

Racks only, no boxes included

PR1MA™ Non-Sterile Racked Pipette Tips

ITEM #	DESCRIPTION	QTY	PRICE
PR-10RK-NS	10 µL, Non-Sterile	96/rk, 10 rks/cs	\$36.51
PR-10XLRK-NS	10 µL XL, Non-Sterile	96/rk, 10 rks/cs	\$36.51
PR-200RK-NS	200 µL, Non-Sterile	96/rk, 10 rks/cs	\$38.60
PR-200XLRK-NS	200 µL XL 5 in. Long, Non-Sterile	96/rk, 8 rks/cs	\$55.27
PR-300RK-NS	300 µL, Non-Sterile	96/rk, 10 rks/cs	\$48.29
PR-1000RK-NS	1000 µL, Non-Sterile	96/rk, 8 rks/cs	\$48.55
PR-1250RK-NS	1250 µL, Non-Sterile	96/rk, 8 rks/cs	\$43.04




PR1MA™ Sterile Racked Pipette Tips

ITEM #	DESCRIPTION	QTY	PRICE
PR-10RK-S	10 µL, Sterile	96/rk, 10 rks/cs	\$41.70
PR-10XLRK-S	10 µL XL, Sterile	96/rk, 10 rks/cs	\$41.70
PR-200RK-S	200 µL, Sterile	96/rk, 10 rks/cs	\$43.75
PR-200XLRK-S	200 µL XL 5 in. Long, Sterile	96/rk, 8 rks/cs	\$61.29
PR-300RK-S	300 µL, Sterile	96/rk, 10 rks/cs	\$53.82
PR-1000RK-S	1000 µL, Sterile	96/rk, 8 rks/cs	\$54.13
PR-1250RK-S	1250 µL, Sterile	96/rk, 8 rks/cs	\$46.27

Certified free of detectable RNase, DNase, DNA, PCR inhibitors and tested pyrogen-free. Sterile product is certified pyrogen-free.

LAB ESSENTIALS

- Compatible with Rainin® single- and multi-channel LiteTouch® pipettors
- A light force is required for adding or ejecting tips, making them easy to use
- LTS tips are made from medical-grade virgin polypropylene
- Racks and grids are made from our exclusive post-consumer recycled polypropylene, which reduces the use of virgin polypropylene by up to 61%



New!

PRIMA[®] LTS[®] Compatible Pipette Tips

ITEM #	DESCRIPTION	QTY	PRICE
PR-LTS20RFL	20 µL LTS-Style Tips, Reloads	96/pk, 10pks/cs	\$31.56
PR-LTS20RFL-S	20 µL LTS-Style Tips, Reloads, Presterile	96/pk, 10pks/cs	\$34.86
PR-LTS20RK-NS	20 µL LTS-Style Tips, Racked	96/pk, 10pks/cs	\$41.36
PR-LTS20RK-S	20 µL LTS-Style Tips, Racked, Sterile	96/pk, 10pks/cs	\$47.25
PR-LTS20RK-FL	20 µL LTS-Style Filter Tips, Racked, Sterile	96/pk, 10pks/cs	\$75.87
PR-LTS200RFL	200 µL LTS-Style Tips, Reloads	96/pk, 10pks/cs	\$30.12
PR-LTS200RFL-S	200 µL LTS-Style Tips, Reloads, Presterile	96/pk, 10pks/cs	\$33.43
PR-LTS200RK-NS	200 µL LTS-Style Tips, Racked	96/pk, 10pks/cs	\$43.73
PR-LTS200RK-S	200 µL LTS-Style Tips, Racked, Sterile	96/pk, 10pks/cs	\$49.57
PR-LTS200RK-FL	200 µL LTS-Style Filter Tips, Racked, Sterile	96/pk, 10pks/cs	\$75.87
PR-LTS300RFL	300 µL LTS-Style Tips, Reloads	96/pk, 10pks/cs	\$31.56
PR-LTS300RFL-S	300 µL LTS-Style Tips, Reloads, Presterile	96/pk, 10pks/cs	\$34.86
PR-LTS300RK-NS	300 µL LTS-Style Tips, Racked	96/pk, 10pks/cs	\$43.73
PR-LTS300RK-S	300 µL LTS-Style Tips, Racked, Sterile	96/pk, 10pks/cs	\$49.57
PR-LTS300RK-FL	300 µL LTS-Style Filter Tips, Racked, Sterile	96/pk, 10pks/cs	\$75.87
PR-LTS1000RFL	1,000 µL XL LTS-Style Tips, Reloads	96/pk, 8pks/cs	\$41.48
PR-LTS1000RFL-S	1,000 µL XL LTS-Style Tips, Reloads, Presterile	96/pk, 8pks/cs	\$48.06
PR-LTS1000RK-NS	1,000 µL XL LTS-Style Tips, Racked	96/pk, 8pks/cs	\$60.96
PR-LTS1000RK-S	1,000 µL XL LTS-Style Tips, Racked, Sterile	96/pk, 8pks/cs	\$65.52
PR-LTS1000RK-FL	1,000 µL XL LTS-Style Filter Tips, Racked, Sterile	96/pk, 8pks/cs	\$62.78

Certified free of detectable RNase, DNase, DNA, PCR inhibitors, and tested pyrogen-free

Rainin, LTS and LiteTouch are registered trademarks of Mettler-Toledo Rainin, LLC.

LAB ESSENTIALS



Aurelia™ Amazing®

ITEM #	SIZE	QTY	PRICE
92885	X-Small	300/bx	\$23.00
92886	Small	300/bx	\$23.00
92887	Medium	300/bx	\$23.00
92888	Large	300/bx	\$23.00
92889	X-Large	270/bx	\$23.00

Non-Sterile, Powder-Free, Ambidextrous, Finger Textured, Transblue Color, 2 mil Palm Thickness, 100% Latex-Free



Aurelia™ Quest®

ITEM #	SIZE	QTY	PRICE
92895	X-Small	100/bx	\$8.65
92896	Small	100/bx	\$8.65
92897	Medium	100/bx	\$8.65
92898	Large	100/bx	\$8.65
92899	X-Large	100/bx	\$8.65

Non-Sterile, Powder-Free, Ambidextrous, Finger Textured, White Color, 3.5 mil Palm Thickness, 100% Latex-Free



Aurelia™ Blush®

ITEM #	SIZE	QTY	PRICE
78885	X-Small	200/bx	\$17.31
78886	Small	200/bx	\$17.31
78887	Medium	200/bx	\$17.31
78888	Large	200/bx	\$17.31
78889	X-Large	200/bx	\$17.31

Non-Sterile, Powder-Free, Ambidextrous, Finger Textured, Pink Color, 2.5 mil Palm Thickness, 100% Latex-Free



Aurelia™ Protégé®

ITEM #	SIZE	QTY	PRICE
93995	X-Small	100/bx	\$8.65
93996	Small	100/bx	\$8.65
93997	Medium	100/bx	\$8.65
93998	Large	100/bx	\$8.65
93999	X-Large	100/bx	\$8.65

Non-Sterile, Powder-Free, Ambidextrous, Finger Textured, Sky Blue Color, 4.0 mil Palm Thickness, 100% Latex-Free



Aurelia™ Perform®

ITEM #	SIZE	QTY	PRICE
92395	X-Small	200/bx	\$17.31
92396	Small	200/bx	\$17.31
92397	Medium	200/bx	\$17.31
92398	Large	200/bx	\$17.31
92399	X-Large	200/bx	\$17.31

Non-Sterile, Powder-Free, Ambidextrous, Finger Textured, Teal Color, 2.5 mil Palm Thickness, 100% Latex-Free



Aurelia™ Bold®

ITEM #	SIZE	QTY	PRICE
73995	X-Small	100/bx	\$9.40
73996	Small	100/bx	\$9.40
73997	Medium	100/bx	\$9.40
73998	Large	100/bx	\$9.40
73999	X-Large	100/bx	\$9.40
73994	XX-Large	100/bx	\$9.40

Non-Sterile, Powder-Free, Ambidextrous, Finger Textured, Black Color, 5 mil Palm Thickness, 100% Latex-Free

Take advantage of starter kit savings!

- Four AlphaPettes™ (10 µL, 20 µL, 200 µL, 1000µL)
- Three racks of **PR1MA** Low-Binding Pipette Tips (1 each 10 µL, 200 µL, 1250 µL tips)
- One backpack
- One 6-place carousel pipettor stand



PR1MA™ Starter Kit

ITEM #	DESCRIPTION	PRICE
PR-LQH01-SK	PR1MA™ Starter Kit, Includes, Backpack, A-10, A-20, A-200, A-1000, 3rks PR1MA Tips, 6-pl Carousel*	\$997.00

*No mix and match or switching pipettor or tip sizes.

PR1MA™ Low Retention Graduated Pipette Tips are certified free of detectable RNase, DNase, DNA, PCR inhibitors, and tested pyrogen-free.



Family owned and trusted since 1983



Stocked products ship in 24 hours



Responses within 24 hours



REAL people answering your questions and placing your orders

We're here to help.
Place your order
with us today!