



Midsci™

EQUIPMENT
EQUIPMENT
EQUIPMENT
EQUIPMENT
EQUIPMENT

GUIDE 2024

**ORDER
TODAY!**

www.midsci.com • custserv@midsci.com
• 800.227.9997 • f: 636.225.9998



TABLE OF CONTENTS

3 •• Liquid Handling	25 •• PCR Equipment
8 •• Mixing	27 •• Electrophoresis Equipment
13 •• Centrifuges	28 •• Cold Storage
19 •• Centrifuge Tubes	33 •• Water Systems
20 •• Incubators	34 •• Balances
22 •• Live Cell Imaging	35 •• Gloves
23 •• Tissue Culture Supplies	36 •• PR1MA Starter Kit

Greetings to our valued customers,

Over the last 40 years MIDSCI™ has established a reputation of being a “go-to” for scientific supplies and equipment. We realize people are our most important asset. We make certain to take care of our MIDSCI™ family, so we can ensure we all work together to take care of you.

We believe in providing the best value to our customers by offering quality products and exceptional service. Our knowledge and expertise have allowed us to earn the trust of our customers daily, and we take pride in maintaining that trust.

In addition, we still have real people answering our phones and providing you with real answers. We not only value our relationship with you but take the responsibility of being the supplier you can count on.

To our long-time customers, we would like to thank you for your incredible support and your belief and trust in us.

To our new customers, we look forward to earning your trust and business well into the future.

Sincerely,
Mike Degenhart
President



If you can't find it here, give us a call at 800-227-9997 or email us at tech@midsci.com



MIDSCI assumes no liability for printing errors, mistakes, or product changes. All prices subject to change without notice.

**All promotions valid through 6/30/24, unless otherwise noted. Promos cannot be combined with any other promotions or discounts.

All terms and conditions apply. Visit our website for more information: <https://midsci.com/terms>

V.04.2024

LIQUID HANDLING


AlphaPette™ by PR1MA™ Single-Channel Pipettors

ITEM #	DESCRIPTION	QTY	PRICE
A-2	0.1 to 2 µL	1/ea	\$273.09
A-10	0.5 to 10 µL	1/ea	\$273.09
A-20	2 to 20 µL	1/ea	\$273.09
A-100	10 to 100 µL	1/ea	\$273.09
A-200	20 to 200 µL	1/ea	\$273.09
A-1000	100 to 1,000 µL	1/ea	\$273.09
A-5000	1,000 to 5,000 µL	1/ea	\$273.09
P4405	Multi-Channel Pipette Racks, 5-Place	1/ea	\$90.72
MID4280-CR	Carousel Pipettor Stand, 6-Place	1/ea	\$98.45



- Universal shaft will accept most pipette tips
- Lower ejection force
- Autoclavable
- Distinguishing color sets it apart for dedicated use
- Lifetime warranty (except calibration and abuse)

- Universal shaft will accept most pipette tips
- Lower ejection force
- 3-year warranty



Alphapette™ by PR1MA™ Multi-Channel Pipettors

ITEM #	DESCRIPTION	QTY	PRICE
AP4508-10	8-Channel, 1 to 10 µL	1/ea	\$710.01
AP4508-50	8-Channel, 5 to 50 µL	1/ea	\$710.01
AP4508-200	8-Channel, 20 to 200 µL	1/ea	\$710.01
AP4508-300	8-Channel, 50 to 300 µL	1/ea	\$710.01
AP4512-10	12-Channel, 1 to 10 µL	1/ea	\$864.95
AP4512-50	12-Channel, 5 to 50 µL	1/ea	\$864.95
AP4512-200	12-Channel, 20 to 200 µL	1/ea	\$864.95
AP4512-300	12-Channel, 50 to 300 µL	1/ea	\$864.95

- Buy one that lasts; 2-year warranty
- 8-hour rechargeable battery life
- QuickStand™ retractable kickstand

PR1MA™ Pipette Controller

ITEM #	DESCRIPTION	QTY	PRICE
MPIP-1	PR1MA™ Pipette Controller	1/ea	\$317.90
MID6950	PR1MA™ Reach Long Neck	1/ea	\$369.50



LIQUID HANDLING

PR1MA™ Low-Retention Pipette Tips



- All PR1MA™ tips are low-retention!
- Don't lose your sample in the tip
- Easier identification with color coded refills & clear racks
- Fits all standard pipettors
- Ideal for anyone that wants more consistent & reproducible data



PR1MA™ Bulk Pipette Tips

ITEM #	DESCRIPTION	QTY	PRICE
PR-10BK	10 µL, Non-Sterile	1,000/pk	\$17.20
PR-10XLBK	10 µL XL, Non-Sterile	1,000/pk	\$17.20
PR-200BK	200 µL, Non-Sterile	1,000/pk	\$13.37
PR-200XLBK	200 µL XL 5" Long, Non-Sterile	1,000/pk	\$43.88
PR-300BK	300 µL, Non-Sterile	1,000/pk	\$21.56
PR-1000BK	1000 µL, Non-Sterile	1,000/pk	\$22.17
PR-1250BK	1250 µL, Non-Sterile	1,000/pk	\$21.11



LIQUID HANDLING



PR1MA™ Pipette Tip Reloads



ITEM #	DESCRIPTION	QTY	PRICE
PR-10RFL	10 µL, Non-Sterile	96/rk, 10 rks/cs	\$27.85
PR-10XLRFL	10 µL XL, Non-Sterile	96/rk, 10 rks/cs	\$27.85
PR-200RFL	200 µL, Non-Sterile	96/rk, 10 rks/cs	\$26.59
PR-300RFL	300 µL, Non-Sterile	96/rk, 10 rks/cs	\$29.58
PR-1250RFL	1250 µL, Non-Sterile	96/rk, 8 rks/cs	\$36.61

Racks only, no boxes included



PR1MA™ Non-Sterile Racked Pipette Tips

ITEM #	DESCRIPTION	QTY	PRICE
PR-10RK-NS	10 µL, Non-Sterile	96/rk, 10 rks/cs	\$36.51
PR-10XLRK-NS	10 µL XL, Non-Sterile	96/rk, 10 rks/cs	\$36.51
PR-200RK-NS	200 µL, Non-Sterile	96/rk, 10 rks/cs	\$38.60
PR-200XLRK-NS	200 µL XL 5 in. Long, Non-Sterile	96/rk, 8 rks/cs	\$55.27
PR-300RK-NS	300 µL, Non-Sterile	96/rk, 10 rks/cs	\$48.29
PR-1000RK-NS	1000 µL, Non-Sterile	96/rk, 8 rks/cs	\$48.55
PR-1250RK-NS	1250 µL, Non-Sterile	96/rk, 8 rks/cs	\$43.04

PR1MA™ Sterile Racked Pipette Tips



ITEM #	DESCRIPTION	QTY	PRICE
PR-10RK-S	10 µL, Sterile	96/rk, 10 rks/cs	\$41.70
PR-10XLRK-S	10 µL XL, Sterile	96/rk, 10 rks/cs	\$41.70
PR-200RK-S	200 µL, Sterile	96/rk, 10 rks/cs	\$43.75
PR-200XLRK-S	200 µL XL 5 in. Long, Sterile	96/rk, 8 rks/cs	\$61.29
PR-300RK-S	300 µL, Sterile	96/rk, 10 rks/cs	\$53.82
PR-1000RK-S	1000 µL, Sterile	96/rk, 8 rks/cs	\$54.13
PR-1250RK-S	1250 µL, Sterile	96/rk, 8 rks/cs	\$46.27

LIQUID HANDLING




Supreme accuracy with the low-retention feature!

- Barrier tips have extra working volume, so liquid doesn't touch the filter
- Great stackable racks
- 10 µL XL tips have 1 µL graduations
- Clear & readable graduations
- Sterile, RNase free and DNase free

PR1MA™ Low-Retention Barrier Pipette Tips

ITEM #	DESCRIPTION	QTY	PRICE
PR-10RK-FL	10 µL, Sterile	96/rk, 10 rks/cs	\$67.16
PR-10XLRK-FL	10 µL XL, Sterile	96/rk, 10 rks/cs	\$66.96
PR-20RK-FL	20 µL, Sterile	96/rk, 10 rks/cs	\$66.96
PR-100RK-FL	100 µL, Sterile	96/rk, 10 rks/cs	\$66.96
PR-200RK-FL	200 µL, Sterile	96/rk, 10 rks/cs	\$66.96
PR-200XLRK-FL	200 µL XL 5 in. Long, Sterile	96/rk, 8 rks/cs	\$84.52
PR-300RK-FL	300 µL, Sterile	96/rk, 10 rks/cs	\$66.96
PR-1000RK-FL	1,000 µL, Sterile	96/rk, 8 rks/cs	\$58.52
PR-1250RK-FL	1,250 µL, Sterile	96/rk, 8 rks/cs	\$55.42

Improve consistency and reproducibility with low-retention tips!




Graduations:

10 µL: 2 & 10 µL	200 µL: 10, 50, & 100 µL
10 µL XL: 1, 2, & 10 µL	300 µL XL: 10, 50, 100, & 200 µL
20 µL: 10 µL	1,000 µL: 100, 500, 1,000 µL
100 µL: 10, 50, & 100 µL	1,250 µL: 100, 500, & 1,000 µL

LIQUID HANDLING

- Compatible with Rainin® single- and multi-channel LiteTouch® pipettors
- A light force is required for adding or ejecting tips, making them easy to use
- LTS tips are made from medical-grade virgin polypropylene
- Racks and grids are made from our exclusive post-consumer recycled polypropylene, which reduces the use of virgin polypropylene by up to 61%



New!

PR1MA[™] LTS[®] Compatible Pipette Tips

ITEM #	DESCRIPTION	QTY	PRICE
PR-LTS20RFL	20 µL LTS-Style Tips, Reloads	96/pk, 10pks/cs	\$31.56
PR-LTS20RFL-S	20 µL LTS-Style Tips, Reloads, Presterile	96/pk, 10pks/cs	\$34.86
PR-LTS20RK-NS	20 µL LTS-Style Tips, Racked	96/pk, 10pks/cs	\$41.36
PR-LTS20RK-S	20 µL LTS-Style Tips, Racked, Sterile	96/pk, 10pks/cs	\$47.25
PR-LTS20RK-FL	20 µL LTS-Style Filter Tips, Racked, Sterile	96/pk, 10pks/cs	\$75.87
PR-LTS200RFL	200 µL LTS-Style Tips, Reloads	96/pk, 10pks/cs	\$30.12
PR-LTS200RFL-S	200 µL LTS-Style Tips, Reloads, Presterile	96/pk, 10pks/cs	\$33.42
PR-LTS200RK-NS	200 µL LTS-Style Tips, Racked	96/pk, 10pks/cs	\$43.73
PR-LTS200RK-S	200 µL LTS-Style Tips, Racked, Sterile	96/pk, 10pks/cs	\$49.57
PR-LTS200RK-FL	200 µL LTS-Style Filter Tips, Racked, Sterile	96/pk, 10pks/cs	\$75.87
PR-LTS300RFL	300 µL LTS-Style Tips, Reloads	96/pk, 10pks/cs	\$31.56
PR-LTS300RFL-S	300 µL LTS-Style Tips, Reloads, Presterile	96/pk, 10pks/cs	\$34.85
PR-LTS300RK-NS	300 µL LTS-Style Tips, Racked	96/pk, 10pks/cs	\$43.73
PR-LTS300RK-S	300 µL LTS-Style Tips, Racked, Sterile	96/pk, 10pks/cs	\$49.57
PR-LTS300RK-FL	300 µL LTS-Style Filter Tips, Racked, Sterile	96/pk, 10pks/cs	\$75.87
PR-LTS1000RFL	1,000 µL XL LTS-Style Tips, Reloads	96/pk, 8pks/cs	\$41.48
PR-LTS1000RFL-S	1,000 µL XL LTS-Style Tips, Reloads, Presterile	96/pk, 8pks/cs	\$48.06
PR-LTS1000RK-NS	1,000 µL XL LTS-Style Tips, Racked	96/pk, 8pks/cs	\$60.96
PR-LTS1000RK-S	1,000 µL XL LTS-Style Tips, Racked, Sterile	96/pk, 8pks/cs	\$65.52
PR-LTS1000RK-FL	1,000 µL XL LTS-Style Filter Tips, Racked, Sterile	96/pk, 8pks/cs	\$62.78

Certified free of detectable RNase, DNase, DNA, PCR inhibitors, and tested pyrogen-free

Rainin, LTS and LiteTouch are registered trademarks of Mettler-Toledo Rainin, LLC.

MIXING

- Powerful vortexing for tubes up to 50 mL
- Pressure activated operation
- "Mini sized", less than 4 inches wide
- 2,800 rpm
- 2-year warranty

PR1MA™ Mini Vortexer



ITEM #	DESCRIPTION	QTY	PRICE
VM1-B	Mini Vortexer, Blue, 4 mm Orbit, Max Speed 2,800 rpm	1/ea	\$229.09
VM1-G	Mini Vortexer, Green, 4 mm Orbit, Max Speed 2,800 rpm	1/ea	\$229.09
VM1-GR	Mini Vortexer, Gray, 4 mm Orbit, Max Speed 2,800 rpm	1/ea	\$229.09
VM1-P	Mini Vortexer, Purple, 4 mm Orbit, Max Speed 2,800 rpm	1/ea	\$229.09
VM1-R	Mini Vortexer, Red, 4 mm Orbit, Max Speed 2,800 rpm	1/ea	\$229.09



- Powerful motor for INSTANT vortexing
- Continuous or "touch" operation
- Q-Drive™ dynamic balancing system
- Wide variety of accessories
- 2-year warranty

PR1MA™ Vortex Mixers

ITEM #	DESCRIPTION	QTY	PRICE
VM-3200	Adjustable Vortex Mixer, 3 mm Orbit, Max Speed 3,200 rpm	1/ea	\$296.62
VM-COMBO	Combo Head Attachment, 1.5 mL, 0.5 mL Microtubes & 1 Microplate	1/ea	\$78.69
VM-FLAT	Flat Head Attachment, 3 in.	1/ea	\$43.88
VM-H15	Horizontal Head, 12 x 1.5 mL Microtubes	1/ea	\$111.76
VM-H150	Horizontal Head, 4 x 15 mL Microtubes	1/ea	\$111.76
VM-H500	Horizontal Head, 2 x 50 mL Microtubes	1/ea	\$111.76

- Blocks available from 96 x 0.2 mL to 6 x 50 mL, custom, or multi-purpose
- Speed 200 to 1,500 rpm
- Heating rate 5°C/min - Cooling rate 2.5°C/min
- 2-year warranty

Multi-Therm™ Incubated Vortexer

ITEM #	DESCRIPTION	QTY	PRICE
MTH-1500	Temp Controlled Vortexer Ambient +5°C to 100°C	1/ea	\$2,425.66
MTHC-1500	Temp Controlled Vortexer Ambient -20°C to 100°C	1/ea	\$3,186.99
MTB-500	Block, 50 mL x 6 Tubes*	1/ea	\$494.19
MTB-5MT	Block, 5 mL x 15 conical Tubes*	1/ea	\$494.19
MTB-02	Block, 0.2 mL x 96 Tubes or 1 PCR Plate*	1/ea	\$416.84
MTB-05	Block, 0.5 mL x 54 Tubes*	1/ea	\$350.82
MTB-15	Block, 1.5 mL x 35 Tubes*	1/ea	\$350.82

*Blocks sold separately



Additional Blocks Available - visit www.midsci.com for more options!

MIXING

- Adjustable speed & tilt for a wide range of applications
- Large 14 x 12 in. work surface
- Dual platform accessory doubles work surface
- Cold room and incubator safe
- 2-year warranty

PR1MA™ Variable Speed 2D Rocker

ITEM #	DESCRIPTION	QTY	PRICE
R2D-30	2D Rocker, 14 x 12 in. w/ Flat Mat	1/ea	\$892.81
BR2000-STACK	14 x 12 in. Stacking Platform w/ Flat Mat	1/ea	\$152.13



- Adjustable speed & tilt for a wide range of applications
- Cold room and incubator safe
- Insta-TILT™ technology for instant adjustment to the tilt angle
- Dual platform accessory doubles work surface
- 2-year warranty


PR1MA™ 3D Nutating Shaker

ITEM #	DESCRIPTION	QTY	PRICE
R3D-30	12 x 12 in. 3D Nutating Shaker, 2 to 30 rpm	1/ea	\$1,120.98
BR1000-STACK	12 x 12 in. Stacking Platform	1/ea	\$139.87
BR1000-DIMPLED	12 x 12 in. Dimpled Mat	1/ea	\$48.72

- Four stainless steel platform options
- Reciprocal motion for staining, washing blots, and general mixing
- Economical analog control with easy-turn knobs
- Safe for cold room and incubator use
- 1-year warranty

Pro 30 Reciprocal Shaker

ITEM #	DESCRIPTION	QTY	PRICE
S2030-RC	Reciprocal Shaker, 6 - 150 Strokes/min, +4° to 70°C, No Platform	1/ea	\$1,636.69
S2030-10	Platform for Four Microplates	1/ea	\$295.40
S2030-12	Platform, w/ Non-Slip Rubber Mat, 12 X 12 in.	1/ea	\$303.00
S2031-12	Platform, w/ Non-Slip Rubber Mat, 30 x 30 in.	1/ea	\$284.87



MIXING

- Digital speed control up to 1,500 rpm
- 3 mm orbit for thorough mixing in microplates
- Includes platform for up to 4 microplates
- Cold room and incubator safe
- Optional adapters for mixing microtubes
- 2-year warranty

PR1MA™ Orbital Microplate Shaker

ITEM #	DESCRIPTION	QTY	PRICE
BT1500	13.1 x 11.7 x 4.7 in. Shaker, Holds Up to 4 Microplates	1/ea	\$1,168.82
BT1500-A05	20 x 0.5 mL Optional Adapter	1/ea	\$95.58
BT1500-A1520	20 x 1.5/2.0 mL Optional Adapter	1/ea	\$95.58



- Gentle orbital motion for gel applications
- Horizontally circular 19 mm orbit for aeration & mixing
- Large 14 x 12 inch work surface
- Variable speed control from 3 to 70 rpm
- Incubator and cold room safe
- 2-year warranty



Orbi-Blotter™

ITEM #	DESCRIPTION	QTY	PRICE
BT30	14 x 12 in. Orbital Shaker, 3 x 2 L or 9 lbs Capacity	1/ea	\$1,014.59
BR2000-STACK	14 x 12 in. Stacking Platform	1/ea	\$152.13
BR2000-FLAT	14 x 12 in. Flat Mat	1/ea	\$52.53


- Optimized speed and tilt for gel blotting, washing & staining
- Tri-directional motion for thorough, gentle mixing
- Optional stacking platforms available
- Cold room and incubator safe
- 2-year warranty

BlotBoy™ Low Speed 3D Rockers

ITEM #	DESCRIPTION	QTY	PRICE
R3D-710	7.5 x 10 in. Low Speed 3D Rocker, 12 rpm	1/ea	\$485.30
R3D-710ST	7.5 x 10 in. Stacking Platform	1/ea	\$139.87
R3D-1212	12 x 12 in. Low Speed 3D Rocker, 12 rpm	1/ea	\$587.40
R3D-1212ST	12 x 12 in. Stacking Platform	1/ea	\$139.87



MIXING



- Adjustable speed & tilt for a wide range of applications
- Insta-TILT™ technology for instant adjustment to the tilt angle
- Dual platform accessory doubles work surface
- Cold room and incubator safe
- 2-year warranty

BenchWaver™ 3-D Rocker

ITEM #	DESCRIPTION	QTY	PRICE
R3D-5000	13.5 X 13 in 3D Rocker, 5 to 105 rpm	1/ea	\$1,548.97
R3D-5000-ST	13 in. Stacking Platform, 5.4" Clearance	1/ea	\$159.14
R3D-5000-DMP	13 in. Dimpled Mat	1/ea	\$105.88

- Complete 360° rotation
- Adjustable 0-90° mixing angle
- Cold room and incubator compatible
- Interchangeable rotisserie options
- 1-year warranty

Labnet Revolver Tube Mixer

ITEM #	DESCRIPTION	QTY	PRICE
H5600	36 x 1.5 mL Mixer, 18 rpm Fixed Speed	1/ea	\$481.96



- Ideal for all tube sizes 1.5 to 50 mL
- Includes grooved mat for tub
- Cold room and incubator safe
- No adapters or mounting hardware required
- 2-year warranty

Ideal for a variety of tube sizes



TubeRocker Mini

ITEM #	DESCRIPTION	QTY	PRICE
M2100	12 x 4.5 in. Rocker w/ Grooved Mat, Fixed 24 rpm, 115V	1/ea	\$427.07

MIXING

- Gentle orbital motion, ideal for gel staining
- 13 x 12 in. non-slip platform included. Clamps sold separately
- 30 to 300 rpm
- Capacity 5 x 1 L or 14 lbs
- 2-year warranty

PRIMA™ Orbi-Shaker (Large)

ITEM #	DESCRIPTION	QTY	PRICE
OS-3000	Orbital Shaker (Large), 3 x 2 L Flasks or 14 lbs Capacity	1/ea	\$1,901.20
OS-PLAT	MAGic Clamp™ Platform for OS-3000	1/ea	\$321.26
H1000-MR-CMB	MAGic Clamp™ Combination Clamp	1/ea	\$52.63



- High-capacity benchtop shaker with 18 x 18 in. platform
- Circular 19 mm orbit
- Recommended to use with MAGic Clamp™ platform
- 2-year warranty



Orbi-Shaker™ XL

ITEM #	DESCRIPTION	QTY	PRICE
OSXL-1818	18 x 18 in. Shaker w/ Touch Screen Display & Rubber Mat Platform	1/ea	\$3,780.10
IS-1010-PLAT	18 x 18 in. MAGic Clamp™ Universal Platform	1/ea	\$413.65
H1000-MR-CMB	Adjustable MAGic Clamp™, Accepts 125, 250, & 500 mL Flasks	1/ea	\$52.63

- Large white ceramic top plate
- Digital speed and temperature control
- Safety LED indicates hot surface
- Available with stir capacities up to 15 L
- 2-year warranty

Digital Hotplates, Stirrers and Hotplate Stirrers

ITEM #	DESCRIPTION	QTY	PRICE
H3710-H	10 x 10 in Digital Hotplate	1/ea	\$558.93
H3710-S	10 x 10 in Magnetic Stirrer	1/ea	\$558.93
H3710-HS	10 x 10 in Digital Hotplate Magnetic Stirrer	1/ea	\$778.18
H3770-H	7 x 7 in Digital Hotplate	1/ea	\$321.24
H3770-S	7 x 7 in Magnetic Stirrer	1/ea	\$321.24
H3770-HS	7 x 7 in Digital Hotplate Magnetic Stirrer	1/ea	\$488.11



CENTRIFUGES

- Eight position microtube rotor
- A second rotor for 0.2 mL PCR tubes and strips is also included
- Storage compartment for PCR rotor
- 2-year warranty


PR1MA™ 8 Mini Centrifuge

ITEM #	DESCRIPTION	QTY	PRICE
MF8-R	8 Mini, Blue, 8 x 1.5/2.0 mL or 32 x 0.2 mL Tubes or 4 x PCR, Max Speed 6,000 rpm	1/ea	\$346.90
MF8-B	8 Mini, Clear, 8 x 1.5/2.0 mL or 32 x 0.2 mL Tubes or 4 x PCR, Max Speed 6,000 rpm	1/ea	\$346.90
MF8-G	8 Mini, Green, 8 x 1.5/2.0 mL or 32 x 0.2 mL Tubes or 4 x PCR, Max Speed 6,000 rpm	1/ea	\$346.90
MF8-P	8 Mini, Purple, 8 x 1.5/2.0 mL or 32 x 0.2 mL Tubes or 4 x PCR, Max Speed 6,000 rpm	1/ea	\$346.90
MF8-C	8 Mini, Red, 8 x 1.5/2.0 mL or 32 x 0.2 mL Tubes or 4 x PCR, Max Speed 6,000 rpm	1/ea	\$346.90

Includes both microtube and PCR strip rotors



Storage compartment for PCR strip rotor



COMBI-Rotor™
accepts tubes and strips!

- Twice the capacity of traditional mini centrifuges
- Starts & stops with closing/opening of the lid
- Unique COMBI-Rotor™ is all that is required for running 12 microtubes and 4 PCR strips simultaneously
- 2-year warranty

PR1MA™ 12 Mini Microcentrifuge

ITEM #	DESCRIPTION	QTY	PRICE
MF12	12 Mini, Capacity 12 x 1.5 mL Tubes, Max Speed 5,500 rpm	1/ea	\$457.47

- Includes rotor for all common centrifuge tubes up to 5mL
- Digital control up to 5,500 rpm / 1,690 xg
- Conserves valuable bench space
- Includes adapters for common blood tubes (12 x 75 mm)
- 2-year warranty

MyFuge™ 5D

ITEM #	DESCRIPTION	QTY	PRICE
C2205	5 Micro w/ Combo Rotor, Capacity 4 x 5 mL & 4 x 1.5/2.0 mL, Max Speed 5,500 rpm	1/ea	\$595.00



When 2 mL is not enough and 15 mL is too much - try 5 mL!

CENTRIFUGES



- Quickly spins down droplets in PCR plates (less than 20 seconds)
- Accepts skirted, non-skirted, standard, and most microtiter PCR plates
- 50% smaller than traditional centrifuges
- Unique rotor design prevents sample spillage
- 2-year warranty

PR1MA™ Microplate Microcentrifuge

ITEM #	DESCRIPTION	QTY	PRICE
PCRCENT	Microplate Centrifuge w/ Rotor & Plate Carriers, Capacity 2 x 2 PCR Plates or 24 x 0.2 mL PCR strips, Max Speed 2,550 rpm	1/ea	\$894.54
PCRCENT-A02	Adapter Plate for 96 x 0.2 mL Tubes or 12 PCR Strips	2/pk	\$67.26
PCRCENT-DWMP	Rotor for Deep Well Plates (Up to 35 mm Tall)	1/ea	\$231.34

- Unique rotor for 12-position PCR strips
- Starts and stops with opening/closing the lid
- Tiny benchtop footprint and affordability
- Place a StripSpin™ at the bench, thermal cycler & water bath for convenience
- 2-year warranty

New! StripSpin™ 12D

ITEM #	DESCRIPTION	QTY	PRICE
C2248	Mini 12-Strip PCR Tubes Centrifuge, Capacity 4 x 12 PCR Strips, Max Speed 5,500 rpm	1/ea	\$595.00



- Full color touch screen
- COMBI-Rotor™ accepts tubes and strips
- Compact, only 9 x 12 x 7.75 inches
- Variable speed control, set in RPM or RCF
- 2-year warranty



COMBI-Rotor™ accepts tubes and strips!

PR1MA™ MC24 Touch Microcentrifuge

ITEM #	DESCRIPTION	QTY	PRICE
MC24	High Speed Microcentrifuge w/ COMBI-Rotor™, Capacity Varies, Max Speed 13,500 rpm	1/ea	\$2,666.57

CENTRIFUGES

- 12 tube capacity / up to 16,100 x g
- Less than 9 inches wide
- Quick deceleration
- Whisper quiet operation
- 2-year warranty

PR1MA™ High Speed Microcentrifuge

ITEM #	DESCRIPTION	QTY	PRICE
MF1612	High Speed 16,100 x g, 12 x 1.5 mL Rotor Microcentrifuge	1/ea	\$1,690.33



- Brushless, maintenance-free motor
- Exceptionally quiet and cool operation
- High maximum speed/g-force
- 1-year warranty

Spectrafuge™ Microcentrifuge

ITEM #	DESCRIPTION	QTY	PRICE
C0160	Micro w/ Rotor, Gray, Capacity 18 x 1.5/2.0 mL, Max Speed 14,000 rpm	1/ea	\$1,484.79
C0160-R	Micro w/ Rotor, Red, Capacity 18 x 1.5/2.0 mL, Max Speed 14,000 rpm	1/ea	\$1,484.79
C0160-P	Micro w/ Rotor, Purple, Capacity 18 x 1.5/2.0 mL, Max Speed 14,000 rpm	1/ea	\$1,484.79
C0160-T	Micro w/ Rotor, Teal, Capacity 18 x 1.5/2.0 mL, Max Speed 14,000 rpm	1/ea	\$1,484.79



- Powerful refrigeration system, Quick-COOL to 4°C in 8 minutes
- Includes easy access rotor (24 x 1.5 mL) up to 17,135 x g
- Exceptionally quiet and compact
- Optional StripSpin™ adapter for 0.2 mL tubes and strips
- 2-year warranty

Best Seller!

Prism™ R Refrigerated Microcentrifuge

ITEM #	DESCRIPTION	QTY	PRICE	PROMO PRICE
C2500-R	Refrigerated Microcentrifuge, Capacity 24 x 1.5/2.0 mL Tubes, Max Speed 13,500 rpm	1/ea	\$5,524.45	\$3,498.00



Offer valid 04/01/24 - 06/30/24. This promo cannot be combined with other promotions or discounts. All terms apply.

CENTRIFUGES



MC-24R Refrigerated Microcentrifuge

ITEM #	DESCRIPTION	QTY	PRICE
C2417-R	Refrigerated High-Speed Microfuge w/ COMBI-rotor	1/ea	\$5,695.99
C2417-R-NR	Refrigerated High-Speed Microfuge, 120V	1/ea	\$5,395.00
C2417-2420H	Hermetically Sealed Rotor	1/ea	\$841.40
C2417-3202	PCR Strip Rotor	1/ea	\$795.00
C2417-RCMB	COMBI-Rotor	1/ea	\$565.00

HERMLE Z216-MK Refrigerated Microcentrifuge

ITEM #	DESCRIPTION	QTY	PROMO PRICE
Z216-K-CMB	High Speed w/ Combi 24 x 1.5 & 4 PCR Strip Rotor, Capacity 44 x 1.5/2.0 mL, Max Speed 15,000 rpm	1/ea	\$6,210.00
Z216-K-BND44	High Speed w/ 44 x 1.5 mL Rotor, Capacity 44 x 1.5/2.0 mL, Max Speed 15,000 rpm	1/ea	\$5,760.00

** For more rotors and accessories visit midsci.com **

Hermle Z216-K-BND44 Bundle!

Bundle Includes Z216-MK with 44 place 1.5/2.0 mL rotor

Hermle Z216-K-CMB Bundle!

Bundle Includes Z216-MK with COMBI-Rotor™ 24 x 1.5/2.0 and 4 x PCR strips

Promo valid through June 30, 2024. This promo cannot be combined with other promotions or discounts. All terms apply.



24 x 1.5/2.0 mL Rotor
Z216-4420



COMBI-Rotor™
Z216-2420S

Hermle Z287-BND-MB Bundle!

Bundle Includes Z287-A Micro bundle with 2 place microplate rotor and 18 place 1.5/2.0 mL rotor

Promo valid through June 30, 2024. All terms apply.

HERMLE Z287-A Microcentrifuge PLUS

ITEM #	DESCRIPTION	QTY	PRICE
Z287-A	Multi-purpose PLUS, Capacity 6 x 50 mL, Max Speed 14,000 rpm	1/ea	\$3,772.29
Z287-02MP	3,800 rpm, Swing-Out Rotor (90°), 2 x Microplate	1/ea	\$2,030.70
Z287-0450H	5,000 rpm, Swing-Out Rotor (90°), 4 x 50 mL	1/ea	\$1,550.12
Z287-1820	14,000 rpm, Angle Rotor (45°), 18 x 1.5/2.0 mL	1/ea	\$708.97
Z287-2420	12,000 rpm, Angle Rotor (45°), 24 x 1.5/2.0 mL	1/ea	\$781.04
Z287-BND-MB	Bundle Items Z287-A, Z287-02MP, Z287-1820	1/ea	\$4,500.00

** Additional rotors and accessories available online! **



Z287-02MP



Z287-2420



Z287-0450H

CENTRIFUGES

- Digital speed control up to 5,000 rpm
- Safety lid lock, prevents opening lid until rotor stops
- Unique airflow design keeps samples cool
- 8 tube capacity, for tubes up to 15 mL
- 2-year warranty



LC-8 & LC-8 Plus Compact Centrifuges

ITEM #	DESCRIPTION	QTY	PRICE
MSLD8	LC-8 w/ Rotor, Capacity 8 x 15 mL, Max Speed 3,500 rpm	1/ea	\$551.98
MSLD8P	LC-8 Plus w/ Rotor, Capacity 8 x 15 mL, Max Speed 5,000 rpm	1/ea	\$797.84



- Compact only 11 inches wide
- Unique EZ-Scroll™ touch pad
- 5-year warranty

Hermle Z207A Compact Centrifuge

ITEM #	DESCRIPTION	QTY	PRICE
Z207-A-CMB	Compact w/ Combination Rotor, Capacity 4 x 15 mL/50 mL, Max Speed 6,800 rpm	1/ea	\$2,368.70

- Swing out plate rotor accepts deepwell plates or can house up to 3 standard plates on each arm
- High speed microtube rotors (up to 18,624 x g)
- Fixed angle rotors for 0.2 mL to 50 mL

Hermle Z307-BND-MP Bundle!

Bundle Includes Z307 Universal Centrifuge, 120V
2 place swing out microplate rotor

Promo valid through June 30, 2024. All terms apply.

HERMLE Z307 Universal Centrifuge

ITEM #	DESCRIPTION	QTY	PROMO PRICE
Z307-BND-MP	Universal Centrifuge w/ Fixed Angle Rotor, Capacity 16 x 15 mL or 8 x 50 mL, Max Speed 14,000 rpm	1/ea	\$5,130.00



CENTRIFUGES

Hermle Z327 Universal Centrifuges

ITEM #	DESCRIPTION	QTY	PRICE
Z327	Universal, Capacity 4 x 200 mL, 16 x 15 mL or 8 x 50 mL, Max Speed 18,000 rpm	1/ea	\$6,653.71
Z327-K	Refrigerated Universal, Capacity 4 x 200 mL, 16 x 15 mL or 8 x 50 mL, Max Speed 18,000 rpm	1/ea	\$9,738.87
Z327-200	Swing-Out Rotor w/ 4 Buckets, 4 x 200 mL (4,500 rpm / 3,350 x g), Requires Inserts	1/ea	\$2,912.00
Z326-02MP	Swing-Out Rotor, 2 x 3 Microplates (4,500 rpm / 2,716 x g)	1/ea	\$2,826.20
Z326-3015	Rotor, 30 x 15 mL (4,500 rpm / 2,830 x g)	1/ea	\$2,084.78
Z326-1215	Rotor, 12 x 15 mL (6,000 rpm / 4,427 x g)	1/ea	\$1,348.26
Z326-0650-LS	Rotor, 6 x 50 mL (6,000 rpm / 4,427 x g)	1/ea	\$1,483.98
Z326-0685-25	Super Rotor, 6 x 85 mL High Speed Rotor (25°), (Z327: 11,000 rpm / 13,932 x g) (-K: 13,500 rpm/20,984 x g)	1/ea	\$2,999.26
Z326-1220H	Rotor, 12 x 1.5/2.0 mL (18,000 rpm/23,542 x g)	1/ea	\$1,789.22



Looking for Hermle bundles? Check out our website for more options!



- Swing-out rotor, capacity up to 16 x 50 mL and 40 x 15 mL
- Advanced programming with auto-rotor recognition
- 5-year warranty

Hermle Z366 Universal Centrifuges

ITEM #	DESCRIPTION	QTY	PRICE
Z366	Universal, Capacity 4 x 250 mL, Max Speed 16,000 rpm	1/ea	\$7,349.99
Z366-K	Refrigerated Universal, Capacity 4 x 250 mL, Max Speed 16,000 rpm, Temp -20°C to 40°C	1/ea	\$11,120.82

Hermle Z496-BND-BND Bundle!

Bundle Includes Centrifuge (230V ONLY), 4 x 1000 mL swing out rotor w/ rectangular buckets, 10 x 50 mL conical insert, & 24 x 15 mL conical insert

Promo valid through June 30, 2024. All terms apply.

Hermle Z496 Ultra High Capacity Centrifuge

ITEM #	DESCRIPTION	QTY	PRICE
Z496-K-BND	4 x 1,000 mL Swing-Out Rotor w/ Rectangular Buckets, 10 x 50 mL Conical Insert 2/pk, & 24 x 15 mL Conical Insert 2/pk	1/ea	\$24,300.00



**** For a complete listing of available rotors & accessories visit midsci.com ****

CENTRIFUGE TUBES



- G-Force: rated to 20,000 x g
- Graduations at every 100 μ L for easy volume estimation
- Low Retention: available with proprietary NoStick[®] resin to ensure maximum sample recovery
- Temperature Range: -80°C to +121°C



PR1MA[™] Premium Microcentrifuge Tubes

ITEM #	DESCRIPTION	QTY	PRICE
PR-MCT6	0.5 mL, Snap Cap, Boil Proof, Clear, Non-Sterile	1,000/pk	\$32.59
PR-MCT6-LR	0.5 mL, Low-Retention, Snap Cap, Boil Proof, Clear, Non-Sterile	500/pk	\$24.13
PR-MCT15	1.5 mL, Snap Cap, Boil Proof, Clear, Non-Sterile	500/pk	\$19.32
PR-MCT15-LR	1.5 mL, Low-Retention, Snap Cap, Boil Proof, Clear, Non-Sterile	250/pk	\$13.78
PR-MCT17	1.5 mL, Large Snap Cap, Boil Resist, Clear, Non-Sterile	500/pk	\$19.32
PR-MCT2	2 mL, Snap Cap, Boil Proof, Clear, Non-Sterile	500/pk	\$21.84

For colors add:

-B (blue), -G (green), -O (orange), -R (red), -P (purple), -Y (yellow), -RB (assorted), or -RBF (assorted neon) to the Item #
Sterile available for clear tubes only. Low-retention available for clear, non-sterile tubes only.



- DNase/RNase and pyrogen-FREE
- Black graduations
- White shaded writing area
- Temp Range: -80°C to +120°C

PR1MA[™] Conical Centrifuge Tubes

ITEM #	DESCRIPTION	QTY	PRICE
C15B	15 mL [PP], Screw Cap, Bagged, Sterile	500/cs	\$155.82
C15R	15 mL [PP], Screw Cap, Racked, Sterile	500/cs	\$142.62
C50B	50 mL [PP], Screw Cap, Bagged, Sterile	500/cs	\$183.52
C50R	50 mL [PP], Screw Cap, Racked, Sterile	500/cs	\$176.27

The 15 mL tubes can withstand up to 12,000 (RCF) g, and the 50 mL can withstand up to 16,000 (RCF) g.

- Outstanding for high centrifugation or ultra low freezing
- Liquid & temp resistant label on screw cap
- Leak proof
- White graduations down to 0.1 mL & 0.5 mL respectively
- Specs: -190°C to 100°C & 15,500 xg [PP]
-30°C to +70°C & 1,700 xg [PS], Sterile

TPP[®] Centrifuge Tubes

ITEM #	DESCRIPTION	QTY	PRICE
TP91015	15 mL, Conical, [PP], Sterile	40/pk, 20 pks/cs	\$653.83
TP91115	15 mL, Conical, [PS], Sterile	40/pk, 20 pks/cs	\$827.31
TP91050	50 mL, Conical, [PP], Sterile	20/pk, 18 pks/cs	\$351.23
TP91051	50 mL, Skirted Conical, [PP], Sterile	20/pk, 16 pks/cs	\$339.49
TP91550	50 mL, Conical, [PP], Sterile, Racked	20/rk, 10 rks/cs	\$413.95



INCUBATORS



- Stainless steel side walls with side capture zones and negative pressure to optimize containment
- Single-piece work tray, spill retaining and easy to clean
- Available in 8" sash height opening
- ULPA Filter, 10x filtration efficiency compared to the HEPA filter

NEW!

Esco® Airstream® Biosafety Cabinet

ITEM #	DESCRIPTION	QTY	PRICE
AR2-3S9-G4	Class II Type A2, 3 ft/0.9m, 8-inch Opening	1/ea	Contact for Quote
AR2-4S9-G4	Class II Type A2, 4 ft/1.2m, 8-inch Opening	1/ea	Contact for Quote
AR2-5S9-G4	Class II Type A2, 5 ft/1.5m, 8-inch Opening	1/ea	Contact for Quote
AR2-6S9-G4	Class II Type A2, 6 ft/1.8m, 8-inch Opening	1/ea	Contact for Quote

- Contamination control with long life ULPA filters and 90°C moist heat decontamination cycle.
- Direct heat & air jacket design with gentle yet effective forced convection.
- Intuitive interfaces and intelligent data and event logger.

CelCulture® CO₂ Incubator

ITEM #	DESCRIPTION	QTY	PRICE
CCL-170B-9-SP	170L CelCulture® CO ₂ Incubator, IR Sensor, CO ₂ Control, ULPA, High Temp Decontamination	1/ea	\$6,363.00



SureTherm™ CO₂ Incubators with IncuView™ LCI



ITEM #	DESCRIPTION	QTY	PRICE
H3565-45	45 L CO ₂ Incubator, w/ Two Shelves	1/ea	\$6,184.74
H3565-180	180 L CO ₂ Incubator, w/ Three Shelves	1/ea	\$7,470.80
H3565-180HD	180 L CO ₂ Incubator, w/ High Heat Decontamination	1/ea	\$9,339.94
H3565-180HDO2	180 L CO ₂ Incubator, w/ High Heat Decontamination, Split Window Door and O ₂ Control	1/ea	\$12,844.62

20 ORDER TODAY!

midsci.com • custserv@midsci.com • 800.227.9997 • f: 636.225.9998

INCUBATORS

- Digital control of temperature and CO₂
- Personal sized, 13 x 14.5 in. footprint (20 L)
- Temp range ambient -10°C to 55°C
- Accepts bottles and flasks up to 2 L
- Includes 2 stainless steel shelves
- Dual beam, infrared sensor
- 2-year warranty

PR1MA™ Mini CO₂ Incubator

ITEM #	DESCRIPTION	QTY	PRICE
INCO2-H	0.8 CuFt Incubator	1/ea	\$4,518.82
H2300-REG	Optional CO ₂ Gas Regulator	1/ea	\$385.53



- Innovative chamber design features smooth rounded internal corners for easy cleaning
- Infra-Red(IR) CO₂ sensor for precise and reliable CO₂ detection
- Stackable design up to 3 units
- 25 mm and 50 mm orbit models available

AEOLUS Compact CO₂ Incubator Shakers

ITEM #	DESCRIPTION	QTY	PRICE
IS-AMA	Incubator Shaker, 7" High-Res Touch Screen, Air Jacket, 25 mm (1.02 in.) Orbit, Dims 970 x 780 x 530 mm	1/ea	\$17,550.00
IS-A5MA	Incubator Shaker, 7" High-Res Touch Screen, Air Jacket, 50 mm (2 in.) Orbit, Dims 970 x 780 x 530 mm	1/ea	\$19,500.00
IS-ACMA	CO ₂ Incubator Shaker, 7" High-Res Touch Screen, Air Jacket, 25 mm (1.02 in.) Orbit, 970 x 780 x 530 mm	1/ea	\$22,725.00
IS-A5CMA	CO ₂ Incubator Shaker, 7" High-Res Touch Screen, Air Jacket, 50 mm (2 in.) Orbit, 970 x 780 x 530 mm	1/ea	\$23,730.00

- 7" High-Res Color Touch Screen with intuitive graphic interface
- Capacity is 6 x 5 L
- Infra-Red(IR) CO₂ sensor for precise and reliable CO₂ detection
- Optional WiFi and RS485 monitoring and control
- Two Sets of UV lights for easy sterilization
- Stackable design up to 3 units.



AEOLUS Stackable 302 L CO₂ Incubator Shaker

ITEM #	DESCRIPTION	QTY	PRICE
IS-18A	302 L Stackable Incu Shaker, 25 mm Orbit	1/ea	\$23,010.00
IS-18CA	302 L Stackable CO ₂ Incu Shaker, 25 mm Orbit	1/ea	\$33,240.00
IS-18B5A	302 L Stackable Incu Shaker, 50 mm Orbit	1/ea	\$35,623.14
IS-18C5A	302 L Stackable CO ₂ Incu Shaker, 50 mm Orbit	1/ea	\$34,240.00

LIVE CELL IMAGING



- Real-time live cell monitoring of bottom layer for multi-layer flasks up to 10 levels.
- Checks cell quality and detects confluency rate while culturing
- Improved reproducibility for mass cell culture and production
- No need to take the cells out of the incubator thus eradicating contamination worries
- Multi-position time-lapse imaging & recording

Curiosis Celloger Stack

ITEM #	DESCRIPTION	QTY	PRICE
CRCLG-SB02	Real-Time Cell Monitoring System	1/ea	\$40,000.00

- Compact design easily fits into standard CO₂ incubators
- Position memory allows multiple area to be set for automatic observation, up to including 96 well plates
- Time-lapse image capturing & video clip
- Software is easy to use
- Provides thermo-hygrometer data in real time
- Allows for imaging of standard plates, dishes and flask.



Curiosis Celloger Mini

ITEM #	DESCRIPTION	QTY	PRICE
CRCLG-MB01	Mini Live Celloger, 4x / 10x	1/ea	\$19,500.00

- One color fluorescence (GFP or RFP) and brightfield imaging technology
- Compact size that easily fits into standard CO₂ incubators
- High quality system yet cost-effective
- Compatible with various cell and tissue culture vessel types
- Intuitive UI/UX and easy to acquire confluency data




Curiosis Celloger Nano

ITEM #	DESCRIPTION	QTY	PRICE
CRCLG-NBG04	Brightfield, Green Fluorescence, 4X	1/ea	\$11,000.00
CRCLG-NBG10	Brightfield, Green Fluorescence, 10X	1/ea	\$12,000.00
CRCLG-NBR04	Brightfield, Red Fluorescence, 4X	1/ea	\$11,000.00
CRCLG-NBR10	Brightfield, Red Fluorescence, 10X	1/ea	\$12,000.00
CRCLG-NB04	Brightfield, 4X	1/ea	\$9,500.00
CRCLG-NB10	Brightfield, 10X	1/ea	\$9,900.00

TISSUE CULTURE SUPPLIES

You can always trust MIDSCI™ for quality TPP® products that help you save money where it matters most, cost per experiment.

- TPP® is Swiss made with excellent flatness and clarity
- Flatness = Less clumping and even growth
- Clarity = Great for microscopy
- Vents = Rapid temperature equilibration for better cell consistency and health
- Stacked dishes, flasks, and plates interlock and create vents between levels
- Interlock makes handling stacks much less prone to dropping



TPP® Tissue Culture Plates

- Minimized edge effect in 96-well plates
- Easy to read yellow marking spot and black alphanumeric

ITEM #	DESCRIPTION	QTY	PRICE
TP92006	6-Well TC Plate [PS], Flat Bottom, Sterile	1/pk, 126 pks/cs	\$711.31
TP92012	12-Well TC Plate [PS], Flat Bottom, Sterile	1/pk, 126 pks/cs	\$711.96
TP92024	24-Well TC Plate [PS], Flat Bottom, Sterile	1/pk, 126 pks/cs	\$710.69
TP92048	48-Well TC Plate [PS], Flat Bottom, Sterile	1/pk, 126 pks/cs	\$711.32
TP92096	96-Well TC Plate [PS], Flat Bottom, Sterile	1/pk, 162 pks/cs	\$877.09
TP92097	96-Well TC Plate [PS], Round Bottom, Sterile	1/pk, 162 pks/cs	\$872.47


Also available in multipacks!

Contact your trusted MIDSCI™ Lab Consultant for more information

- Easy to read painted yellow marking spot
- Ridged grip around base greatly improves handling. No more dropped dishes

TPP® Tissue Culture Dishes

ITEM #	DESCRIPTION	QTY	PRICE
TP93040	40 x 11 mm TC Dish [PS], Sterile	20/pk, 45 pks/cs	\$701.57
TP93060	60 x 16 mm TC Dish [PS], Sterile	14/pk, 60 pks/cs	\$725.52
TP93100	96 x 21 mm TC Dish [PS], Sterile	10/pk, 24 pks/cs	\$393.88
TP93150	146 x 21 mm TC Dish [PS], Sterile	5/pk, 20 pks/cs	\$466.96



TISSUE CULTURE SUPPLIES

- Wider neck allows for better access for 100% cell retrieval when scraping
- Short access to cell floor reduces chance of contamination



TPP® Tissue Culture Flasks

ITEM #	DESCRIPTION	QTY	PRICE
TP90025	25 cm ² (92 x 51 x 29 mm) TC Flask [PS], Solid Vent Cap [PE], Sterile	10/pk, 36 pks/cs	\$794.00
TP90026	25 cm ² (92 x 51 x 29 mm) TC Flask [PS], Filter Cap [PE], Sterile	10/pk, 36 pks/cs	\$867.03
TP90075	75 cm ² (155 x 87 x 40 mm) TC Flask [PS], Solid Vent Cap [PE], Sterile	5/pk, 20 pks/cs	\$451.70
TP90076	75 cm ² (155 x 87 x 40 mm) TC Flask [PS], Filter Cap [PE], Sterile	5/pk, 20 pks/cs	\$485.57
TP90150	150 cm ² (210 x 122 x 50 mm) TC Flask [PS], Solid Vent Cap [PE], Sterile	3/pk, 13 pks/cs	\$364.88
TP90151	150 cm ² (210 x 122 x 50 mm) TC Flask [PS], Filter Cap [PE], Sterile	3/pk, 13 pks/cs	\$388.30
TP90300	300 cm ² (275 x 170 x 50 mm) TC Flask [PS], Solid Vent Cap [PE], Sterile	3/pk, 6 pks/cs	\$389.53
TP90301	300 cm ² (275 x 170 x 50 mm) TC Flask [PS], Filter Cap [PE], Sterile	3/pk, 6 pks/cs	\$408.34

Filter caps come with a 0.22 µm PTFE Membrane

- Ideal for suspension cultures
- Scale up or scale down ability to do large in parallel studies
- Better viability & proliferation than flasks or bottles
- Better gas exchange control & cyclonic stirring than flasks or bottles

TPP® TubeSpin® Bioreactor

ITEM #	DESCRIPTION	QTY	PRICE
TP87015	15 mL Conical [PP], Sterile	40/pk, 20 pks/cs	\$1,111.48
TP87017	15 mL Round [PP], Sterile	40/pk, 20 pks/cs	\$1,665.50
TP87017-200	15 mL Round [PP], Sterile	40/pk, 5 pks/cs	\$414.91
TP87017-400	15 mL Round [PP], Sterile	40/pk, 10 pks/cs	\$834.22
TP87050	50 mL Conical [PP], Sterile	20/pk, 9 pks/cs	\$484.26
TP87050-60	50 mL Conical [PP], Sterile	20/pk, 3 pks/cs	\$177.91
TP87056	50 mL Round [PP], Sterile	20/pk, 9 pks/cs	\$755.42
TP87450	450 mL Conical [PP], Sterile	1/pk, 32 pks/cs	\$643.86
TP87600	600 mL Conical [PP], Sterile	1/pk, 26 pks/cs	\$1,015.84
TP87600-13	600 mL Conical [PP], Sterile	1/pk, 13 pks/cs	\$506.13
TP87600-3	600 mL Conical [PP], Sterile	1/pk, 3 pks/cs	\$117.33



PCR EQUIPMENT

- Gradient capability
- Precision thermal control
- 6-segment block temperature control
- TM calculator
- 3-yr warranty

Best Seller!

MultiGene™ Optimax Thermal Cycler

ITEM #	DESCRIPTION	QTY	PRICE	PROMO PRICE
TC9610	96 Sample Thermal Cycler, 115V Model	1/ea	\$4,241.54	\$2,849.00



Offer valid 04/01/24 - 06/30/24. This promo cannot be combined with other promotions or discounts. All terms apply.



- Run multiple templates and programs at the same time!
- Dual or triple block models are available, both offering extremely fast cycling profiles
- Large touch screen interfaces
- Self-adjusting heated lids, gradient capabilities
- Easy-to-use program wizard
- 2-yr warranty

MultiCycler™ Multi-block Thermal Cycler

ITEM #	DESCRIPTION	QTY	PRICE
T3-0248G	MultiCycler™ 248, 2 x 48 x 0.2 mL Blocks	1/ea	\$8,868.90
T3-0332G	MultiCycler™ 332, 3 x 32 x 0.2 mL Blocks	1/ea	\$10,194.90

The Azure Cielo is a 96-well real-time PCR machine that features innovative and high-performance optical technology with broad application compatibility.

- Unparalleled data uniformity
- Data emailed after run completion
- Shorter run times
- Multiplex qPCR of up to six-targets

qPCR Azure Systems

ITEM #	DESCRIPTION	QTY	PRICE	PROMO PRICE
AIQ030	Cielo 3 qPCR System, 0.1 mL x 96	1/ea	\$27,850.00	
AIQ060	Cielo 6 qPCR System, 0.1 mL x 96*	1/ea	\$34,850.00	\$16,895
AIQ1000	Cielo 3 Software Plate Setup, PCR Compatible	1/ea	\$1,600.00	



SPECIAL PROMOTIONAL PRICE:

*Includes: Azure Cielo 6, software, and 18 months of warranty.
Offer valid 04/01/24 - 06/30/24

PCR EQUIPMENT

- Variable speed has mix/pause/rock function
- Includes tube holders for 11 x 1.5/2.0 mL, 5 x 10-15 mL, 5 x 5-7 mL, and 6 x 50 mL

- Light weight (3 lbs)
- 2-year warranty

**Need Temperature Control?
Contact us!**

PR1MA™ Roto-Mini Series

ITEM #	DESCRIPTION	QTY	PRICE
RFP2020	10.8 x 6 x 6.5 in., 24 rpm, Fixed Speed	1/ea	\$441.50
RMP2024	10.8 x 6 x 6.5 in., 5 to 70 rpm, Variable Speed	1/ea	\$530.71



- Gradient capabilities- ideal for protocol optimization
- Intuitive touchscreen control- Works even when wearing gloves!
- No need for interchangeable blocks - multi-format accepts:
 - 0.2 mL tubes or strips (96)
 - 0.5 mL tubes (39)
 - 96-well PCR plate

PR1MA™ TC9639 Thermal Cycler

ITEM #	DESCRIPTION	QTY	PRICE
BSTC-96	10.2 x 14.2 in. Cycler w/ Combination Block, 96 x 0.2 mL	1/ea	\$5,995.00
BSTC-384	10.2 x 14.2 in. Cycler w/ 384-Well Block	1/ea	\$6,979.04
BSTC-AIS	Combo Adapter (4 slides)	1/ea	\$587.20

- Complete system, includes gel box, power supply and casting sets
- Quick and easy gel casting, no tape, clamps or spacers required
- Compact design, footprint: 7 x 5.1 in.

PR1MA™ Mini Gel Electrophoresis System

ITEM #	DESCRIPTION	QTY	PRICE
EMINIGEL	Horizontal System, One 10.5 x 6 cm or two 5 x 6 cm Gels, Variable Teeth	1/ea	\$636.38
E1101-COMB1	22/12 Reversible Comb	2/pk	\$9.50
E1101-COMB2	18/10 Reversible Comb	2/pk	\$12.14
E1101-CS1	10.5 x 6 cm Gel Casting Stand	1/ea	\$49.67
E1101-CS2	5 x 6 cm Gel Casting Stand	1/ea	\$78.74



ELECTROPHORESIS EQUIPMENT



- Molded, durable construction
- Rubber gates for easy, tapeless gel casting
- Adjustable height combs
- Loading guides aid in well visualization



Enduro™ Horizontal Electrophoresis System

ITEM #	DESCRIPTION	QTY	PRICE
E1007-7	7 x 7 cm (Gel Tank), Buffer Vol 225 mL, Max Sample Cap 32	1/ea	\$609.47
E1007-10	7 x 10 cm (Gel Tank), Buffer Vol 225 mL, Max Sample Cap 64	1/ea	\$634.73
E1010-10	10 x 10 cm (Gel Tank), Buffer Vol 300 mL, Max Sample Cap 100	1/ea	\$669.04
E1015-10	15 x 10 cm (Gel Tank), Buffer Vol 500 mL, Max Sample Cap 140	1/ea	\$775.78
E1015-15	15 x 15 cm (Gel Tank), Buffer Vol 500 mL, Max Sample Cap 210	1/ea	\$1,003.29
E1020-20	20 x 20 cm (Gel Tank), Buffer Vol 1,200 mL, Max Sample Cap 450	1/ea	\$1,048.84

- Compact and simple with a 2-year warranty
- Choose constant 100 V or 200 V output , 200 mA
- Two pairs of output terminals

MYVOLT™ Mini Power Supply

ITEM #	DESCRIPTION	PRICE
PS2100	Mini Power Supply	\$212.42



- Touch screen control
- Bright color display
- Dual orientation, auto rotating display
- 2-year warranty

PR1MA™ Touch Power Supply

ITEM #	DESCRIPTION	QTY	PRICE
PS2300	Compact Electrophoresis Power Supply, 100 W	1/ea	\$566.32

Output Settings	Constant voltage or constant current
Voltage Output	10-300 V (1 V increments)
Current Output	4-400 mA (1 mA increments)
Max Wattage	100 W
Timer	1 min. to 20 hours or continuous
Output Terminals	4 pairs (wired in parallel)

Safety	No load detection, overload detection with power limiting function
Display	2 x 2 in. color touchscreen
Dimensions	8.5 x 6.3 x 4.0 inches
Weight	2 lbs
Electrical	100 V to 230 V, 50-60 Hz (Universal Input)

COLD STORAGE

- Microprocessor-controlled system designed for controllable range of -40°C to -86°C for cabinet space with 1°C increment
- Large LED display for cabinet temperature, set temperature, ambient temperature, and input voltage
- Automatic clean-filter alarm and sensor error alert

Haier -86°C ULT Freezers

ITEM #	DESCRIPTION	PRICE
DW-86L338J	11.9 CuFt, -40°C to -86°C, Upright, 115 V/60 HZ	Contact for Quote
DW-86L578J	20.4 CuFt, -40°C to -86°C, Upright, 115 V/60 HZ	Contact for Quote
DW-86L728J	25.7 CuFt, -40°C to -86°C, Upright, 115 V/60 HZ	Contact for Quote
DW-86L828JA	29.2 CuFt, -40°C to -86°C, Upright, 220 V/60 HZ	Contact for Quote

Energy Star® Certified



- Advanced innovative design delivers excellent energy savings
- Optional barcode scanner for simple, effortless and precise identification
- Standard equipped with key lock, padlock and electromagnetic lock with optional fingerprint lock, providing multiple safeguards for sample safety

Haier -80°C Smart Frequency ULT Freezers

ITEM #	DESCRIPTION	PRICE
DW-86L579BPT	20.5 CuFt, -86°C, Upright, 110-230 V/60 Hz	Contact for Quote
DW-86L729BPT	26 CuFt, -86°C, Upright, 110-230 V/60 Hz	Contact for Quote
DW-86L829BPT	29 CuFt, -86°C, Upright, 208-230 V/60 Hz	Contact for Quote

Energy Star® Certified

- Dual independent refrigeration systems for maximum sample safety
- High speed refrigeration system for faster pull down and temperature recovery after door opening
- World leading energy saving refrigeration technology
- Smart full-size touchscreen and wireless monitoring connectivity

Haier -86°C TwinCool ULT Freezers

ITEM #	DESCRIPTION	PRICE
DW-86L578SAT	20.4 CuFt, -40°C to -86°C, Upright, 220 V/60 HZ	Contact for Quote
DW-86L578ST	20.4 CuFt, -40°C to -86°C, Upright, 115 V/60 HZ	Contact for Quote
DW-86L728SAT	25.7 CuFt, -40°C to -86°C, Upright, 220 V/60 HZ	Contact for Quote
DW-86L728ST	25.7 CuFt, -40°C to -86°C, Upright, 115 V/60 HZ	Contact for Quote

Energy Star® Certified



For a limited time, Purchase any Haier ULT series with MIDSCI™
(TwinCool, Smart Frequency, Energy Savings Salvum J Series)

And get: One Shelf of Slider Racks and Dock to Dock Shipping for Contiguous USA For FREE!

Promo valid through June 30, 2024. This promo cannot be combined with other promotions or discounts. All terms apply.

COLD STORAGE

Organize your PCR tubes with our
PCR boxes!



PR1MA™ Cardboard Freezer Boxes

ITEM #	DESCRIPTION	PRICE
SB2CC-(xx)*	Standard White Cardboard 5.25 x 5.25 x 2 in Box with Divider	\$5.00

* Please specify (xx) as 16, 25, 49, 64, 81, 100, or 144 when ordering

Not sure what to buy?
Scan the code to check out
our freezer rack configurator



Freezer Inventory Racks

ITEM #	TYPE	FREEZER STYLE	D X H X W	BOX HEIGHT	RACK CAPACITY (DXH)	PRICE
UFE-332	Eco	Upright	16 5/16 x 6 5/8 x 5 1/2 in	2 in	3 x 3 = 9	\$68.39
UFE-342	Eco	Upright	16 5/16 x 8 7/8 x 5 1/2 in	2 in	3 x 4 = 12	\$78.42
UFE-352	Eco	Upright	16 5/16 x 11 x 5 1/2 in	2 in	3 x 5 = 15	\$92.48
UFE-362	Eco	Upright	16 5/16 x 12 5/8 x 5 1/2 in	2 in	3 x 6 = 18	\$119.74
UFE-432	Eco	Upright	21 5/16 x 6 5/8 x 5 1/2 in	2 in	4 x 3 = 12	\$83.02
UFE-442	Eco	Upright	21 5/16 x 8 7/8 x 5 1/2 in	2 in	4 x 4 = 16	\$122.88
UFE-452	Eco	Upright	21 5/8 x 11 x 5 1/2 in	2 in	4 x 5 = 20	\$139.60
UFE-462	Eco	Upright	21 5/8 x 12 5/8 x 5 1/2 in	2 in	4 x 6 = 24	\$153.62
UFE-532	Eco	Upright	26 3/4 x 6 5/8 x 5 1/2 in	2 in	5 x 3 = 15	\$107.81
UFE-542	Eco	Upright	26 3/4 x 8 7/8 x 5 1/2 in	2 in	5 x 4 = 20	\$170.64
UFE-552	Eco	Upright	26 3/4 x 11 x 5 1/2 in	2 in	5 x 5 = 25	\$175.07
UF-332	Std	Upright	16 5/16 x 6 5/8 x 5 1/2 in	2 in	3 x 3 = 9	\$75.97
UF-342	Std	Upright	16 5/16 x 6 5/8 x 5 1/2 in	2 in	3 x 4 = 12	\$97.52
UF-352	Std	Upright	16 5/16 x 11 x 5 1/2 in	2 in	3 x 5 = 15	\$115.56
UF-442	Std	Upright	21 5/8 x 8 7/8 x 5 1/2 in	2 in	4 x 4 = 16	\$108.97
UF-452	Std	Upright	21 5/8 x 11 x 5 1/2 in	2 in	4 x 5 = 20	\$141.09
UF-462	Std	Upright	21 5/8 x 12 5/8 x 5 1/2 in	2 in	4 x 6 = 24	\$192.04
UF-542	Std	Upright	26 3/4 x 8 7/8 x 5 1/2 in	3 in	5 x 4 = 20	\$155.04
UF-433	Std	Upright	21 5/8 x 9 3/8 x 5 1/2 in	2 in	4 x 3 = 12	\$96.27
CF-10-2	Std	Chest	22 x 5 5/8 x 5 1/2 in	2 in	1 x 10 = 10	\$76.98
CF-12-2	Std	Chest	26 3/8 x 5 5/8 x 5 1/2 in	2 in	1 x 12 = 12	\$92.24

More racks available at midsci.com.

ORDER TODAY!

midsci.com • custserv@midsci.com • 800.227.9997 • f: 636.225.9998 29

COLD STORAGE



- Cycle defrost
- Microprocessor temperature controller
- Password-protected control panel
- Touch screen digital temperature display and alarm module with battery back-up to support temperature monitoring and alarms

Templog Premier Solid Door Lab Refrigerator

ITEM #	DESCRIPTION	PRICE
ABT-HC-12S-TS	12 CuFt, 1°C to 10°C, 29.75 x 25 x 65.75 in	Contact for Quote
ABT-HC-16S-TS	16 CuFt, 1°C to 10°C, 29.75 x 25 x 79 in.	Contact for Quote
ABT-HC-23S-TS	23 CuFt, 1°C to 10°C, 34.75 x 26.875 x 81.75 in.	Contact for Quote
ABT-HC-26S-TS	26 CuFt, 1°C to 10°C, 36.50 x 28.375 x 78.25 in.	Contact for Quote
ABT-HC-36S-TS	36 CuFt, 1°C to 10°C, 34.75 x 39.625 x 81.75 in.	Contact for Quote

- Microprocessor temperature controller
- Cycle defrost
- Touch screen digital temperature display and alarm module
- Password-protected control panel

Templog Premier Glass Door Lab Refrigerator

ITEM #	DESCRIPTION	PRICE
ABT-HC-LP-12-TS	12 CuFt, 1°C to 10°C, 29.75 x 25 x 62.25 in.	Contact for Quote
ABT-HC-LP-16-TS	16 CuFt, 1°C to 10°C, 29.75 x 25 x 79 in.	Contact for Quote
ABT-HC-LP-23-TS	23 CuFt, 1°C to 10°C, 34.75 x 26.875 x 81.75 in.	Contact for Quote
ABT-HC-LP-26-TS	26 CuFt, 1°C to 10°C, 36.50 x 28.375 x 78.25 in.	Contact for Quote
ABT-HC-LP-33-TS	33 CuFt, 1°C to 10°C, 33.25 x 39.625 x 81.75 in.	Contact for Quote
ABT-HC-LP-36-TS	36 CuFt, 1°C to 10°C, 34.75 x 39.625 x 81.75 in.	Contact for Quote
ABT-HC-LP-47-TS	47 CuFt, 1°C to 10°C, 33.25 x 54.125 x 81.75 in.	Contact for Quote



- 2 Probes; 1 in air and 1 in glycol bottle
- Six fixed shelves
- Keyed door lock
- Magnetic door gasket for positive seal
- Includes monitoring system

-20°C Lab Freezers Standard Manual Defrost

ITEM #	DESCRIPTION	PRICE
ABT-HC-MFP-14	14 CuFt, -15°C to -25°C, w/ Natural Refrigerants	Contact for Quote
ABT-HC-MFP-17	17 CuFt, -15°C to -25°C, w/ Natural Refrigerants	Contact for Quote
ABT-HC-MFP-20	20 CuFt, -15°C to -25°C, w/ Natural Refrigerants	Contact for Quote

Don't see what you need?
Call us, we have many more options available for you!

COLD STORAGE



- Rugged and durable aluminium construction
- Vapour phase cryogenic storage
- Straw and cryovial storage options
- Designed to ensure no LN₂ spillage
- 3-year vacuum warranty

Haier LN₂ Dryshipper Series for Transportation

ITEM #	DESCRIPTION	PRICE
YDH-3	3 L, 132 Straws (0.5 mL), Includes Lockable Cap, Protective Jacket	Contact for Quote
YDH-6-80	6 L, 374 Straws (0.5 mL), Includes Lockable Cap, Protective Jacket	Contact for Quote
YDH-10-125-F	10 L, 100 Vials (2 mL), Includes Rack, Cryoboxes, Cap, Jacket	Contact for Quote
YDH-15-216-F	15 L, 300 Vials (2 mL), Includes Rack, Cryoboxes, Cap	Contact for Quote
YDH-25-216-F	25 L, 500 Vials (2 mL), Includes Rack, Cryoboxes, Cap	Contact for Quote

- Durable aluminum construction
- Temperature monitoring available
- Liquid or vapor phase option available
- High thermal efficiency
- 5-year vacuum warranty

Haier LN₂ Medium Sized Storage Series

ITEM #	DESCRIPTION	PRICE
YDS-65-216	65 L, 2,400 Vials (2 mL), Includes Lockable Cap	Contact for Quote
YDS-95-216	95 L, 3,000 Vials (2 mL), Includes Lockable Cap	Contact for Quote
YDS-115-216	115 L, 3,600 Vials (2 mL), Includes Lockable Cap	Contact for Quote
YDS-140-216	140 L, 4,800 Vials (2 mL), Includes Lockable Cap	Contact for Quote
YDS-175-216	175 L, 6,000 Vials (2 mL), Includes Lockable Cap	Contact for Quote



- LN₂ splash proof
- Data logging
- Hot gas bypass
- 5-year vacuum warranty



Haier Smart Type LN₂ Biobank Series

ITEM #	DESCRIPTION	PRICE
YDD-350-326	350 L, 13,000 Vials, w/ Vacuum Insulated Hose	Contact for Quote
YDD-550-445	587 L, 27,000 Vials, w/ Vacuum Insulated Hose	Contact for Quote
YDD-750-445	783 L, 37,800 Vials, w/ Vacuum Insulated Hose	Contact for Quote
YDD-850-465	890 L, 42,900 Vials, w/ Vacuum Insulated Hose	Contact for Quote
YDD-1300-635	1340 L, 58,500 Vials, w/ Vacuum Insulated Hose	Contact for Quote
YDD-1800-635	1880 L, 94,875 Vials, w/ Vacuum Insulated Hose	Contact for Quote

**Don't see what you need?
Contact us for other cryo
storage solutions!**

COLD STORAGE

- Product Advantages
 - Easy to access
 - Brand-new intelligent monitoring system
 - Reduces frosting and freezing
- Key Design Features
 - High degree of safety
 - Ergonomic design

Haier Zhiling LN Container

ITEM #	DESCRIPTION	PRICE
YDD-200-610	200 L Volume, Holds 11,000 or 13,000 Vials	\$31,000.00
YDD-340-790	340 L Volume, Holds 20,350 or 24,050 Vials	\$35,000.00
YDD-550-1000	550 L Volume, Holds 33,550 or 39,650 Vials	\$39,000.00



- LCD touch operation
- USB interface and supports USB data export
- Monitors the temperature and liquid level in real time and displays the expected remaining service time

Haier LN₂ Transport Trolley

ITEM #	DESCRIPTION	PRICE
YDC-3000H	50.98 x 20.59 x 43.11 in. Trolley, -196°C~Ambient Temp Range	\$16,890.00

- Microprocessor control with large LED display
- Hermetically sealed, optimized single stage refrigeration system
- User programmable high and low temperature alarm set points, power failure, sensor error, extreme voltage, hot condenser, and extremely high ambient

Haier -86° ULT Energy Saving Chest Freezer

ITEM #	DESCRIPTION	PRICE
DW-86W420JA	14.8 CuFt Chest Freezer, 65.7 x 33.5 x 41.7 in. 220V/60Hz	\$12,995.00



WATER SYSTEMS

Upgrade your lab water systems with Rephile's innovative and high-quality products, available at MIDSCI. Don't let high prices rain on your parade - choose Rephile with confidence.



The PURIST® Ultrapure Lab Water System produces Type 1 ultrapure water from reverse osmosis, deionized or distilled water.

- Volumetric dispensing from 0.1 L to 99 L
- 2-year warranty
- Simplified consumable replacement – no tools required

Rephile PURIST® UV Ultrapure

ITEM #	DESCRIPTION	PRICE
RS2200QVVK	PURIST® UV Ultrapure Water System Set, Includes Dispenser	\$6,310.00 <i>Includes shipping!</i>

The NuZar Q is a newly developed intelligent ultrapure water system with advanced remote control and diagnosis functions, taking advantage of modern information technology.

It produces Type I ultrapure water from RO, deionized or distilled water.

- Volumetric dispensing from 0.1 L to 99 L
- Simplified consumable replacement – no tools required
- RFID tag technology
- 2-year warranty

Rephile NuZar Q

ITEM #	DESCRIPTION	PRICE
RN0Q0000K	NuZar Q Water System Set	\$7,853.00 <i>Includes shipping!</i>



The Genie PURIST® ultrapure water system produces Type I ultrapure water from RO, deionized or distilled water and is designed for critical laboratory applications.

Water quality, operation parameters, system status, components and peripheral devices are stored and displayed on the large color touch screen monitor. The system can be linked to multiple dispensers via wired or wireless protocols.



- Volumetric dispensing from 0.1 L to 99 L
- Monitors incoming feed water
- 2-year warranty
- Simplified consumable replacement – no tools required
- RFID tag technology

Rephile Genie PURIST®

ITEM #	DESCRIPTION	SALE PRICE
RG0S0000K	Genie PURIST® Lab Water System Set	\$10,290.00 <i>Includes shipping!</i>

BALANCES

- External calibration
- Ideal for bulk chemicals, liquids and animals
- Easy in-lab calibration feature
- RS232 port for data output to a computer or printer
- Readability from 0.001 g to 0.01 g

Accuris™ Precision Balances

ITEM #	DESCRIPTION	PRICE
W3200-120	Precision Balance 120 Max Capacity, 0.001 g Readability	\$1,138.52
W3200-320	Precision Balance 320 Max Capacity, 0.001 g Readability	\$1,220.23
W3200-500	Precision Balance 500 Max Capacity, 0.001 g Readability	\$1,412.90
W3200-1200	Precision Balance 1200 Max Capacity, 0.01 g Readability	\$1,144.34
W3200-3200	Precision Balance 3200 Max Capacity, 0.01 g Readability	\$1,220.23
W3200-5K	Precision Balance 5000 Max Capacity, 0.01 g Readability	\$1,681.47



- Readability from 0.1 to 0.01 g
- Ideal for liquids, animals, general weighing and transportability
- Easy in-lab calibration feature
- RS232 port for data output to a computer or printer

Accuris™ Compact Top Loading Balances

ITEM #	DESCRIPTION	PRICE
W3300-120	Compact Balance, 120 Max Capacity, 0.01 g Readability	\$327.22
W3300-300	Compact Balance, 300 Max Capacity, 0.01 g Readability	\$381.68
W3300-500	Compact Balance, 500 Max Capacity, 0.01 g Readability	\$381.68
W3300-1200	Compact Balance, 1200 Max Capacity, 0.1 g Readability	\$424.12
W3300-5K	Compact Balance, 3200 Max Capacity, 0.1 g Readability	\$424.12
W3300-10K	Compact Balance, 5000 Max Capacity, 0.1 g Readability	\$569.57

- Ideal for reagents and basic analytical balance needs
- Multiple mode options (weighting, count mode, percent)
- RS232 port for data output to a computer or printer

Accuris™ Analytical Balances

ITEM #	DESCRIPTION	PRICE
W3100-120	Analytical Balance 120 Max Capacity, 0.0001 g Readability, External Calibration	\$2,143.02
W3100-210	Analytical Balance 220 Max Capacity, 0.0001 g Readability, External Calibration	\$2,488.30



GLOVES

Non-Sterile, Powder-Free,
Ambidextrous, Finger
Textured, Transblue Color,
2 Mil*, 100% Latex-Free



Aurelia™ Amazing®

ITEM #	SIZE	QTY	PRICE
92885	X-Small	300/bx	\$23.00
92886	Small	300/bx	\$23.00
92887	Medium	300/bx	\$23.00
92888	Large	300/bx	\$23.00
92889	X-Large	270/bx	\$23.00

Non-Sterile, Powder-Free,
Ambidextrous, Finger
Textured, Pink Color,
2.5 Mil*, 100% Latex-Free



Aurelia™ Blush®

ITEM #	SIZE	QTY	PRICE
78885	X-Small	200/bx	\$17.31
78886	Small	200/bx	\$17.31
78887	Medium	200/bx	\$17.31
78888	Large	200/bx	\$17.31
78889	X-Large	200/bx	\$17.31

Non-Sterile, Powder-Free,
Ambidextrous, Finger
Textured, Teal Color,
2.5 mil*, 100% Latex-Free



Aurelia™ Perform®

ITEM #	SIZE	QTY	PRICE
92395	X-Small	200/bx	\$18.00
92396	Small	200/bx	\$18.00
92397	Medium	200/bx	\$18.00
92398	Large	200/bx	\$18.00
92399	X-Large	200/bx	\$18.00

Non-Sterile, Powder-Free,
Ambidextrous, Finger
Textured, White Color,
3.5 mil*, 100% Latex-Free



Aurelia™ Quest®

ITEM #	SIZE	QTY	PRICE
92895	X-Small	100/bx	\$8.65
92896	Small	100/bx	\$8.65
92897	Medium	100/bx	\$8.65
92898	Large	100/bx	\$8.65
92899	X-Large	100/bx	\$8.65

Non-Sterile, Powder-Free,
Ambidextrous, Finger
Textured, Sky Blue Color,
4.5 Mil*, 100% Latex-Free



Aurelia™ Protégé®

ITEM #	SIZE	QTY	PRICE
93995	X-Small	100/bx	\$12.00
93996	Small	100/bx	\$12.00
93997	Medium	100/bx	\$12.00
93998	Large	100/bx	\$12.00
93999	X-Large	100/bx	\$12.00

Non-Sterile, Powder-Free,
Ambidextrous, Finger
Textured, Black Color,
5.5 Mil*, 100% Latex-Free



Aurelia™ Bold®

ITEM #	SIZE	QTY	PRICE
73995	X-Small	100/bx	\$9.40
73996	Small	100/bx	\$9.40
73997	Medium	100/bx	\$9.40
73998	Large	100/bx	\$9.40
73999	X-Large	100/bx	\$9.40
73994	XX-Large	100/bx	\$9.40

*Palm Thickness.

Take advantage of starter kit savings!

- Four AlphaPettes™ (10 µL, 20 µL, 200 µL, 1000µL)
- Three racks of **PR1MA™** Low-Binding Pipette Tips (1 each 10 µL, 200 µL, 1250 µL tips)
- One backpack
- One 6-place carousel pipettor stand



PR1MA™ Starter Kit

ITEM #	DESCRIPTION	PRICE
PR-LQH01-SK	PR1MA™ Starter Kit, Includes, Backpack, A-10, A-20, A-200, A-1000, 3rks PR1MA Tips, 6-pl Carousel*	\$997.00

*No mix and match or switching pipettor or tip sizes.

PR1MA™ Low Retention Graduated Pipette Tips are certified free of detectable RNase, DNase, DNA, PCR inhibitors, and tested pyrogen-free.



Family owned and trusted since 1983



Stocked products ship in 24 hours



Responses within 24 hours



REAL people answering your questions and placing your orders

We're here to help.
Place your order
with us today!

ORDER TODAY! midsci.com • custserv@midsci.com • 800.227.9997 • f: 636.225.9998