

JULY 27 – DECEMBER 31, 2026

Labnet LabValues 2026

Save big with special offers for laboratory equipment and consumables

Labnet™ Enduro™ Gel XL

The compact design of Enduro Gel XL results in a much smaller footprint than that of traditional gel boxes and power supplies. The 150V power supply connects directly into the gel box, eliminating tangled leads.

The molded gel tank provides leak-proof performance while the vented lid aids in dispersing heat. For safety, a magnetic sensor recognizes the presence or absence of the lid and allows current to flow accordingly. Operating parameters are set digitally, including amperage and voltage with automatic crossover. A run can be paused at any time and the parameters changed.

- ▶ Casting gels is quick and easy – place the gel tray in the stand, slide the comb into the slots, and pour
- ▶ Multi-channel compatible combs and UV-transparent gel trays
- ▶ Covered by 1-year warranty



2LV-1.26

Cat. No.	Description	List Price	Special Price
E0160	Enduro Gel XL, with standard casting set (E0168), gel tank with safety lid and power supply, 120V	\$821.41	\$574.99



2LV-2.26

Labnet MPS 1001 Mini Plate Spinner

Centrifuges that accept microplates are typically larger units that take up valuable space on the bench. The unique design of the Labnet MPS 1001 Mini Plate Spinner allows for a compact centrifuge that has an extremely small footprint. Two standard microplates are loaded through a slot in the top of the unit. Sealed microplates are loaded vertically. Surface tension keeps samples in place so there is no danger of spilling.

Upon closing the lid of the unit, the rotor quickly accelerates the microplates to 2,500 rpm. Opening the lid engages the brake and brings the rotor to a stop. For temperature-sensitive applications, the Mini Plate Spinner may be used in a cold room.

- ▶ Capacity for two standard plates
- ▶ Top loading design
- ▶ Exceptionally small footprint
- ▶ Starts and stops in seconds

Cat. No.	Description	List Price	Special Price
C1001	MPS 1001 Mini Plate Spinner, 120V	\$706.47	\$481.89



2LV-3.26

-30%

Labnet™ Mini Microcentrifuge

Simple and convenient to operate, the Labnet Mini Microcentrifuges feature greater capacity (8-place microcentrifuge tubes and 4-place PCR strip rotors), a quick-release rotor system, and electronic brake. Once samples are loaded and the lid is closed, the rotor rapidly accelerates to 6,000 rpm, making it an ideal centrifuge for bringing small droplets to the bottom of tubes.

- ▶ Quick release rotor system – no tools needed to exchange rotors
- ▶ Electronic braking – allows for more rapid deceleration, less waiting

Cat. No.	Description	List Price	Special Price
C1601-G	Mini centrifuge (grey lid), 100-240V	\$230.41	\$161.29
C1601-P	Mini centrifuge (purple lid), 100-240V	\$230.41	\$161.29
C1601-R	Mini centrifuge (red lid), 100-240V	\$230.41	\$161.29

Labnet Spectrafuge™ 6C Centrifuge

The Spectrafuge 6C has been designed for quick production of PPP (platelet-poor plasma) or PRP (platelet-rich plasma) as well as other applications in the clinical and research laboratory.

The fixed angle rotor of the Spectrafuge 6C holds up to 6 x 15 mL conical tubes of maximum tube height of 120 mm, and with adapters will accommodate most common blood draw tubes. Variable speed allows for a variety of applications to be performed.

- ▶ Large digital displays for time and speed
- ▶ Small footprint
- ▶ Cold room compatible

Cat. No.	Description	List Price	Special Price
C0060	Spectrafuge 6C centrifuge, with 6 x 10/15 mL rotor, 120V	\$1,214.28	\$850.00



-30%

2LV-4.26

-40%



2LV-5.26

Labnet Spectrafuge 16M Microcentrifuge

The Spectrafuge 16M is the perfect high speed compact microcentrifuge. Producing a maximum of 16,000 x g, it is the ideal unit for DNA applications. A continuous flow of ambient air through the rotor chamber keeps the rotor and samples cool. A brushless, maintenance-free motor provides extremely fast acceleration and deceleration of the 18-place rotor and makes the Spectrafuge 16M one of the quietest microcentrifuge in its class.

- ▶ Extremely small footprint

Cat. No.	Description	List Price	Special Price
C0160	Spectrafuge 16M microcentrifuge, standard gray, 120V	\$1,499.79	\$899.87
C0160-R	Spectrafuge 16M microcentrifuge, red, 120V	\$1,499.79	\$899.87

-40%



2LV-6.26



Labnet™ Excel™ Electronic Pipettors

The Excel is an ergonomic, intuitive electronic pipettor for use in most laboratory applications. Its high precision stepping motor controls piston action, resulting in exceptional accuracy and precision.

- Factory calibrated to strict ISO8655 standard
- Large LCD screen and simple interface
- Ergonomically designed handle and finger rest

Cat. No.	Description	List Price	Special Price
P3600L-10	Excel single-channel electronic pipettor, 0.5-10 µL	\$406.20	\$243.72
P3600L-20	Excel single-channel electronic pipettor, 2-20 µL	\$406.20	\$243.72
P3600L-200	Excel single-channel electronic pipettor, 20-200 µL	\$406.20	\$243.72
P3600L-1200	Excel single-channel electronic pipettor, 100-1200 µL	\$406.20	\$243.72
P3608L-10	Excel 8-channel electronic pipettor, 0.5-10 µL	\$829.44	\$497.66
P3608L-20	Excel 8-channel electronic pipettor, 2-20 µL	\$829.44	\$497.66
P3608L-200	Excel 8-channel electronic pipettor, 20-200 µL	\$829.44	\$497.66
P3608L-1200	Excel 8-channel electronic pipettor, 100-1200 µL	\$940.27	\$564.16
P3612L-10	Excel 12-channel electronic pipettor, 0.5-10 µL	\$968.18	\$580.91
P3612L-20	Excel 12-channel electronic pipettor, 2-20 µL	\$968.18	\$580.91
P3612L-200	Excel 12-channel electronic pipettor, 20-200 µL	\$968.18	\$580.91
P3612L-1200	Excel 12-channel electronic pipettor, 100-1200 µL	\$1,106.31	\$663.79

Labnet FastPette™ Pro Pipet Controller

The FastPette Pro is a lightweight and ergonomically designed automatic pipet controller ideal for liquid handling with 0.5 to 100 mL glass or plastic serological pipets.

- Triggers convenient for one or two finger control of speed and mode
- LCD display for convenient visualization of speed, modes, and battery status
- Li-ion battery for fast charging and longer life

Cat. No.	Description	List Price	Special Price
P2004	FastPette Pro is supplied with a universal power supply*, an additional autoclavable 0.2 µm filter, and a wall hanger	\$449.12	\$179.65

*Universal power supply comes with US, EU, UK and AU plugs.



2LV-7.26

\$179.65

Special offers valid
July 27 – December 31, 2026