

CRYOGENIC TUBES



**LIMITED
TIME
DEAL!**



**Buy 2 Cases of PR1MA™ Cryo Tubes and
Get 1 CryoFreeze® Rack FREE!**

Buy any 2 cases of PR1MA™ Cryogenic Tubes, Get 1
CryoFreeze® GripRack™ absolutely FREE!

For free item, add item# 5160-0G-F to your order.

PR1MA™ Cryogenic Tubes

Cryo Tubes with Internally-Threaded Caps

ITEM #	DESCRIPTION	QTY	PRICE
PR-CVIT10-SB	1.0 mL, Cryo Tubes, Int. Thread, Star Base, Sterile, Clear	100/pk	\$30.66
PR-CVIT18-RB	1.8 mL, Cryo Tubes, Int. Thread, Round Base, Sterile, Clear	100/pk	\$30.66
PR-CVIT18-SB	1.8 mL, Cryo Tubes, Int. Thread, Star Base, Sterile, Clear	100/pk	\$30.66
PR-CVIT36-RB	3.6 mL, Cryo Tubes, Int. Thread, Round Base, Sterile, Clear	100/pk	\$40.22
PR-CVIT36-SB	3.6 mL, Cryo Tubes, Int. Thread, Star Base, Sterile, Clear	100/pk	\$40.22
PR-CVIT45-RB	4.5 mL, Cryo Tubes, Int. Thread, Round Base, Sterile, Clear	100/pk	\$41.13
PR-CVIT45-SB	4.5 mL, Cryo Tubes, Int. Thread, Star Base, Sterile, Clear	100/pk	\$41.13

Note: not suitable for immersion in liquid nitrogen

Cryo Tubes with Externally-Threaded Caps

ITEM #	DESCRIPTION	QTY	PRICE
PR-CVET10-SB	1.0 mL, Cryo Tubes, Ext. Thread, Star Base, Sterile, Clear	100/pk	\$26.97
PR-CVET18-RB ¹	1.8 mL, Cryo Tubes, Ext. Thread, Round Base, Sterile, Clear	100/pk	\$26.97
PR-CVET18-SB ²	1.8 mL, Cryo Tubes, Ext. Thread, Star Base, Sterile, Clear	100/pk	\$27.56
PR-CVET45-SB	4.5 mL, Cryo Tubes, Ext. Thread, Star Base, Sterile, Clear	100/pk	\$27.56


¹Cap and tube have same outer diameter; ideal for cane storage.

ITEM #	DESCRIPTION	QTY	PRICE
PR-CVCI	Cap Inserts Cryo Tubes, Assorted	500/pk	\$31.65

²Fits all caps except PR-CVET18-RB

Offer cannot be combined with other discounts or institutional pricing. Customer must pay the listed price for the 2 cases to receive the free item. Not valid on prior purchases. All terms and conditions apply. Promo valid until 12/31/26 or while supplies last.

Internally-Threaded Caps



Externally-Threaded Caps



Assorted Colors Include:

Blue Green Orange
Red Violet Yellow



PR-CVCI

Consumables

Equipment & Instruments