

- Polypropylene housing with nylon mesh filter eliminated the need to fabricate filters from gauze, paper, nylon mesh, or other materials
- Extended gripping tab on housing for easy handling
- The elongated lip on the housing extends beyond the edge of the centrifuge tube enabling aseptic handling with forceps
- 5.5 inch pestle length for convenient hand manipulation
- Convex pestle head includes molded textured surface to better manipulate slippery material in the cell strainer



ORDERING INFORMATION

ITEM #	DESCRIPTION	QUANTITY
PR-40BCS-50	40 µm, Blue, Bulk Packed, Sterile	1/pk, 50 pks/cs
PR-40ICS-50	40 µm, Blue, Indiv. Wrapped, Sterile	1/pk, 50 pks/cs
PR-70BCS-50	70 µm, White, Bulk Packed, Sterile	1/pk, 50 pks/cs
PR-70ICS-50	70 µm, White, Indiv. Wrapped, Sterile	1/pk, 50 pks/cs
PR-100BCS-50	100 µm, Yellow, Bulk Packed, Sterile	1/pk, 50 pks/cs
PR-100ICS-50	100 µm, Yellow, Indiv. Wrapped, Sterile	1/pk, 50 pks/cs
PR-SG-PEST-B	5.5 in. Pestle, Blue, Indiv Wrapped, Sterile	1/pk, 50 pks/cs