

PR1MA™ Non-Treated Roller Bottles



- Available in two sizes:
1000 mL Vol with 490 cm² of surface area
and 2000 mL Vol with 850 cm² of surface
area, treated or non-treated versions
- Large, knurled cap for easy grip
- Vented or non-vented cap for open or
closed systems
- Manufactured from medical grade virgin
polystyrene
- Tissue Culture Treated Roller Bottles
have uniform hydrophilic surface ideal for
consistent cell attachment
- Easy-to-read black printed graduations
- Packaged sterile (gamma irradiated) in
individual bags



SPECIFICATIONS

STERILE	Yes
NON-PYROGENIC	Yes
PRIMARY MATERIAL	Bottle: Polystyrene
SECONDARY MATERIAL	Cap: HDPE
GRADUATED	Yes
GRADUATION INCREMENTS (SUBDIVISIONS)	50 mL
TISSUE CULTURE TREATED	No
RNASE/DNASE FREE	Yes

ORDERING INFORMATION

ITEM #	DESCRIPTION	QUANTITY
PR-RBNVC-1K	1,000 mL, Suspension Culture, Printed Graduations, Non-Vented Cap	1/pk, 24/cs
PR-RBNVC-2K	2,000 mL, Suspension Culture, Printed Graduations, Non-Vented Cap	1/pk, 12/cs



ORDER TODAY! midsci.com • custserv@midsci.com • 800.227.9997 • f: 636.225.9998