

- Stronger, heavier construction
- Surface is hydrophobic and does not facilitate cell attachment
- 35 mm dish has an off-set bottom for improved handling
- Lids contain molded spacers for gas exchange
- Ventilation ribs reduce condensation and allow air circulation
- Consistent flatness allows for even media distribution
- Stackable/Slideable – Features partial edge ring on the base which allows a single dish to be slid from a large stack and enhances stability when stacked
- Slippable – Designed for automated filling and streaking systems, features side arrows and ISO Mark Target for automation
- Packaged sterile (gamma irradiated)



ORDERING INFORMATION

ITEM #	DESCRIPTION	QUANTITY
PR-PD-60STK	60 x 15 mm, Sterile	25/pk, 20 pks/cs
PR-PD100SLD	100 x 15 mm, Sterile, Slideable	25/pk, 20 pks/cs
PR-PD-100STK	100 x 15 mm, Sterile, Stackable	25/pk, 20 pks/cs
PR-PD-100SLP	100 x 15 mm, Sterile, Slippable	25/pk, 20 pks/cs
PR-PD-1002	100 x 15 mm, Sterile, 2 Compartment	25/pk, 20 pks/cs
PR-PD-1003	100 x 15 mm, Sterile, 3 Compartment	25/pk, 20 pks/cs
PR-PD-1004	100 x 15 mm, Sterile, 4 Compartment	25/pk, 20 pks/cs