

PR1MA™ Treated Roller Bottles



- Expanded surface roller bottle increases surface area for cell culture
- Large knurled cap for easy grip
- Vented (HRBV) or non-vented cap (HRBNV) for open or closed systems
- Tissue Culture Treated Roller Bottles have uniform hydrophilic surface ideal for consistent cell attachment
- Vent caps (RBCAP) feature 0.22µm hydrophobic membrane
- Easy-to-read black printed graduations
- Manufactured from medical grade virgin polystyrene
- Packaged sterile in individual bags



SPECIFICATIONS

STERILE	Yes
NON-PYROGENIC	Yes
PRIMARY MATERIAL	Bottle: Polystyrene
SECONDARY MATERIAL	Cap: HDPE
GRADUATED	Yes
GRADUATION INCREMENTS (SUBDIVISIONS)	50 mL
TISSUE CULTURE TREATED	No
CELL GROWTH AREA	490 cm ² , 850 cm ² , 1900 cm ²
RNASE/DNASE FREE	Yes

ORDERING INFORMATION

ITEM #	DESCRIPTION	QUANTITY
PR-RBVC-490	490 cm ² , Tissue Culture, Vent Cap, Sterile	1/pk, 24/cs
PR-RBNVC-490	490 cm ² , Tissue Culture, Non-Vent Cap, Sterile	1/pk, 24/cs
PR-RBVC-850-40	850 cm ² , Tissue Culture, Vent Cap, Sterile	1/pk, 12/cs
PR-RBNVC-850	850 cm ² , Tissue Culture, Non-Vent Cap, Sterile	1/pk, 12/cs
PR-RBVC-1900	1900 cm ² , Tissue Culture, Vent Cap, Sterile	1/pk, 12/cs
PR-RBNVC-1900	1900 cm ² , Tissue Culture, Non-Vent Cap, Sterile	1/pk, 12/cs



ORDER TODAY! midsci.com • custserv@midsci.com • 800.227.9997 • f: 636.225.9998