

PR1MA™ Untreated Multi-Well Plates



- Available in 6-, 12-, 24-, 48-, and 96-well format
- Polystyrene, NOT Tissue Culture Treated (suspension culture)
- Plates are labeled with alpha-numeric markings for easy identification of individual wells
- Each well has a raised rim that nests into the rings on the single-direction lid to reduce evaporation and limit the risk of cross-contamination
- Lids have a stacking feature for added stability when handling tall stacks of plates
- Sterile (gamma irradiated)



SPECIFICATIONS

STERILE	Yes
NON-PYROGENIC	Yes
COLOR	Clear
PRIMARY MATERIAL	Polystyrene
TISSUE CULTURE TREATED	No
RNASE/DNASE FREE	Yes

ORDERING INFORMATION

ITEM #	DESCRIPTION	QUANTITY
PR-NTPL-06	6-Well, Non-Treated, w/ Lid, Sterile	1/pk, 100/cs
PR-NTPL-12	12-Well, Non-Treated, w/ Lid, Sterile	1/pk, 100/cs
PR-NTPL-24	24-Well, Non-Treated, w/ Lid, Sterile	1/pk, 100/cs
PR-NTPL-48	48-Well, Non-Treated, w/ Lid, Sterile	1/pk, 100/cs
PR-NTPRBL-96	96-Well, Round Bottom, Non-Treated, w/ Lid, Sterile	1/pk, 100/cs
PR-NTPRB-96	96-Well, Round Bottom, Non-Treated, No Lid, Sterile	1/pk, 100/cs
PR-NTPL-96	96-Well, Non-Treated, w/ Lid, Sterile	1/pk, 100/cs

