Benchmark Scientific | MIDSCI™

Fall into SAVINGS

BLACK FRIDAY SALE

Offers expire 12/04/2026
Use promo code: TURKEY26



PR1MA™ 2D Microplate Centrifuge

FREE PR1MA™ 8 MINI With Purchase!
$347 VALUE

- *Digital control, adjustable rpm and time*
- *Swing out rotor prevents sample spillage (no plate seal required)*
- *Personal sized, 9 x 12" footprint*
- *Accepts all popular types of PCR plates (skirted, semi, and non) and most microtiter plates*

PR-PCRCENT-2D	*PR1MA™ 2D Microplate Centrifuge*	**$939.00**
MF8-B-VF	*PR1MA™ 8 Mini Microcentrifuge, Blue*	**FREE**

Add item at the listed price, and one MF8-B-VF for free item.



MF8-B



ACCURIS INSTRUMENTS COMPACT BALANCE

- *User friendly interface*
- *Easy in-lab calibration feature*
- *RS232 port for data output*
- *500g capacity / 0.01g readability*

W3300-500	*Accuris™ Compact Balance.*	**$412.92**
MF8-B-VF	*PR1MA™ 8 Mini Microcentrifuge, Blue*	**FREE**

Add item at the listed price, and one MF8-B-VF for free item.

BUY ONE GET ONE *FREE!*

PR1MA™ Roto-Mini™ Series

- *Included tube holders:*
  - *11 x 1.5/0.2ml (pack of 2)*
  - *5 x 10-15ml (pack of 2)*
  - *5 x 5-7ml (pack of 2)*
  - *6 x 50ml (pack of 1)*
- *Variable speed control, 5 to 70 rpm*
- *Advanced shaking/mixing modes: Mix (oscillating), Pause, Rock*



RMP2024	*PR1MA™ Roto-Mini™ Series*	**$541.32**
RMP2024-VF	*FREE! PR1MA™ Roto-Mini™ Series*	**FREE**

Add item at the listed price and one RMP2024-VF for free item.



PR1MA™ V2 Vortex Mixer

- *Reversible top platform, cup or flat-dimpled surface*
- *Powerful motor for INSTANT vortexing*
- *Continuous or "touch" operation*
- *Q-Drive™ dynamic balancing system*
- *Variable speed from 200 to 3200 rpm*
- *Wide variety of accessories, for tubes up to 50ml*

BUY ONE GET ONE *FREE!*

VM-1003	*PR1MA™ V2 Vortex Mixer with flip top cup head*	**$342.08**
VM-1003-VF	*FREE! V2 Vortex Mixer*	**FREE**

Add item at the listed price and one VM-1003-VF for free item

ORDER TODAY! midsci.com • custserv@midsci.com • 800.227.9997 • f: 636.225.9998





myFUGE™ 12 Plus

- *COMBI-Rotor for tubes **and** strips*
- *Digital control up to 7,000 rpm / 2,739 xg!*
- *Twice the capacity of traditional mini centrifuges*
- *Near silent operation*
- *Timer up to 99 mins.*


BUY ONE GET ONE FREE!

C2212	*MyFuge™ 12 Plus Digital Mini Centrifuge*	**$708.72**
C2212-VF	*FREE! MyFuge™ 12 Plus*	**FREE**

Add item at the listed price and C2212-VF for free item.

FREE 500g of Agarose!





Includes two casting sets

PR1MA™
Mini Gel Electrophoresis System

- *Complete system, includes gel box, power supply and casting sets*
- *Quick and easy gel casting, no tape, clamps or spacers required*
- *Compact design, footprint: 7 x 5.1 in.*

EMINIGEL	*PR1MA™ Mini Gel Electrophoresis System*	**$691.86**
A1705-VF	*Benchmark Agarose LE, 500g*	**FREE**

Add item EMINIGEL then A1705-VF for free item



A1705





BEADBUG™ 3
MICROTUBE HOMOGENIZER

- *Extremely powerful mixing for lysis, grinding or homogenization*
- *Faster and more effective than tissue grinders*
- *Extremely compact, less than 7" wide*
- *Eliminates the cross contamination associated with handheld homogenizers*
- *Full range of beads in prefilled tubes and bulk available*

D1030	*BeadBug™ 3 Microtube homogenizer, 115V*	**$1,022.95**
D1030-VF	*FREE! BeadBug™ 3*	**FREE**

Add item at the listed price then D1030-VF for free item.



Standard Glass Beads (Acid Washed & Heat Treated) | Triple-Pure™ (Molecular Biology Grade) Zirconium Beads | Stainless Steel Beads (Alcohol Washed & Heat Treated)

A variety of bead sizes to support different sample types and applications are available. For customized recommendations or specific application support, contact our team directly at **tech@midsci.com**.

Buy One Get One Free!

Benchmark Scientific | MIDSCI™

Free offers require purchase of designated items at listed price. ***Offers expire 12/04/2026.*** Not valid on previous purchases. Special price applies only to item specified. Use promo code: ***TURKEY26***

BLACK FRIDAY

30% SALE



SMARTWASHER™ 96

- ***Precision control of aspiration needle position***
  - ***Accurate to 0.1mm***
  - ***Minimal residual liquid, ≤1µl/well***
- ***Auto rinse function eliminates cross contamination***
- ***Up to 100 wash protocols can be customized and stored***
- ***Liquid level sensors in wash and waste bottles***
- ***Wash protocols can be transferred via flash drive***

More Great Deals For Your Lab This Season!

MW9600	*SmartWasher™ 96 Plate Washer*	~~$5609.72~~	**$3926.80**
MW9600-12CM	*12 Channel Manifold for MW9600*	**$404.48**	
MW9600-8CM	*8 Channel Manifold for MW9600*	**$367.08**	

Over 30% OFF!

PR1MA™ Mini HC Shaking Incubator




- ***Heating and Cooling, -8°C from ambient to 50°C***
- ***Extremely compact, just 11" wide***
- ***TOUCH Screen control over time, speed and temperature***
- ***Eco-friendly HFC free cooling***
- ***Instantly exchange flask clamps with MAGic Clamp™ accessories (PATENTED)***

ECO FRIENDLY

PR-ISHC-300	*PR1MA™ Mini HC Shaking Incubator*	~~$3,996.00~~	**$ 2797.20**
IS-PLAT	*MAGic Clamp universal platform, (9.5 x 11.5")*	**$326.70**	
H1000-MR-CMB	*MAGic Clamp™ COMBO, for 125ml, 250ml, 500ml flasks*	**$58.64**	

MC-24R™ Refrigerated Microcentrifuge





COMBI-Rotor (Included)

- ***Included COMBI-Rotor™ accepts tubes and strips (lid included)***
- ***Powerful refrigeration, as low as -20°C***
- ***Touch Operation with speed control to 16,800 xg***
- ***Status indicator light***

C2417-R	*MC-24R™ Refrigerated High Speed Microcentrifuge*	~~$6,087.87~~	**$4261.51**

BOGO

PR1MA™ qMAX™ Green qPCR Mix With Blue Tracking Dye



- ***Superior sensitivity and fast cycling with exceptional results***
- ***Includes Accuris Hot Start Taq for greater specificity and accuracy***
- ***Ideal for low copy number templates***
- ***2X master mix with inert blue tracking dye***

PR2000-N-100	*PR1MA™ qMAX Green qPCR, No ROX, 100 rxns*	**$83.67**	**FOR TWO!**
PR2000-L-100	*PR1MA™ qMAX Green qPCR, Low ROX, 100 rxns*	**$83.67**	**FOR TWO!**
PR2000-H-100	*PR1MA™ qMAX Green qPCR, High ROX, 100 rxns*	**$83.67**	**FOR TWO!**

Add item at the listed price and the same item plus "-VF" for free item